

Cassette Storage

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. 8342..... \$121.00
Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks in 4 drawers, each divided into 12 rows. Cabinets can be stacked up to 8 high. Dimensions: 49.5 x 38 x 21 cm. (19½ x 15 x 8¼")



Tissue Cassette Cabinet, Steel, 4-drawer. 485400 \$315.00

Store up to 2100 cassettes in 6 drawers. Each drawer includes sponge block spacers and an identification label. Cabinets can be stacked up to 4 high. Dimensions: 44.5 x 24.5 x 39.5cm(H) (17½ x 9½ x 15½")




Tissue Cassette Cabinet, ABS. 97-0600 \$162.00



Store 1500 cassettes (250/drawer) or 990 embedding rings (165/drawer) in 6 stackable drawers. High-impact polystyrene resists temperatures from -80°C to +80°C. Drawer dimensions: 40.5 x 23 x 5.1 cm.

Tissue Cassette Drawers, Polystyrene, Pk/6. M495-6..... \$163.00

Brushes




Brushes

Flat Bristles (paint style)
Camel Hair, 1" wide, 6½" total length. **LB-37** \$6.25
Round Bristles, camel hair
3/16" Diameter (artist style), 7½" total length. **LB-31** \$1.80
3/4" Diameter, 8½" total length. **LB-30**..... \$10.45

Freeze Spray

Ultra pure and nonflammable; chills to -62°F.

Freeze Spray, 10 oz Can. 49185..... \$13.00



Microtome Blades

Low Profile, 80 x 8.3 mm, Pk/50. **4689** \$138.00
High Profile, 75.7 x 14 mm, Pk/50. **4685** \$143.00
DURAEDGE®, PTFE coated, **USA**
Low Profile, 77 x 8 mm, Pk/50. **7223** \$64.80
High Profile, 77 x 14 mm, Pk/50. **7310** \$64.80
Encore™, PTFE coated, **USA**
Low Profile, 83 x 8 mm, Pk/50. **7380** \$79.80
High Profile, 83 x 14 mm, Pk/50. **7390** \$79.80

Remove paraffin, ice, etc., safely from your laboratory equipment with these plastic scrapers.



Scrapers, plastic

5 cm (2") Wide Blade. **311-245** \$3.70
10 cm (4") Wide Blade. **311-250** \$4.70
15 cm (6") Wide Blade. **311-255** \$5.70



Tissue Flotation Bath. XH-1001..... \$235.00