

# GENERAL LAB SUPPLIES & EQUIPMENT


## Applicators



### Applicators, 6" long

Nonsterile	
Bg/100. <b>95-8702</b> .....	\$1.55
Bx/1000, 10 bg/bx. <b>95-8702</b> .....	\$12.45
Sterile Pk/1, Bx/100. <b>44024-106-1S</b> .....	\$7.85
Sterile Pk/2, Bx/200. <b>44024-106-2S</b> .....	\$9.80
Plain, wood, Bg/1000. <b>95-8701</b> .....	\$7.55

## Aprons



Protect yourself with this disposable 0.8-mil thick apron.

### Aprons, disposable, polyethylene, 46"(L) x 28"(W)

0.8 mil Thick (0.02 mm), Bx/100. <b>105DA2846</b> .....	\$14.35
1.75 mil Thick (0.044 mm), Bx/100. <b>59010485</b> .....	\$38.00

## Autoclaves



B4000-M



B4000-16


### Autoclaves, benchtop, Benchmark (1-year warranty)

8 L, Bioclave™ Mini. <b>B4000-M</b> .....	\$4,039.00
Chamber Dimensions: 6.7" ID x 12.5" (D)	
18 L, Bioclave™. <b>B4000-18</b> .....	\$5,549.00
Chamber Dimensions: 9" ID x 13.75" (D)	

## Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



HS10044H

<b>Bags, Biohazard, Benchtop, Bx/500. HS10320</b> .....	\$76.00
Biohazard Bag Holder, Benchtop. <b>HS10044H</b> .....	\$23.00
Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.	



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

### Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick

14 x 19" (36 x 48 cm), Bx/200. <b>131641419</b> .....	\$66.00
19 x 23" (48 x 58 cm), Bx/200. <b>131641923</b> .....	\$102.00
25 x 35" (64 x 89 cm), Bx/200. <b>131642535</b> .....	\$243.00

### Bags, Biohazard, Heavy-Duty, 2 mil (0.05 mm) thick

14 x 19" (36 x 48 cm), Bx/200. <b>131651419</b> .....	\$214.00
19 x 23" (48 x 58 cm), Bx/200. <b>131651923</b> .....	\$312.00
25 x 35" (64 x 89 cm), Bx/200. <b>131652535</b> .....	\$530.00
31 x 38" (79 x 96 cm), Bx/200. <b>131653138</b> .....	\$665.00

Safely transport patient specimens with these OSHA and NCCLS compliant bags. Each resealable, 2 mil thick polyethylene bag has an attached document pouch and imprinted biohazard symbol.



<b>Bag, Specimen Transport, 6 x 9", Bx/1000. 4919</b> .....	\$76.30
---	---------



### Bags, Waste Disposal, HDPE, Cleanware™

8 x 10", Pk/50. <b>131741008</b> .....	\$22.70
12 x 16", Pk/50. <b>131741012</b> .....	\$38.10



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

<b>Bag Sealer, electric, 12" (30.5 cm) Wide. HS28272</b> .....	\$200.00
Bag Sealer Replacement Part Kit. <b>HS28280</b> .....	\$31.00