

Beakers



Polypropylene Beakers



Premiu

Beakers

Glass (borosilicate), Griffin low form, single scale, economy

10 mL, no graduations, each. 6102022910	\$1.85
50 mL, each. 7-440050	\$3.10
100 mL, each. 61020229100	\$3.40
150 mL	
Pk/12. 61020229150	\$46.30
Cs/48. 61020229150	\$156.00
250 mL, each. 61020229250	\$3.50
400 mL, each. 61020229400	\$4.40
600 mL, each. 61020229600	\$5.15
800 mL	
Pk/6. 600BG800	\$65.80
Cs/24. 600BG800	\$236.00
1000 mL, each. 600BG1000	\$9.10
1000 mL, Heavy-Duty, each. 600BG1003	\$14.75
2000 mL, each. 600BG2000	\$16.65
4000 mL, each. 600BG4000	\$46.25

Polypropylene

Economy, Triple-Spout

50 mL, 5 mL graduations, Cs/100. 3640	\$25.90
100 mL, 10 mL graduations, Cs/100. 80093	\$39.40
250 mL, 20 mL graduations, Cs/100. 3642	\$58.20
1000 mL, 50 mL graduations, Cs/100. 80098	\$129.00

Premium, Griffin low form, single scale

50 mL, 10 mL graduations, each. 262100000	\$3.20
100 mL, 25 mL graduations, each. 262110000	\$3.05
250 mL, 50 mL graduations, each. 262130000	\$5.40
600 mL, 100 mL graduations, each. 262150000	\$8.50
1000 mL, 200 mL graduations, each. 262160000	\$8.60
2000 mL, 500 mL graduations, each. 262170000	\$23.15
4000 mL, 1000 mL graduations, each. 262190000	\$45.10

Bench Protectors



Bench Liner, Roll

Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

Bench Liners, white

Roll, 50' x 20", Labmat™. 246750000	\$70.40
Sheet, Ahlstrom, 46 x 57 cm (18 x 22½"), Pk/50. 4002-4657	\$86.90

Bottles



Store immersion oil or mounting media in this borosilicate glass bottle with an applicator rod and ground glass cap.

Bottle, Balsam, glass, 45 mL (empty), **USA. 616**\$18.75




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. 116370001	\$19.35
60 mL (2 oz), Bg/12. 116370002	\$22.50
125 mL (4 oz), Bg/12. 116370004	\$37.40
185 mL (6 oz), Bg/12. 116370006	\$37.40

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. 18145	\$7.15
---------------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. 600BDPSC30	\$14.50
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. 3104053D	\$9.95
60 mL (2 oz), Bx/12. 600BDPSC60	\$16.70
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. 3104054D	\$9.85

French Square



Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. 600BDFSQ15	\$17.55
Bottle Dropper Assembly (replacement)	
20/400, 7 x 60 mm Pipet, Pk/12. 3104074	\$10.85
30 mL (1 oz), Bx/12. 600BDFSQ30	\$19.00
Bottle Dropper Assembly (replacement)	
24/400, 7 x 75 mm Pipet, Pk/12. 3104075	\$13.70

Continued