



Boxes, Glass Disposal, cardboard, polyethylene bag included
 Bench Model, 8 x 8 x 10", Pk/6. **246530002** \$90.45
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** \$122.00



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
H-28/334B2	3 ¾ x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 ½ x 1 ¾ x ½"	1	Forceps & Scissors	\$1.80

Brushes



Brushes

Bench, horse hair, 10" long. **BRS-969.10** \$9.20
 Bristle head: 4¾"(L) x 1½"(W); bristle length: 1".
 Flat Bristles, paint style
 Camel Hair, 1" wide, 6½" total length. **LB-37** \$6.25
 Horsehair, heavy-duty bristles, metal handle. **33761880** \$0.85
 Round Bristles, camel hair, 8½" long. **LB-30** \$10.45
 Swivel-Head, black nylon bristles. **LB-14** \$4.30
 Bristle head: 4¼"(L) x 1½"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



Pen Style (scratch), retractable, 11.5 cm (4½") long, Germany
 Brass. **BRS-295.00** \$12.20
 Nylon. **BRS-293.00** \$9.70
 Steel. **BRS-292.00** \$12.80

Burners



Burner, Alcohol, 120 mL, glass body, metal cap. **97-5313** \$6.95
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** \$5.10



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6 1/8")

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8 1/2")



Burners, Bunsen, natural gas, Humboldt Mfg, **USA**

Standard, flame stabilizer. **H-5870** \$29.15
 High Temperature, Meker-type. **H-5550** \$98.80

Continued