



**Boxes, Glass Disposal**, cardboard, polyethylene bag included  
 Bench Model, 8 x 8 x 10", Pk/6. **246530002** ..... \$90.45  
 Floor Model, 12 x 12 x 27", Pk/6. **246530001** ..... \$122.00



**Boxes, Plastic**, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
<b>H-1</b>	1 x 1 x ¼"	-	Reticles, Pointers	\$0.40
<b>H-1/11B</b>	1 x 1 x ¼"	1	Reticles, Pointers	\$0.50
<b>H-54</b>	2 9/16 x 2 9/16 x ½"	-	Filters	\$1.20
<b>H-54/2916G</b>	2 9/16 x 2 9/16 x ½"	1	Filters	\$1.35
<b>H-28/334B2</b>	3 ¾ x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
<b>H-10</b>	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
<b>H-21/512B</b>	5 ½ x 1 ¾ x ½"	1	Forceps & Scissors	\$1.80

## Brushes



### Brushes

Bench, horse hair, 10" long. **BRS-969.10** ..... \$9.20  
 Bristle head: 4¾"(L) x 1½"(W); bristle length: 1".  
 Flat Bristles, paint style  
 Camel Hair, 1" wide, 6½" total length. **LB-37** ..... \$6.25  
 Horsehair, heavy-duty bristles, metal handle. **33761880** .... \$0.85  
 Round Bristles, camel hair, 8½" long. **LB-30** ..... \$10.45  
 Swivel-Head, black nylon bristles. **LB-14** ..... \$4.30  
 Bristle head: 4¼"(L) x 1½"(W). Total length: 12".

Remove debris, rust, or corrosion in hard to reach places.



Pen Style (scratch), retractable, 11.5 cm (4½") long, Germany  
 Brass. **BRS-295.00** ..... \$12.20  
 Nylon. **BRS-293.00** ..... \$9.70  
 Steel. **BRS-292.00** ..... \$12.80

## Burners



**Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313** ..... \$6.95  
 Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310** ..... \$5.10



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6 1/8")

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8 1/2")



**Burners, Bunsen**, natural gas, Humboldt Mfg, **USA**

Standard, flame stabilizer. **H-5870** ..... \$29.15  
 High Temperature, Meker-type. **H-5550** ..... \$98.80

Continued



**Burner Lighter**, replaceable single flint, round file. **SOL-290.00** ... \$3.40  
 Flint(s) (replacement), Pk/5. **SOL-290.01** ..... \$1.85

### Carboys



**Carboys with Spigot**, polypropylene, round  
 10 L (2.5 gal). **118460025** ..... \$72.85  
 20 L (5 gal). **60034052** ..... \$102.00

### Centrifuges



C1012



C1008-C




MLX-110

**Centrifuges, Mini**  
 5,500 RPM (2000 x g), 12-place, MyFuge™12. **C1012** ..... \$408.00  
 (1-year warranty)  
 6,000 RPM (2000 x g), 8-place, MyFuge™Mini. **C1008-\*** \$317.00  
 \*Specify clear (-C), green (-G), or purple (-P) for lid color.  
 (1-year warranty)  
 10,000 RPM (5000 x g), 6-place, blue body. **MLX-110** ..... \$180.00 ↓  
 Holds (16) 0.2 mL tubes/strips, (6) 0.5 mL tubes,  
 or (6) 1.5/2.0 mL tubes with included accessories.  
 (2-year warranty)

### Clamps & Supports

Rectangular  
Steel Base




Triangular  
Cast-Iron Base

#### Support (Ring) Stands, finished metal base/plated support rod

Rectangular Steel Base  
 6 x 4" base, 5/16 x 18" rod. **3104SSS66** ..... \$13.90  
 8 x 5" base, 3/8 x 20" rod. **3104SSS67** ..... \$17.05  
 9 x 6" base, 1/2 x 24" rod. **3104SSS68** ..... \$20.80

Rectangular Cast-Iron Base  
 Provides greater stability than steel base.  
 6 x 4" base, 5/16 x 18" rod. **3104SSC14** ..... \$15.65  
 8 x 5" base, 3/8 x 20" rod. **3104SSC15** ..... \$18.30  
 9 x 6" base, 1/2 x 24" rod. **3104SSC16** ..... \$25.10

Triangular Cast-Iron Base  
 5" Legs, 1/2 x 24" rod. **1257G2324** ..... \$25.00  
 6" Legs, 1/2 x 36" rod. **1257G23** ..... \$30.65




3104BC37




600THCL01



600CLBURDD



97-4023



1257G49

#### Support Clamps

Burette, plastic coated round jaws, 360° rotatable. **3104BC37** \$7.45  
 Burette, Double, deluxe. **600CLBURDD** ..... \$19.75  
 Plastic covered spring jaws on a metal frame securely grip  
 burettes up to 30 mm OD. Fits rods up to 3/4" (19 mm) OD.  
 Hook Collar, nickel plated. **1257G49** ..... \$5.80  
 Fits rods up to 1/2" (13 mm) in diameter.  
 Right-Angle, aluminum. **97-4023** ..... \$6.10  
 Thermometer, aluminum, rotates 360°. **600THCL01** ..... \$8.05  
 Fits rods up to 1/2" (13 mm) in diameter.



#### Support Rings, plated cast-iron

3" OD (2 1/2" ID). **600SRC103** ..... \$7.25  
 4" OD (3 1/2" ID). **600SRC104** ..... \$9.45  
 5" OD (4 1/2" ID). **600SRC105** ..... \$9.95