

## Desiccators



**Non-Vacuum**

Glazed porcelain plate has 5 mm perforations with a 25 mm center hole.



**Vacuum**

### Desiccators

Non-Vacuum, borosilicate glass	
8" (20 cm) ID. <b>600DSGL210</b> .....	\$123.00
Vacuum, polycarbonate top/polypropylene bottom	
6" (15 cm), 13.8 cm plate diameter. <b>60055204</b> .....	\$58.35
8" (20 cm), 19 cm plate diameter . <b>60055205</b> .....	\$79.25
10" (25 cm), 23.8 cm plate diameter . <b>60055206</b> .....	\$100.00

## Drying Racks



188050000



HS23243A

### Drying Racks

Non-Electric	
Bel-Art. <b>188050000</b> .....	\$132.00
Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).	
Heathrow. <b>HS23243A</b> .....	\$127.00
Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).	

## Eyewear

Absorbs  
99.9% UV



**Glasses, Safety, integral side shield. 43574002**..... \$5.60



Flexible frame and adjustable elastic strap provides a tight seal.

**Goggles, Safety, vented (chemical splash). 97-2001**..... \$3.85

## Filter & Blotting Papers



Ahlstrom - Munksjo  
Filter Paper

### Qualitative, Standard Grades (for gravity filtration)

#### Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. <b>6100-0900</b> .....	\$10.50
12.5 cm, Pk/100. <b>6100-1250</b> .....	\$13.30
15 cm, Pk/100. <b>6100-1500</b> .....	\$16.00
20.5 cm, Pk/100. <b>6100-2050</b> .....	\$25.30
24 cm, Pk/100. <b>6100-2400</b> .....	\$35.90

#### Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. <b>6420-0900</b> .....	\$8.20
15 cm, Pk/100. <b>6420-1500</b> .....	\$14.85

#### Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. <b>6010-0425</b> .....	\$6.35
9 cm, Pk/100. <b>6010-0900</b> .....	\$8.45
12.5 cm, Pk/100. <b>6010-1250</b> .....	\$12.25
15 cm, Pk/100. <b>6010-1500</b> .....	\$14.95
18.5 cm, Pk/100. <b>6010-1850</b> .....	\$21.60
20.5 cm, Pk/100. <b>6010-2050</b> .....	\$23.35
33 cm, Pk/100. <b>6010-3300</b> .....	\$63.00

#### Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. <b>6130-0900</b> .....	\$7.20
15 cm, Pk/100. <b>6130-1500</b> .....	\$12.65
20.5 cm, Pk/100. <b>6130-2050</b> .....	\$19.00
33 cm, Pk/100. <b>6130-3300</b> .....	\$42.20

#### Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. <b>6310-2050</b> .....	\$23.50
33 cm, Pk/100. <b>6310-3300</b> .....	\$63.65

#### Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. <b>6150-0900</b> .....	\$7.45
12.5 cm, Pk/100. <b>6150-1250</b> .....	\$10.85
15 cm, Pk/100. <b>6150-1500</b> .....	\$13.05
20.5 cm, Pk/100. <b>6150-2050</b> .....	\$19.55
33 cm, Pk/100. <b>6150-3300</b> .....	\$43.45

#### Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. <b>6170-0900</b> .....	\$6.75
12.5 cm, Pk/50. <b>6170-1250</b> .....	\$9.10
15 cm, Pk/50. <b>6170-1500</b> .....	\$10.60
20.5 cm, Pk/50. <b>6170-2050</b> .....	\$16.15
33 cm, Pk/50. <b>6170-3300</b> .....	\$34.65

### Qualitative, Pleated (Fluted)

#### Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. <b>5420-1250</b> .....	\$72.85
18.5 cm, Pk/100. <b>5420-1850</b> .....	\$82.45

Continued