

Desiccators



Non-Vacuum

Glazed porcelain plate has 5 mm perforations with a 25 mm center hole.



Vacuum

Desiccators

Non-Vacuum, borosilicate glass	
8" (20 cm) ID. 600DSGL210	\$123.00
Vacuum, polycarbonate top/polypropylene bottom	
6" (15 cm), 13.8 cm plate diameter. 60055204	\$58.35
8" (20 cm), 19 cm plate diameter . 60055205	\$79.25
10" (25 cm), 23.8 cm plate diameter . 60055206	\$100.00

Drying Racks



188050000



HS23243A

Drying Racks

Non-Electric	
Bel-Art. 188050000	\$132.00
Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).	
Heathrow. HS23243A	\$127.00
Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).	

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. 43574002..... \$5.60



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). 97-2001..... \$3.85

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. 6100-0900	\$10.50
12.5 cm, Pk/100. 6100-1250	\$13.30
15 cm, Pk/100. 6100-1500	\$16.00
20.5 cm, Pk/100. 6100-2050	\$25.30
24 cm, Pk/100. 6100-2400	\$35.90

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. 6420-0900	\$8.20
15 cm, Pk/100. 6420-1500	\$14.85

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. 6010-0425	\$6.35
9 cm, Pk/100. 6010-0900	\$8.45
12.5 cm, Pk/100. 6010-1250	\$12.25
15 cm, Pk/100. 6010-1500	\$14.95
18.5 cm, Pk/100. 6010-1850	\$21.60
20.5 cm, Pk/100. 6010-2050	\$23.35
33 cm, Pk/100. 6010-3300	\$63.00

Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. 6130-0900	\$7.20
15 cm, Pk/100. 6130-1500	\$12.65
20.5 cm, Pk/100. 6130-2050	\$19.00
33 cm, Pk/100. 6130-3300	\$42.20

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. 6310-2050	\$23.50
33 cm, Pk/100. 6310-3300	\$63.65

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. 6150-0900	\$7.45
12.5 cm, Pk/100. 6150-1250	\$10.85
15 cm, Pk/100. 6150-1500	\$13.05
20.5 cm, Pk/100. 6150-2050	\$19.55
33 cm, Pk/100. 6150-3300	\$43.45

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. 6170-0900	\$6.75
12.5 cm, Pk/50. 6170-1250	\$9.10
15 cm, Pk/50. 6170-1500	\$10.60
20.5 cm, Pk/50. 6170-2050	\$16.15
33 cm, Pk/50. 6170-3300	\$34.65

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. 5420-1250	\$72.85
18.5 cm, Pk/100. 5420-1850	\$82.45

Continued

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200** \$68.60
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100** \$10.20

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900** \$29.45
 15 cm, Pk/100. **0940-1500** \$50.70
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900** \$29.45
 15 cm, Pk/100. **0740-1500** \$50.70
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900** \$33.30
 15 cm, Pk/100. **0550-1500** \$53.45

Blotting/Chromatography/Electrophoresis Papers

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657** \$122.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598, 2316)
 7 x 10 cm, Pk/100. **2388-0710** \$28.30
 20 x 20 cm, Pk/100. **2388-2020** \$63.80
 46 x 57 cm, Pk/100. **2388-4657** \$249.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114** \$27.90
 20 x 20 cm, Pk/100. **2438-2020** \$54.55
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114** \$52.30
 20 x 20 cm, Pk/50. **3208-2020** \$74.15



Smooth tips & polypropylene grips


- Forceps, Filter, 11.5 cm. 300-025** \$4.40




Vinyl-coated lead wire is 300 mm long.

- Flask Weight, lead wire, 114 g. 183240000** \$36.10

Freezer Racks



175HS-9-3



175HD-442



175V-13-2L

Freezer Racks, stainless steel, USA

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

	1-9	≥10
Horizontal, Holds 2" Boxes		
9-Place (3L x 3H). 175HS-9-2	\$59.00	53.00
Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)		
12-Place (3L x 4H). 175HS-12-2	\$78.00	71.00
Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)		
15-Place (3L x 5H). 175HS-15-2	\$98.00	89.00
Dimensions: 16.31(L) x 5.5(W) x 11"(H)		
16-Place (4L x 4H). 175HL-16-2	\$105.00	95.00
Dimensions: 22(L) x 5.5(W) x 8.75"(H)		
20-Place (4L x 5H). 175HL-20-2	\$132.00	119.00
Dimensions: 22(L) x 5.5(W) x 11"(H)		
Horizontal, Holds 3" Boxes		
9-Place (3L x 3H). 175HS-9-3	\$61.00	55.00
Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)		
12-Place (4L x 3H). 175HL-12-3	\$80.00	73.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		



Horizontal Drawer (fits upright freezers)

	1-9	≥10
Horizontal Drawer, stainless steel drawers		
Holds 2" Boxes		
9-Place (3L x 3H). 175HD-332	\$85.00	77.00
Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)		
12-Place (3L x 4H). 175HD-342	\$111.00	101.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
15-Place (3L x 5H). 175HD-352	\$141.00	128.00
Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)		
16-Place (4L x 4H). 175HD-442	\$151.00	137.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		
20-Place (5L x 4H). 175HD-542	\$188.00	171.00
Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)		
Horizontal Drawer		
Holds 5" Boxes		
6-Place (3L x 2H). 175RD-325	\$74.00	67.00
Dimensions: 18.125(L) x 6(W) x 10.75"(H)		
8-Place (4L x 2H). 175RD-425	\$99.00	90.00
Dimensions: 24(L) x 6(W) x 10.75"(H)		
Holds 15 mL tubes		
160-Place (80L x 2H). 175HD-32ML-15	\$99.00	90.00
Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)		
220-Place (110L x 2H). 175HD-42ML-15	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		
Holds 50 mL tubes		
84-Place (42L x 2H). 175HD-42ML-50	\$129.00	118.00
Dimensions: 22(L) x 5.5(W) x 9.5"(H)		

Continued