

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200** \$68.60
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100** \$10.20

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900** \$29.45
 15 cm, Pk/100. **0940-1500** \$50.70
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900** \$29.45
 15 cm, Pk/100. **0740-1500** \$50.70
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900** \$33.30
 15 cm, Pk/100. **0550-1500** \$53.45

Blotting/Chromatography/Electrophoresis Papers

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657** \$122.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598,2316)
 7 x 10 cm, Pk/100. **2388-0710** \$28.30
 20 x 20 cm, Pk/100. **2388-2020** \$63.80
 46 x 57 cm, Pk/100. **2388-4657** \$249.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114** \$27.90
 20 x 20 cm, Pk/100. **2438-2020** \$54.55
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114** \$52.30
 20 x 20 cm, Pk/50. **3208-2020** \$74.15



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025** \$4.40



Vinyl-coated lead wire is 300 mm long.

- Flask Weight, lead wire, 114 g. 183240000** \$36.10

Freezer Racks



175HS-9-3



175HD-442



175V-13-2L

Freezer Racks, stainless steel, USA

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes	1-9	≥10
9-Place (3L x 3H). 175HS-9-2 Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)	\$59.00	53.00
12-Place (3L x 4H). 175HS-12-2 Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)	\$78.00	71.00
15-Place (3L x 5H). 175HS-15-2 Dimensions: 16.31(L) x 5.5(W) x 11"(H)	\$98.00	89.00
16-Place (4L x 4H). 175HL-16-2 Dimensions: 22(L) x 5.5(W) x 8.75"(H)	\$105.00	95.00
20-Place (4L x 5H). 175HL-20-2 Dimensions: 22(L) x 5.5(W) x 11"(H)	\$132.00	119.00
Horizontal, Holds 3" Boxes		
9-Place (3L x 3H). 175HS-9-3 Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)	\$61.00	55.00
12-Place (4L x 3H). 175HL-12-3 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$80.00	73.00



Horizontal Drawer (fits upright freezers)

Horizontal Drawer, stainless steel drawers	1-9	≥10
Holds 2" Boxes		
9-Place (3L x 3H). 175HD-332 Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)	\$85.00	77.00
12-Place (3L x 4H). 175HD-342 Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)	\$111.00	101.00
15-Place (3L x 5H). 175HD-352 Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)	\$141.00	128.00
16-Place (4L x 4H). 175HD-442 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$151.00	137.00
20-Place (5L x 4H). 175HD-542 Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)	\$188.00	171.00
Horizontal Drawer		
Holds 5" Boxes		
6-Place (3L x 2H). 175RD-325 Dimensions: 18.125(L) x 6(W) x 10.75"(H)	\$74.00	67.00
8-Place (4L x 2H). 175RD-425 Dimensions: 24(L) x 6(W) x 10.75"(H)	\$99.00	90.00
Holds 15 mL tubes		
160-Place (80L x 2H). 175HD-32ML-15 Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)	\$99.00	90.00
220-Place (110L x 2H). 175HD-42ML-15 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$129.00	118.00
Holds 50 mL tubes		
84-Place (42L x 2H). 175HD-42ML-50 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$129.00	118.00

Continued