

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200** \$68.60
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100** \$10.20

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900** \$29.45
 15 cm, Pk/100. **0940-1500** \$50.70
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900** \$29.45
 15 cm, Pk/100. **0740-1500** \$50.70
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900** \$33.30
 15 cm, Pk/100. **0550-1500** \$53.45

Blotting/Chromatography/Electrophoresis Papers

- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657** \$122.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598, 2316)
 7 x 10 cm, Pk/100. **2388-0710** \$28.30
 20 x 20 cm, Pk/100. **2388-2020** \$63.80
 46 x 57 cm, Pk/100. **2388-4657** \$249.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114** \$27.90
 20 x 20 cm, Pk/100. **2438-2020** \$54.55
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114** \$52.30
 20 x 20 cm, Pk/50. **3208-2020** \$74.15



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025** \$4.40




Vinyl-coated lead wire is 300 mm long.

- Flask Weight, lead wire, 114 g. 183240000** \$36.10

Freezer Racks



175HS-9-3



175HD-442



175V-13-2L

Freezer Racks, stainless steel, USA

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

Horizontal, Holds 2" Boxes	1-9	≥10
9-Place (3L x 3H). 175HS-9-2 Dimensions: 16.31(L) x 5.5(W) x 6.5"(H)	\$59.00	53.00
12-Place (3L x 4H). 175HS-12-2 Dimensions: 16.31(L) x 5.5(W) x 8.75"(H)	\$78.00	71.00
15-Place (3L x 5H). 175HS-15-2 Dimensions: 16.31(L) x 5.5(W) x 11"(H)	\$98.00	89.00
16-Place (4L x 4H). 175HL-16-2 Dimensions: 22(L) x 5.5(W) x 8.75"(H)	\$105.00	95.00
20-Place (4L x 5H). 175HL-20-2 Dimensions: 22(L) x 5.5(W) x 11"(H)	\$132.00	119.00
Horizontal, Holds 3" Boxes		
9-Place (3L x 3H). 175HS-9-3 Dimensions: 16.31(L) x 5.5(W) x 9.5"(H)	\$61.00	55.00
12-Place (4L x 3H). 175HL-12-3 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$80.00	73.00



Horizontal Drawer (fits upright freezers)

Horizontal Drawer, stainless steel drawers	1-9	≥10
Holds 2" Boxes		
9-Place (3L x 3H). 175HD-332 Dimensions: 16.5(L) x 5.5(W) x 7.125"(H)	\$85.00	77.00
12-Place (3L x 4H). 175HD-342 Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)	\$111.00	101.00
15-Place (3L x 5H). 175HD-352 Dimensions: 16.5(L) x 5.5(W) x 11.812"(H)	\$141.00	128.00
16-Place (4L x 4H). 175HD-442 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$151.00	137.00
20-Place (5L x 4H). 175HD-542 Dimensions: 27.25(L) x 5.5(W) x 9.5"(H)	\$188.00	171.00
Horizontal Drawer		
Holds 5" Boxes		
6-Place (3L x 2H). 175RD-325 Dimensions: 18.125(L) x 6(W) x 10.75"(H)	\$74.00	67.00
8-Place (4L x 2H). 175RD-425 Dimensions: 24(L) x 6(W) x 10.75"(H)	\$99.00	90.00
Holds 15 mL tubes		
160-Place (80L x 2H). 175HD-32ML-15 Dimensions: 16.5(L) x 5.5(W) x 9.5"(H)	\$99.00	90.00
220-Place (110L x 2H). 175HD-42ML-15 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$129.00	118.00
Holds 50 mL tubes		
84-Place (42L x 2H). 175HD-42ML-50 Dimensions: 22(L) x 5.5(W) x 9.5"(H)	\$129.00	118.00

Continued

Horizontal Drawer-Tray (fits upright freezers)


Drawer-Tray, aluminum drawer-trays		
Holds 2" Boxes	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75"(H)		

Vertical (fits chest freezers)

Vertical, Holds 2" Boxes		
13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$36.40	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$36.40	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$38.10	
65 mm Top ID, Pk/6. 600GF614065	\$45.80	
75 mm Top ID, Pk/6. 600GF614075	\$47.25	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$20.75	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$34.05	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$27.90	
35 mm Top ID, 8.1 mL capacity, Pk/12. 146850035	\$25.70	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$59.35	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$15.90	↓
100 mm Top ID, Pk/6. 146600100	\$15.40	↓
145 mm Top ID, Pk/6. 60057224	\$19.60	
Tapered Stem		
101 mm Top ID, each. 146810000	\$30.75	
152 mm Top ID, each. 146820000	\$34.80	

Continued

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art		
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.		
55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75	
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60	
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65	
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60	
Standard Stem, Economy, polypropylene		
50 mm Top ID, each. 310FUNPP50	\$0.80	
62 mm Top ID, each. 310FUNPP62	\$0.95	
75 mm Top ID, each. 310FUNPP75	\$1.15	
100 mm Top ID, each. 310FUNPP100	\$1.65	

Gauze Pads


These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads , 10 pk/200, Cs/2000. 44122000	\$95.00
---	---------

Gloves & Hand Protection

Cotton Gloves



Heat-Resistant Gloves



Gloves

Cotton		
Protect optics and photographs from fingerprints. Can also be used as glove liners.		
Small, Pk/12. GLV-190.10	\$4.95	
Large, Pk/12. GLV-190.20	\$5.25	
Heat-Resistant		
Terry cloth gloves protect hands from high temperatures.		
33 cm (13") Long, Pair. 132010000	\$36.50	
47 cm (18.5") Long, Pair. 132010001	\$49.70	



Safely grip hot or cold labware with this silicone rubber hand protector partially lined with mini suction cups. Withstands temperatures from 260°C (500°F) to -57°C (-70°F).

Hand Protector , red. 380000000	\$18.35
---	---------