


## Homogenizers & Tissue Grinder



Rapidly homogenize, emulsify, suspend and disrupt 0.1 mL to 250 mL samples. Includes a 5 x 50 mm flat-bottom probe & 7 x 50 mm saw tooth probe for use with micro tubes.

**Homogenizer, Handheld. D1000**..... \$1,099.00



Homogenize a variety of sample types, from soft tissue to bone, in 35 seconds or less.  
Benchmark BeadBlaster™ 24

**Homogenizer, Benchtop, Microtube. D2400**..... \$7,549.00

- Holds (24) 2.0 mL tubes, or (6) 5 mL tubes w/optional tube holder
- Speed adjusts from 4.00 to 7.00 m/s in 0.05 m/s increments
- Stores 50 programs in memory, each with a maximum of 10 cycles
- 2-year warranty

### Homogenizer Prefilled Tubes

2.0 mL screw-capped tubes are prefilled with your choice of the following beads:

Glass, acid washed

0.1 mm, Pk/50. **D1031-01**.....\$119.00

0.5 mm, Pk/50. **D1031-05**.....\$119.00


Zirconium Beads, acid washed/heat treated/nuclease-free

1.5 mm, Pk/50. **D1032-15**.....\$166.00

3.0 mm, Pk/50. **D1032-30**.....\$166.00

Stainless Steel Beads, acid washed

2.8 mm, Pk/50. **D1033-28**.....\$166.00



### Tissue Grinder, Dual® Glass Pestle

1 mL Capacity, **LG10640100** .....\$95.00

3 mL Capacity, **LG10640102** .....\$107.00

## Hot Plate, Magnetic Stirrer & Hot Plate Stirrers



310HOTPLT

H4000-H

### Hot Plates, analog

6" Diameter Cast Iron Burner, 1000W. **310HOTPLT** ..... \$29.30

7.5 x 7.5" Ceramic Platform, 80-380°C. **H4000-H**..... \$259.00



S1005

H4000-S



### Magnetic Stirrers, analog, Benchmark

3.4 x 3.4" Platform, 100-1500 rpm, MiniMag™. **S1005** ..... \$122.00  
(1-year warranty)

7.1 x 7.1" Ceramic Platform, 60-1500 rpm. **H4000-S** ..... \$259.00  
(2-year warranty)

**Hot Plate Stirrer, 80-380°C, Benchmark. H4000-HS**..... \$389.00

## Ice Buckets & Pans



Ice Buckets



Ice Pans

Temperature Range: -196 to 100°C

### Ice Buckets with lids, expanded PVC foam, Bel-Art, USA

2.5 L Capacity (Handy)

Black. **188482002** ..... \$54.55

Blue. **188482001** ..... \$54.55

4 L Capacity (N'Icer)

Black. **188484002** ..... \$66.50

Blue. **188484001** ..... \$66.50

### Ice Pans w/o lids, expanded PVC foam, Bel-Art, USA

1 L Capacity (Mini), 15.5 x 15.5 x 8.5 cm(H)

Black. **188481102** ..... \$40.65

Blue. **188481101** ..... \$40.65

Purple. **188481106** ..... \$40.65

4 L Capacity (Midi), 33 x 24.5 x 12.5 cm(H)

Black. **188484102** ..... \$64.80

Blue. **188484101** ..... \$64.80

Purple. **188484106** ..... \$64.80

9 L Capacity (Maxi), 44 x 35 x 13 cm(H)

Black. **188489102** ..... \$76.20

Blue. **188489101** ..... \$76.20

Purple. **188489106** ..... \$76.20