

Stereo Microscopes



420T-430PLL-10

National

420-430PLL-10 Stereoscope, Binocular, 10-40X Zoom \$649.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 10X-40X total magnification; working distance of 79 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (LED with rheostat control and directional beam), bottom only (LED), or both top and bottom
- Post stand with locking support collar allows viewing specimens up to 93 mm tall
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Lifetime mechanical and optical/3-year electrical warranty

420T-430PLL-10 Stereoscope, Trinocular, 10-40X Zoom \$758.00

- Trinocular version of 420-430PLL-10
- 2X SLR adapter & 0.4X C-mount adapter included



425-445PLL-10

National

425-445PLL-10 Stereoscope, Binocular, 7.5-50X Zoom \$914.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-50X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warranty



Meiji Access www.meijitechno.com

Fiber Optic Illuminator

Please call regarding availability of items not listed



Fiber optic light source with dual gooseneck and two focusing lenses is shown. Items are sold separately.

Fiber Optics, Fiber-Lite™, Dolan-Jenner

Light Source, 150W EKE lamp, rheostat. **FL150/115** \$392.00
(MI-150)


Gooseneck, Dual, chrome plated, 23" long. **FL150/05**..... \$230.00
(EEG2823M)

Focusing Lens, each. **FL150/50** \$39.00
(LH-759)

Accessories & Parts

Counting Chambers

Bright-Line®
Hemocytometer



Hemocytometer, bright-lined, V-slash, Germany

Chamber with 2 Coverslips, Marienfeld. **0650030** \$134.00

Hemocytometers, Bright-Line®, molded borosilicate glass, **USA**

Coverslips, 20 x 26 x 0.4 mm, Erie, Pk/2. **2026H** \$11.05

Reduced Thickness/phase counting chamber (3.5 mm thick)

Chamber with 2 Coverslips, Hausser. **1475** \$305.00

Coverslips, 20 x 26 x 0.4 mm, Erie, Pk/2. **2026H** \$11.05

Howard Mold Chamber, Hausser, **USA**

Requires Howard Disc - see RETICLES & STAGE MICROMETERS

Chamber with 2 Coverslips (5080 & 5090). **3820** \$383.00

Coverslip only, 28 x 33 x 0.5 mm (thin), each. **5080** \$46.80

Petroff-Hausser Bacteria Counter, Hausser, **USA**

Coverslip only, 20 x 26 x 0.2 mm, each. **5051**..... \$46.80

England Field Finder, 1 mm thick, **USA**. **KR-834** \$212.00