


glass w/rubber bulb; 75 mm (3")

Medicine (Eye) Dropper, Pk/12. 97-4301 \$3.60



Pipet Fillers


- 3-Valve, rubber, black. **207B** \$15.05
- Wheel-Type, Pipette Pump™
 - Fits Pipets Up To 2 mL, blue. **378970000** \$27.80
 - Fits Pipets Up To 10 mL, green. **378980000** \$29.00
 - Fits Pipets Up To 25 mL, red. **378990000** \$39.55
- Wheel-Type, Pipette Pump™ II, quick-release ea pk/12
 - Fits Pipets Up To 2 mL, blue. **379111002** \$35.15 373.00
 - Fits Pipets Up To 10 mL, green. **379111010** \$35.15 373.00
 - Fits Pipets Up To 25 mL, red. **379111025** \$46.70 496.00



Package includes: Diamond SeroLogic portable pipette controller; two rechargeable batteries; battery charger; three PTFE membrane filters – (1) 0.20 µm & (2) 0.45 µm (one installed); desktop stand and wall mounting bracket.

Pipet Fillers

- Motorized, rechargeable Diamond SeroLogic (1-year warranty). **3380** \$298.00



Dimensions: 29 x 8.8 x 40 cm

Pipet Rack, magnetized, blue, Heathrow. HS20612M \$79.15

Pipets, Mohr (Measuring), Class A, individually certified


- 5 mL, (±0.020 mL), Pk/5. **600PM7058H** \$71.55
- 10 mL, (±0.030 mL), Pk/5. **600PM7058I** \$71.55
- 25 mL, (±0.050 mL), Pk/5. **600PM7058J** \$104.00



Dispense fluids accurately and safely with these color-coded and plugged pipets. Negative graduations provide extra volume (not available on 25 mL). RNase-/DNase-/Pyrogen-free, and noncytotoxic. Calibrated "To Deliver" (TD).

Pipets, Serological, polystyrene, UniPlast™


- Nonsterile
 - 5 mL, blue band, pk/25, Cs/375. **1755** \$109.00
 - 10 mL, orange band, pk/25, Cs/250. **1775** \$76.40
 - 25 mL, red band, pk/10, Cs/100. **1784** \$71.80
- Sterile Pk/1 (individually paper/plastic wrapped)
 - 1 mL, yellow band, Cs/500. **1700** \$109.00
 - 2 mL, green band, Cs/500. **1720** \$117.00
 - 5 mL, blue band, reverse grad., Cs/250. **1740** \$86.55
 - 10 mL, orange band, reverse grad., Cs/200. **1760** \$77.15
 - 25 mL, red band, reverse grad., Cs/100. **1780** \$83.05
 - 50 mL, black band, reverse grad., Cs/100. **1820** \$186.00
- Sterile Bg/25
 - 1 mL, yellow band, Cs/1000. **1710** \$180.00
 - 10 mL, orange band, Cs/250. **1770** \$88.75
 - 25 mL, red band, Cs/100. **1782** \$73.95



Store up to 4 sterilization cans with this epoxy-coated wire rack. (Pipet cans not included)



Pipet Sterilizing Box Rack, Bel-Art. 189420000 \$68.20



- 5 mL Capacity, 1 mL graduated
- 7 mL Capacity, 3 mL graduated

Pipets, Transfer, LDPE (low-density polyethylene)

- Graduated, Nonsterile, 15 cm long
 - 5 mL capacity, 1 mL graduated, 25 drops/mL
 - Bx/500. **PTP-02** \$18.75
 - Cs/10 Bx. **PTP-02** \$178.00
 - 5 mL capacity, thick-walled
 - Cs/10 Bx. **188BP1200** \$147.00
 - 7 mL capacity, 3 mL graduated, 14 drops/mL
 - Bx/500. **PTP-01** \$18.75
 - Cs/10 Bx. **PTP-01** \$178.00
- Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
 - 5 mL capacity, 1 mL graduated, Simport
 - Bx/400. **P200-521S** \$46.75
 - Cs/10 Bx. **P200-521S** \$433.00
- Graduated, Sterile Pk/20, 15 cm long
 - 5 mL capacity, 1 mL graduated, Bx/500. **HS206371B** \$40.30
- Long Stem, Sterile Pk/1 (individually wrapped)
 - 22.5 cm (9") long, 6 mL capacity, Cs/400. **139038** \$76.30
- Long Stem, Nonsterile
 - 30 cm (12") long, 23 mL capacity, Bx/100. **139050** \$30.55

Continued