

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



- Pipets, Volumetric, Class A, individually certified**
- 2 mL (± 0.006 mL), orange band, Pk/5. **600PT7100B** \$75.20
 - 5 mL (± 0.010 mL), white band, Pk/5. **600PT7100C** \$79.00
 - 10 mL (± 0.020 mL), red band, Pk/5. **600PT7100D** \$90.25
 - 20 mL (± 0.030 mL), yellow band, Pk/5. **600PT7100E** \$129.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

- Pipettors, Single Channel**
Fixed Volume, pocket-size, Diamond Jr.
- 5 μ L, red. **3322-5** \$26.35
 - 10 μ L, pink. **3322-10** \$26.35
 - 20 μ L, green. **3322-20** \$26.35
 - 50 μ L, purple. **3322-50** \$26.35
 - 100 μ L, orange. **3322-100** \$26.35



- Adjustable Volume, Diamond Advance™ (3-year warranty)
- 0.5-10 μ L, 0.02 μ L increments, red. **3341-10** \$242.00
 - 5-50 μ L, 0.1 μ L increments, yellow. **3341-50** \$242.00
 - 10-100 μ L, 0.1 μ L increments, blue. **3341-100** \$242.00
 - 100-1000 μ L, 1 μ L increments, orange. **3341-1000** \$242.00




- Pipettor Holders, Carousel, Flip & Grip™, Heathrow**
- Blue. **HS120026** \$83.95
 - Lime. **HS120056** \$83.95



- Pipettor Tips, Filtered (Barrier), sterile** *Inquire*
- 1000 μ L, 35 racks of 96, Cs/3360. **BP71008B** \$247.00
 - Length: 102 mm

- Pipettor Tips, Standard, graduated**
- Nonsterile, natural
 - 0.1-20 μ L, 45 mm, Bg/1000. **151145** \$27.85
 - 50-1250 μ L, 40 racks of 96, Cs/3840. **BP71017** \$197.00
 - Length: 102 mm
 - Sterile, blue
 - 100-1250 μ L, 6 racks of 96, Bx/576. **151153BRS-96** \$56.15



- Reagent Reservoir, 50 mL, v-bottom, Heathrow**
Sterile Bg/1, white polystyrene, Pk/100. **120644** \$99.40

Pitchers




Polypropylene



Stainless Steel

- Pitchers**
- Polypropylene, graduated (mL), tall form, tapered
 - 1000 mL, 50 mL graduations, each. **289910000** \$6.35
 - 2000 mL, 100 mL graduations, each. **289920000** \$7.30
 - 3000 mL, 100 mL graduations, each. **289930000** \$13.80
 - Stainless Steel, graduated (mL/oz), tapered, Polar Ware
 - 1000 mL/32 oz, each. **T1063** \$73.20

Rockers, Shakers & Mixers



- Rockers, Fixed Speed, Benchmark (2-year warranty)**
- 2-D (seesaw), 12 x 12" platform w/flat mat
 - BenchBlotter™. **BR1000** \$499.00
 - Fixed speed of 12 rocks/minute
 - Cold room & incubator safe (4-45°C)
 - Load capacity is 2.3 kg (5 lbs)
 - 3-D (nutating)
 - Small, 10.5 x 7.5" platform w/dimpled mat
 - MiniMixer™. **B3D1020** \$429.00
 - Large, 12 x 12" platform w/dimpled mat
 - BioMixer™. **B3D1320** \$516.00
 - Fixed speed of 24 rocks/minute
 - Platform angle is $\pm 20^\circ$
 - Cold room & incubator safe (4-45°C)
 - Load capacity is 1.8 kg (4 lbs)

Continued