

Made of borosilicate glass in compliance with ASTM specification E-969, Class A requirements. Individually serialized and certified. Calibrated "To Deliver" (TD).



- Pipets, Volumetric, Class A, individually certified**
- 2 mL ( $\pm 0.006$  mL), orange band, Pk/5. **600PT7100B** ..... \$75.20
  - 5 mL ( $\pm 0.010$  mL), white band, Pk/5. **600PT7100C** ..... \$79.00
  - 10 mL ( $\pm 0.020$  mL), red band, Pk/5. **600PT7100D** ..... \$90.25
  - 20 mL ( $\pm 0.030$  mL), yellow band, Pk/5. **600PT7100E** ..... \$129.00



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

- Pipettors, Single Channel**
- Fixed Volume, pocket-size, Diamond Jr.
- 5  $\mu$ L, red. **3322-5** ..... \$26.35
  - 10  $\mu$ L, pink. **3322-10** ..... \$26.35
  - 20  $\mu$ L, green. **3322-20** ..... \$26.35
  - 50  $\mu$ L, purple. **3322-50** ..... \$26.35
  - 100  $\mu$ L, orange. **3322-100** ..... \$26.35



- Adjustable Volume, Diamond Advance™ (3-year warranty)
- 0.5-10  $\mu$ L, 0.02  $\mu$ L increments, red. **3341-10** ..... \$242.00
  - 5-50  $\mu$ L, 0.1  $\mu$ L increments, yellow. **3341-50** ..... \$242.00
  - 10-100  $\mu$ L, 0.1  $\mu$ L increments, blue. **3341-100** ..... \$242.00
  - 100-1000  $\mu$ L, 1  $\mu$ L increments, orange. **3341-1000** ..... \$242.00




- Pipettor Holders, Carousel, Flip & Grip™, Heathrow**
- Blue. **HS120026** ..... \$83.95
  - Lime. **HS120056** ..... \$83.95




- Pipettor Tips, Filtered (Barrier), sterile** ..... *Inquire*
- 1000  $\mu$ L, 35 racks of 96, Cs/3360. **BP71008B** ..... \$247.00
  - Length: 102 mm

- Pipettor Tips, Standard, graduated**
- Nonsterile, natural
- 0.1-20  $\mu$ L, 45 mm, Bg/1000. **151145** ..... \$27.85
  - 50-1250  $\mu$ L, 40 racks of 96, Cs/3840. **BP71017** ..... \$197.00
  - Length: 102 mm
- Sterile, blue
- 100-1250  $\mu$ L, 6 racks of 96, Bx/576. **151153BRS-96** ..... \$56.15



- Reagent Reservoir, 50 mL, v-bottom, Heathrow**
- Sterile Bg/1, white polystyrene, Pk/100. **120644** ..... \$99.40

## Pitchers




Polypropylene



Stainless Steel

- Pitchers**
- Polypropylene, graduated (mL), tall form, tapered
- 1000 mL, 50 mL graduations, each. **289910000** ..... \$6.35
  - 2000 mL, 100 mL graduations, each. **289920000** ..... \$7.30
  - 3000 mL, 100 mL graduations, each. **289930000** ..... \$13.80
- Stainless Steel, graduated (mL/oz), tapered, Polar Ware
- 1000 mL/32 oz, each. **T1063** ..... \$73.20

## Rockers, Shakers & Mixers



- Rockers, Fixed Speed, Benchmark (2-year warranty)**
- 2-D (seesaw), 12 x 12" platform w/flat mat
- BenchBlotter™. **BR1000** ..... \$499.00
  - Fixed speed of 12 rocks/minute
  - Cold room & incubator safe (4-45°C)
  - Load capacity is 2.3 kg (5 lbs)
- 3-D (nutating)
- Small, 10.5 x 7.5" platform w/dimpled mat
- MiniMixer™. **B3D1020** ..... \$429.00
- Large, 12 x 12" platform w/dimpled mat
- BioMixer™. **B3D1320** ..... \$516.00
  - Fixed speed of 24 rocks/minute
  - Platform angle is  $\pm 20^\circ$
  - Cold room & incubator safe (4-45°C)
  - Load capacity is 1.8 kg (4 lbs)

Continued



### Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000**.....\$776.00

Stacking Platform. **BR2000-STACK**.....\$135.00

14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**..... \$48.40

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300**.....\$939.00

Stacking Platform. **BR1000-STACK**.....\$124.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**..... \$48.40

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



### Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**.....\$603.00

Double Platform, 3/4" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**.....\$647.00

Double Platform, 3/4" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



**Rocker, Tube**, 12 x 4.5" grooved mat, Benchmark. **M2100** .....\$379.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs )
- 2-year warranty

### Orbital Shakers



BT30

Aerate and mix with these low to high capacity orbital shakers.



BT302

### Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**..... \$856.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302** .....\$1,009.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



BT3001

Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**.....\$1,619.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

*Optional Magic Clamp™ platform & accessories sold separately - specify configuration*

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**... \$286.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** ..... \$40.30

500 mL Flask Magnetic Clamp. **H1000-MR-500** ..... \$40.30

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** ..... \$61.60