



Rockers, Variable Speed

Benchmark (2-year warranty)

2-D (seesaw), 14 x 12" platform w/flat mat

BenchRocker™ 2D. **BR2000**.....\$776.00

Stacking Platform. **BR2000-STACK**.....\$135.00

14 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**..... \$48.40

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

3-D (nutating), 12 x 12" platform w/dimpled and flat mats

BenchRocker™ 3D. **B3D2300**.....\$939.00

Stacking Platform. **BR1000-STACK**.....\$124.00

12 x 12" platform w/flat mat has 3" of clearance.

Stacking Platform Extenders, Pk/4. **BR2000-SP**..... \$48.40

Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

2-D (seesaw)

Standard, 11 x 14" platform w/rubber mat

Single Platform. **55S1114**.....\$603.00

Double Platform, 3/4" apart. **55D1114**.....\$665.00

Large, 12 x 16" platform w/rubber mat

Single Platform. **55S1216**.....\$647.00

Double Platform, 3/4" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100**\$379.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.



Orbital Shakers, Benchmark (2-year warranty)

Low Speed, 14 x 12" platform w/flat rubber mat

Orbi-Blotter™. **BT30**..... \$856.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

Small, 9.5 x 8" platform w/flat rubber mat

Orbi-Shaker™ JR. **BT302**\$1,009.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



Large, 13 x 12" platform w/flat rubber mat

Orbi-Shaker™. **BT3001**.....\$1,619.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**... \$286.00

250 mL Flask Magnetic Clamp. **H1000-MR-250** \$40.30

500 mL Flask Magnetic Clamp. **H1000-MR-500** \$40.30

1000 mL Flask Magnetic Clamp. **H1000-MR-1000** \$61.60

Continued