


- Tube Rotator**, programmable, Benchmark (2-year warranty)
- Roto-Bot Mini. **R4040** \$599.00
- Roto-Bot. **R4045** \$899.00
- Speed adjusts from 1-80 rpm
 - Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
 - Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



- Vortex Mixers**, Benchmark (2-year warranty)
- MiniMixer, fixed 3000 rpm, blue cup. **VM-03UA** \$138.00
- Vornado™, mini vortexer, fixed 2800 rpm. **BV101**__* \$209.00
- *Specify Gray, Blue or Purple cup head
- BenchMixer™, single cup head, 200-3200 rpm. **BV1000** \$269.00
- Flat Head Attachment, 3". **BV1000-FLAT** \$39.40

Rulers, Meter Stick & Tape Measure



- Rulers**, 15 cm/6" long, scaled in both metric and English units
- Vinyl, flexible
- Transparent. **4461** \$0.95
- White. **4460** \$0.55
- Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065** \$5.70

- Meter Stick**, lacquered wood, plain ends. **605-050** \$4.90
- One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

- Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036** \$3.70
- One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



- Sharps Containers**, red polypropylene
- 1 Quart, translucent locking lid. **440065** \$3.80
- Dimensions: 4.25" x 4.25" x 6.25"(H)
- 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269** \$9.90
- Dimensions: 12" x 4.75" x 11"(H)
- 12 Gallon, hinged vertical entry lid, each. **440321391** \$28.95
- Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers



- Specimen Containers**, polypropylene (PP)
- Screw-Cap
- 90 mL, sterile, Cs/400. **840175** \$110.00
- 480 mL (16 oz), Cs/24. **59014510** \$71.95
- 83 mm (3 1/4") opening
- 960 mL (32 oz), Cs/12. **59014511** \$59.05
- 111 mm (4 3/8") opening
- Snap-Cap, white HDPE bottom/white LDPE lid
- 480 mL (16 oz), Cs/100. **271016W** \$72.50
- 9.9 cm opening x 11.1 cm high
- 1000 mL (34 oz), Cs/100. **271032W** \$109.00
- 11.7 cm opening x 13.7 cm high



- Spot Plates**, glazed white porcelain, 21 mm diameter wells
- 6-Well, each. **LPC-161** \$2.35
- 12-Well, 115 x 95 x 7 mm, each. **3104187701** \$4.15