

Tape

Kraft paper with rubber resin adhesive.



Tape, Label, Write-On™, Bel-Art


1" Core, 15 yd (540") roll

1/2" Wide, each	\$6.00
Green. 134623005	Yellow. 134622005
3/4" Wide, Pk/4	\$25.90
Green. 134820075	Yellow. 134810075
White. 134800075	

3" Core, 40 yd roll

1/2" Wide, each	\$10.95
Green. 134633005	White. 134630005
3/4" Wide, each	\$14.55
Green. 134633075	Yellow. 134632075
3/4" Wide, White, Pk/4. 134830075	\$58.10
3/4" Wide, Mixed, Pk/6. 134630600	\$88.35

blue, green, orange, red, white & yellow




Label Tape Dispenser, multi-roll. 134730000..... \$46.40

Thermal Cycler



Versatile multi-format block holds (96) 0.2 mL tubes/strips, (39) 0.5 mL tubes, or (1) nonskirted 96-well PCR plate. In-situ adapter sold separately. Touchscreen control and height adjustable heated lid. 2-year warranty.



Thermal Cycler, multi-format block, Benchmark. T5000-96 .. \$4,874.00
In-Situ Adapter, holds 4 slides. T5000-A-IS\$489.00

Trays



Trays

Cafeteria Style, polypropylene, 12 x 16 1/4", NSF, **USA**

Orange. 105CT1216O	\$5.95
Red. 105CT1216R	\$5.95
Black. 105CT1216B	\$5.95



Contain spills or transport lab products with this rimmed, durable, large white tray. Holds 5 gallons (18.9 L). *Not autoclavable.*

HDPE (high-density polyethylene), **USA. 1053506**..... \$35.50
 Dimensions: 26 x 18 x 3 1/2" (66 x 45.7 x 8.9 cm)



Polypropylene, autoclavable

12 x 7 3/4 x 5 1/8", Bel-Art, each. 162600000	\$33.70
15 x 14 x 5", Pk/6. 60081721	\$137.50
20 x 17 x 5", Pk/6. 60081722	\$144.50