



Tubes & Vials



- Tubes, Centrifuge, Conical** bottom, sterile, white writing patch
RNase-/DNase-/Pyrogen-free
- 15 mL, PP (polypropylene) tube/HDPE cap, natural
Bag/25, 20 bags, Cs/500. **BP1502** \$90.10
Maximum RCF is 12,000 x g. Frosted writing patch.
 - Rack/25, 20 racks, Cs/500. **BP1501** \$98.45
Maximum RCF is 12,000 x g. Frosted writing patch.
 - 15 mL, PS (polystyrene)
Rack/50, 10 racks, Cs/500. **6274**  \$220.00
Maximum RCF is 3,000 x g.
 - 50 mL, natural, polypropylene (PP) tube/HDPE cap
Bag/25, 20 bags, Cs/500. **BP5002** \$116.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - Bag/25, 20 bags, Cs/500. **6297** \$191.00
Maximum RCF is 20,000 x g
 - Rack/25, 20 racks, Cs/500. **BP5001** \$124.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - 50 mL, Black, polypropylene (PP) tube/HDPE cap
Bag/50, 10 bags, Cs/500. **HS4429** \$230.00
Maximum RCF is 9,500 x g




Self-
Standing



- Tubes, Centrifuge, Self-Standing** (skirted)
- 30 mL
Nonsterile, PP tube/PP cap, Cs/500. **109220** \$143.00
Dimensions: 23(OD) x 90 mm
 - Sterile, PP tube/HDPE cap, Cs/500. **C572-1** \$145.00
Molded-in graduations from 10 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with
dual-position caps.



- Tubes, Culture,**
- 10 x 75 mm (3 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1503** 58.10
 - 12 x 75 mm (5 mL)
Polystyrene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110428** \$76.40
Polystyrene, Nonsterile, Cs/1000. **110410** \$53.70
Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110438** \$76.60
 - 13 x 100 mm (10 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510** \$81.65
 - 17 x 100 mm (16 mL)
Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110178** \$121.00
 - 16 x 125 mm (19 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515** \$135.00
 - 16 x 150 mm (23 mL)
 Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517** \$132.00



Natural



Black

RNase/DNase-free

- Tubes, Microcentrifuge, Economy**, polypropylene
- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** \$35.00
 - 1.5 mL, 0.25 mL graduations, Heathrow
Natural, max RCF is 16,400 x g
Pk/500. **HS4323** \$17.00
Pk/1000. **HS4325** \$33.00
Black, max RCF is 18,000 x g, Pk/500. **HS4323K** \$19.00
 - 5 mL, natural, max RCF is 25,000 x g, Benchmark
Nonsterile, 4 bags/50, Pk/200. **C1005T5** \$62.95
Sterile, 10 bags/20, Pk/200. **C1005T5ST** \$69.95