



**Tubes & Vials**



- Tubes, Centrifuge, Conical** bottom, sterile, white writing patch  
RNase-/DNase-/Pyrogen-free
- 15 mL, PP (polypropylene) tube/HDPE cap, natural  
Bag/25, 20 bags, Cs/500. **BP1502** ..... \$90.10  
Maximum RCF is 12,000 x g. Frosted writing patch.
  - Rack/25, 20 racks, Cs/500. **BP1501** ..... \$98.45  
Maximum RCF is 12,000 x g. Frosted writing patch.
  - 15 mL, PS (polystyrene)  
Rack/50, 10 racks, Cs/500. **6274** .....  \$220.00  
Maximum RCF is 3,000 x g.
  - 50 mL, natural, polypropylene (PP) tube/HDPE cap  
Bag/25, 20 bags, Cs/500. **BP5002** ..... \$116.00  
Maximum RCF is 12,000 x g. Frosted writing patch.
  - Bag/25, 20 bags, Cs/500. **6297** ..... \$191.00  
Maximum RCF is 20,000 x g
  - Rack/25, 20 racks, Cs/500. **BP5001** ..... \$124.00  
Maximum RCF is 12,000 x g. Frosted writing patch.
  - 50 mL, Black, polypropylene (PP) tube/HDPE cap  
Bag/50, 10 bags, Cs/500. **HS4429** ..... \$230.00  
Maximum RCF is 9,500 x g




Self-  
Standing



- Tubes, Centrifuge, Self-Standing** (skirted)
- 30 mL  
Nonsterile, PP tube/PP cap, Cs/500. **109220** ..... \$143.00  
Dimensions: 23(OD) x 90 mm
  - Sterile, PP tube/HDPE cap, Cs/500. **C572-1** ..... \$145.00  
Molded-in graduations from 10 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with  
dual-position caps.



- Tubes, Culture,**
- 10 x 75 mm (3 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1503** ..... 58.10
  - 12 x 75 mm (5 mL)  
Polystyrene, Sterile, dual-position caps  
Cs/500, 20 sterile bags/25. **110428** ..... \$76.40  
Polystyrene, Nonsterile, Cs/1000. **110410** ..... \$53.70  
Polypropylene, Sterile, dual-position caps  
Cs/500, 20 sterile bags/25. **110438** ..... \$76.60
  - 13 x 100 mm (10 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1510** ..... \$81.65
  - 17 x 100 mm (16 mL)  
Polypropylene, Sterile, dual-position caps  
Cs/500, 20 sterile bags/25. **110178** ..... \$121.00
  - 16 x 125 mm (19 mL)  
Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1515** ..... \$135.00
  - 16 x 150 mm (23 mL)  
 Glass (Type 1, Class A borosilicate), rimless  
Cs/1000. **1517** ..... \$132.00




Natural

Black

RNase/DNase-free



- Tubes, Microcentrifuge, Economy**, polypropylene
- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422** ..... \$35.00
  - 1.5 mL, 0.25 mL graduations, Heathrow  
Natural, max RCF is 16,400 x g  
Pk/500. **HS4323** ..... \$17.00  
Pk/1000. **HS4325** ..... \$33.00
  - Black, max RCF is 18,000 x g, Pk/500. **HS4323K** ..... \$19.00
  - 5 mL, natural, max RCF is 25,000 x g, Benchmark  
Nonsterile, 4 bags/50, Pk/200. **C1005T5** ..... \$62.95  
Sterile, 10 bags/20, Pk/200. **C1005T5ST** ..... \$69.95



Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

**Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA**

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500 ..... \$24.40  
Amber. **4048** Lavender. **4047** Natural. **4040**  
Mixed (blue, green, orange, pink, & yellow). **4049** ..... \$29.40
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500..... \$24.40  
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**  
Mixed (blue, green, orange, pink, & yellow). **4039** ..... \$29.40
- 2.0 mL, 0.5 mL graduations, frosted area on side, Pk/500..... \$26.20  
Amber. **4058** Blue. **4053** Natural. **4050**  
Mixed (blue, green, orange, pink, & yellow). **4059** ..... \$31.20

**Tubes, Microcentrifuge, Siliconized, G-Tube™, Bio Plas, USA**

- 1.5 mL, natural color, Pk/500. **4165SL**..... \$48.95  
*Siliconized products should NOT be autoclaved*



**Tubes, PCR, 0.2 mL thin-wall, polypropylene**

- Individual, flat caps, Simport  
Blue, Pk/1000. **T325-2VB** ..... \$47.10  
Natural, Pk/1000. **T325-2VN** ..... \$47.10
- Strips, Amplitube™, Low Profile, individually attached caps  
8-Tube, Pk/125 (1000 tubes). **T320-2LPN** ..... \$96.65  
Holds 0.1 mL. Tubes have flat caps & natural color.
- Strips, Bio Plas, *Caps are sold separately; available in natural only.*  
8-Tube, Pk/125 (1000 tubes)..... \$71.80  
Blue. **5010-4** Natural. **5010-1**  
8-Tube Cap Strips, domed, Pk/300 (2400 caps). **5030-1** .. \$79.20  
12-Tube, Pk/100 (1200 tubes)..... \$84.85  
Blue. **5020-4** Natural. **5020-1**  
12-Tube Cap Strips, domed, Pk/200 (2400 caps). **5040-1** \$79.20



- Vials, Cryogenic, round bottom, internal thread**  
4 mL, Cs/1000. **T301-4**..... \$287.00



Vinyl grip, **USA**


- Vial Decapping Pliers, for 13 mm Crimp Seals. D-13** ..... \$72.00

**Tube Storage**



**Baskets, test tube, aluminum, Heathrow**

- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A** ..... \$27.00
- Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B** ..... \$30.00



Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.

**Tube Cooler, benchtop, 0°C, CoolCube™, Benchmark**

- Microtube and PCR Plate. **R1000** ..... \$132.00

**Cryogenic Storage Boxes**



SB2CC-64 Cryogenic Tube Storage Box

**Tube Storage Boxes, Cardboard, vented**

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC** ..... \$1.65
- 2" Box w/dividers, white  
100-Cell (10 x 10), Pk/12. **SB2CC-100** ..... \$26.65
- ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81** ..... \$26.65
- 64-Cell (8 x 8), Pk/12. **SB2CC-64** ..... \$26.65 ↓
- 2" Mini box w/dividers, 3 x 3 x 2"  
25-Cell (5 x 5), Pk/12. **MB2C-25** ..... \$29.55  
Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C** ..... \$1.65 ↓
- 3" Box w/dividers  
100-Cell (10 x 10), Pk/12. **SB3C-100** ..... \$26.65 ↓
- 81-Cell (9 x 9), Pk/12. **SB3C-81** ..... \$26.65 ↓
- 64-Cell (8 x 8), Pk/12. **SB3C-64** ..... \$26.65 ↓
- Dividers only, Fit 2 & 3" Boxes  
100-Cell (10 x 10), Pk/12. **CD-100** ..... \$11.55  
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
- 81-Cell (9 x 9), Pk/12. **CD-81** ..... \$11.55  
Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 64-Cell (8 x 8), Pk/12. **CD-64** ..... \$11.55  
Holds 14 mm vials.
- 49-Cell (7 x 7), Pk/12. **CD-49** ..... \$11.55  
Holds 14 mm vials.
- 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25** ..... \$11.55
- 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16** ..... \$11.55

Continued