

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500 \$24.40
Amber. **4048** Lavender. **4047** Natural. **4040**
Mixed (blue, green, orange, pink, & yellow). **4049** \$29.40
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500..... \$24.40
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**
Mixed (blue, green, orange, pink, & yellow). **4039** \$29.40
- 2.0 mL, 0.5 mL graduations, frosted area on side, Pk/500..... \$26.20
Amber. **4058** Blue. **4053** Natural. **4050**
Mixed (blue, green, orange, pink, & yellow). **4059** \$31.20

Tubes, Microcentrifuge, Siliconized, G-Tube™, Bio Plas, USA

- 1.5 mL, natural color, Pk/500. **4165SL**..... \$48.95
Siliconized products should NOT be autoclaved



Tubes, PCR, 0.2 mL thin-wall, polypropylene

- Individual, flat caps, Simpourt
Blue, Pk/1000. **T325-2VB** \$47.10
Natural, Pk/1000. **T325-2VN** \$47.10
- Strips, Amplitube™, Low Profile, individually attached caps
8-Tube, Pk/125 (1000 tubes). **T320-2LPN** \$96.65
Holds 0.1 mL. Tubes have flat caps & natural color.
- Strips, Bio Plas, *Caps are sold separately; available in natural only.*
8-Tube, Pk/125 (1000 tubes)..... \$71.80
Blue. **5010-4** Natural. **5010-1**
8-Tube Cap Strips, domed, Pk/300 (2400 caps). **5030-1** .. \$79.20
12-Tube, Pk/100 (1200 tubes)..... \$84.85
Blue. **5020-4** Natural. **5020-1**
12-Tube Cap Strips, domed, Pk/200 (2400 caps). **5040-1** \$79.20



- Vials, Cryogenic, round bottom, internal thread**
4 mL, Cs/1000. **T301-4**..... \$287.00



Vinyl grip, **USA**


- Vial Decapping Pliers, for 13 mm Crimp Seals. D-13** \$72.00

Tube Storage



Baskets, test tube, aluminum, Heathrow

- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A** \$27.00
- Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B** \$30.00



Each cube holds (36) 1.5/2.0 mL tubes, or (96) 0.2 mL tubes, or (12) 0.2 mL strips, or (1) 96-well PCR plate and (10) 1.5/2.0 mL tubes. Maintains 0°C for 4 hours.

Tube Cooler, benchtop, 0°C, CoolCube™, Benchmark

- Microtube and PCR Plate. **R1000** \$132.00

Cryogenic Storage Boxes



SB2CC-64 Cryogenic Tube Storage Box

Tube Storage Boxes, Cardboard, vented

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC** \$1.65
- 2" Box w/dividers, white
100-Cell (10 x 10), Pk/12. **SB2CC-100** \$26.65
- ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81** \$26.65
- 64-Cell (8 x 8), Pk/12. **SB2CC-64** \$26.65 ↓
- 2" Mini box w/dividers, 3 x 3 x 2"
25-Cell (5 x 5), Pk/12. **MB2C-25** \$29.55
Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C** \$1.65 ↓
- 3" Box w/dividers
100-Cell (10 x 10), Pk/12. **SB3C-100** \$26.65 ↓
- 81-Cell (9 x 9), Pk/12. **SB3C-81** \$26.65 ↓
- 64-Cell (8 x 8), Pk/12. **SB3C-64** \$26.65 ↓
- Dividers only, Fit 2 & 3" Boxes
100-Cell (10 x 10), Pk/12. **CD-100** \$11.55
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
- 81-Cell (9 x 9), Pk/12. **CD-81** \$11.55
Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 64-Cell (8 x 8), Pk/12. **CD-64** \$11.55
Holds 14 mm vials.
- 49-Cell (7 x 7), Pk/12. **CD-49** \$11.55
Holds 14 mm vials.
- 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25** \$11.55
- 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16** \$11.55

Continued

5" Box Only, 5 ³ / ₄ x 5 ³ / ₄ x 5", each. SB4C	\$3.20
5" Box w/dividers	
64-Cell (8 x 8), Pk/12. SB4C-64	\$42.45
36-Cell (6 x 6), 15 mL tubes, Pk/12. SB4C-36	\$42.45
16-Cell (4 x 4), 50 mL tubes, Pk/12. SB4C-16	\$42.45
Dividers only, Fit 5" Box	
64-Cell (8 x 8), Pk/12. CD4-64	\$11.55
36-Cell (6 x 6), Pk/12. CD4-36	\$11.55
16-Cell (4 x 4), Pk/12. CD4-16	\$11.55



Holds 1.0-2.0 mL tubes and cryovials. Features a locking hinged lid, and base with molded alphanumeric coordinates. Dimensions: 12.8 x 12.8 x 5.4 cm.

Tube Storage Boxes, Polypropylene

81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. 60066501	\$28.00
---	---------

Holds 1.5-2.0 mL tubes in a two-piece rack.



Tube Storage Boxes, Polystyrene Foam (styrofoam), USA


50-Well, Pk/10. HS4283P	\$44.50
100-Well, Pk/10. HS4282	\$55.65

Tube Racks




Tube Racks, Grip, polypropylene, No-Wire™, Bel-Art

15-16 mm Tubes, 60-place, each. 187490001	\$28.75
18-20 mm Tubes, 40-place, each. 187490002	\$28.75



Tube Racks, Interlocking, polypropylene, Heathrow

Cube, mixed colors, Pk/5. HS29050A	\$44.45
Holds (4) 50 mL tubes, (10) 15 mL tubes, (12) 12 x 75 or 12 x 100 tubes, or (16) 1.5-2.0 mL tubes.	
4-Way, mixed colors, Pk/5. HS29022G	\$40.85
Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes.	
4-Way Rotatable, Rota-Rack™, four connected modules	
Large Tubes, 22.4 x 11 x 11 cm, each. HS29040A	\$43.05
Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.	



Tube Racks, Mega Rack, polypropylene, Heathrow

10-13 mm Tubes	
216-Place, 33.7 x 26.7 x 9.1 cm(H)	
Blue, each. HS24322A	\$38.00
White, each. HS24322E	\$38.00
432-Place, 67.4 x 26.7 x 9.1 cm(H)	
Blue, each. HS120113	\$70.00
White, each. HS120112	\$70.00



Tube Racks, Microtube, polypropylene, Heathrow

80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm	
Holds 1.5-2.0 mL tubes in conical wells.	
Blue, Pk/5. HS29025B	\$26.10
Pink, Pk/5. HS29025D	\$26.10
Mixed, Pk/5. HS29025G	\$26.10
96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm	
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.	
Blue, Pk/5. HS2345C	\$51.15
Green, Pk/5. HS2345D	\$51.15
Pink, Pk/5. HS2345E	\$51.15