





Microscope Filters, glass

- 29.8 mm Diameter
 - Blue, clear. **MA289/05** \$35.00
 - Blue Frosted. **MA289** \$35.00
 - Diffusion (frosted). **MA288**..... \$35.00
- 32 mm Diameter
 - Blue, clear. **923-160** \$15.00
 - Blue Frosted. **3236** \$7.50
 - Diffusion (frosted). **3235-W**..... \$7.50
 - Green, clear. **CM-0321-G**..... \$7.00
- 40 mm Diameter
 - Blue Frosted. **MA563/05**..... \$50.00
 - Diffusion (frosted). **MA566** \$50.00
- 45 mm Diameter
 - Blue, clear. **14671**..... \$10.40
 - Color Balancing, NCB-11, Nikon. **MBN11710** \$77.00
 - Diffusion (frosted). **14673** \$35.00
 - Green Interference, 546/60, Nikon. **MBN11200** \$52.00
 - Neutral Density, Nikon
 - 50% Transmittance, ND2A. **MBN21802**..... \$42.00
 - 12.5% Transmittance, ND8A. **MBN21808**..... \$42.00



- Microscope Stage Laminate, USA. 16251** \$8.00
- 6.25 x 7.75" sheet provided with a #10 scalpel.
- White Teflon® (PTFE) on self-adhesive vinyl is 0.010" thick.



- Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420**..... \$0.50



Mark slides, dishes, flasks, and chambers with a 1.8 mm ring for later referral. Works on upright and inverted microscopes.

- Object Marker, RMS thread (20.32 mm Dia.)**
- Refill Ink, black, 10 mL, Shachihata. **22112** \$7.00
- Adapters
 - RMS/45 to CFI60, 15 mm long. **14111** \$43.00
 - RMS ID to 25 mm OD Thread, brass. **14115** \$52.00


Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.
Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.
Microscopes with 160 mm tube length objectives do not require an adapter.



Objectives, 160 mm tube length, DIN, Jenco, Meiji, National

- 4X, N.A. 0.10. **704-160** \$32.00
- 10X, N.A. 0.25. **710-160** \$33.00
- 40X, N.A. 0.65. **740-160** \$42.00
- 100X Oil, N.A. 1.25. **799-160**..... \$69.00
- Super High Contrast Achromats, National
 - 4X, N.A. 0.10. **704-160ASC** \$31.00
 - 10X, N.A. 0.25. **710-160ASC** \$39.00
 - 40X, N.A. 0.65. **740-160ASC** \$72.00
 - 100X Oil, N.A. 1.25. **799-160ASC** \$107.00
- Semi-Plan Achromats, National (unless noted)
 - 2X, N.A. 0.05, Accu-Scope. **3193SP-00**..... \$70.00
 - 4X, N.A. 0.10. **704-160SP** \$35.00
 - 10X, N.A. 0.25. **710-160SP** \$36.00
 - 40X, N.A. 0.65. **740-160SP** \$46.00
 - 50X Oil, N.A. 0.95, Meiji. **MA694** \$295.00
 - 60X, Strain Free, N.A. 0.85. **MA645** \$315.00
 - 100X Oil, N.A. 1.25. **799-160SP** \$75.00
- Plan Achromats, National (unless noted)
 - 4X, N.A. 0.10. **704-160P**..... \$36.00
 - 10X, N.A. 0.25. **710-160P**..... \$60.00
 - 20X, N.A. 0.40. **720-160P**..... \$135.00
 - 40X, N.A. 0.65. **740-160P**..... \$113.00
 - 50X Oil, N.A. 0.95, Accu-Scope. **3190PL-00**..... \$315.00
 - 60X, N.A. 0.85. **760-160P**..... \$214.00
 - 100X Oil, N.A. 1.25. **799-160P** \$181.00

Protect spare microscope objectives.



- Objective Case, DIN, polystyrene. 15601** \$0.50



Extender tubes do not contain lenses.

- Objective Extender Tubes, RMS thread, brass, USA**
- 10 mm Long. **14105** \$44.20
- 18 mm Long. **14106** \$42.75