



Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene

10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95

18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**..... \$49.95 ↓



Tools not included.

Dissecting Cases

6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30

12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40

13-Loop with zipper, brown leatherette. **36900**..... \$7.80

Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

Inquire about customizing a kit for your surgical procedure.

Kit components are discounted 10%.

Basic kit. 38501..... \$5.30

Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551..... \$17.85

Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90

25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic

Black Enameled, Pk/100..... \$7.00

#000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**

#00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**

#0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**

Stainless Steel, Pk/100..... \$8.75

#000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**

#00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**

#0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**..... \$43.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1 1/2". **35-406**..... \$26.60

Large, 10 x 4 x 2". **35-408**..... \$35.00

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461**..... \$0.95

White. **4460**..... \$0.55

Stainless Steel, 12 mm wide, **USA. 605-065**..... \$5.50




Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450..... \$419.00

Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$25.00

Continued



- Tube Rotator**, programmable, Benchmark (2-year warranty)
- Roto-Bot Mini. **R4040**\$647.00
- Roto-Bot. **R4045**.....\$972.00
- Speed adjusts from 1-80 rpm
 - Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
 - Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



- Vortex Mixers**, Benchmark (2-year warranty)
- MiniMixer, fixed 3000 rpm, blue cup. **VM-03UA**\$143.00
- Vornado™, mini vortexer, fixed 2800 rpm. **BV101**__*.....\$224.00
- *Specify **Gray**, **Blue** or **Purple** cup head
- BenchMixer™, single cup head, 200-3200 rpm. **BV1000**\$290.00
- Flat Head Attachment, 3". **BV1000-FLAT**\$42.60

Rulers, Meter Stick & Tape Measure



- Rulers**, 15 cm/6" long, scaled in both metric and English units
- Vinyl, flexible
- Transparent. **4461**\$0.95
- White. **4460**\$0.55
- Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065**\$5.70

- Meter Stick**, lacquered wood, plain ends. **605-050**\$4.90
- One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

- Tape Measure**, 152.4 cm/60", reinforced fiberglass. **11036**\$3.70
- One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



- Sharps Containers**, red polypropylene
- 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269**\$9.90
- Dimensions: 12" x 4.75" x 11"(H)
- 12 Gallon, hinged vertical entry lid, each. **440321391**\$30.75
- Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers



- Specimen Containers**, polypropylene (PP)
- Screw-Cap
- 90 mL, sterile, Cs/400. **840175**.....\$110.00
- 480 mL (16 oz), Cs/24. **59014510**\$71.95
- 83 mm (3 1/4") opening
- 960 mL (32 oz), Cs/12. **59014511**\$59.05
- 111 mm (4 3/8") opening
- Snap-Cap, white HDPE bottom/white LDPE lid
- 480 mL (16 oz), Cs/100. **271016W**.....\$82.45
- 9.9 cm opening x 11.1 cm high
- 1000 mL (34 oz), Cs/100. **271032W**.....\$124.00
- 11.7 cm opening x 13.7 cm high

Spot Plates



- Spot Plates**, glazed white porcelain, 21 mm diameter wells
- 6-Well, each. **LPC-161**\$2.35
- 12-Well, 115 x 95 x 7 mm, each. **3104187701**\$4.30

Rocks & Minerals



Specimens are in a package of ten pieces, each 3-4 cm in size.
All specimens are shipped in a labeled, resealable bag.

Minimum order of \$100 required for rocks and minerals (any combination). Please call about availability of rocks not listed.




Minerals

Albite with twinning, Pk/10. 310560501	\$8.90
Augite, cleavable, Pk/10. 310562802	\$12.85
Calcite, cleavages, white-brown, Pk/10. 310564701	\$7.90
Fluorite, cleavable, Pk/10. 310569402	\$11.90
Graphite, pure, Pk/10. 310570601	\$16.85
Magnetite, black, coarse crystalline, Pk/10. 310573301	\$5.95
Muscovite, cleavable, Pk/10. 310574601	\$11.40
Olivine, massive, medium-grained, Pk/10. 3105751	\$6.55
Plagioclase Feldspar, Pk/10. 3105725	\$7.90
Pyrite, pure, Pk/10. 310576902	\$17.80
Quartz, Rose, massive, Pk/10. 310578012	\$5.95
Quartz, Silicified Wood, Pk/10. 310578014	\$6.95
Rhodonite, pink, massive, Pk/10. 310577801	\$15.85
Talc, schist, Pk/10. 3105801	\$6.95



Igneous Rocks

Anorthosite, gray, coarse-grained, Pk/10. 3105006	\$5.95
Gabbro, brown-black, coarse-grained, Pk/10. 3105024	\$6.55
Granite, red, coarse-grained, Pk/10. 310502701	\$4.95
Pegmatite, very coarse-grained, Pk/10. 3105044	\$4.95
Pumice, white-gray, Pk/10. 310504801	\$7.90
Volcanic Ash, light gray, Pk/10. 3105069	\$11.85



Sedimentary Rocks

Arkose, red, coarse-grained, Pk/10. 3105103	\$4.95
Breccia, crystalline in rock, Pk/10. 310510803	\$5.45
Coal, bituminous, Pk/10. 3105112	\$5.95
Conglomerate, quartz pebbles, Pk/10. 310511601	\$4.95
Limestone, oolitic, granular, Pk/10. 310512002	\$5.95
Shale, gray, calcereous, Pk/10. 310514501	\$5.95




Metamorphic Rocks

Gneiss, coarse-grained, banded, Pk/10. 310521601	\$5.95
Marble, white, coarse-grained, Pk/10. 310523101	\$4.95
Phyllite, gray-green, fine texture, Pk/10. 3105235	\$4.95
Quartzite, blue-gray, fine-grained, Pk/10. 310523801	\$5.95
Slate, black-gray, Pk/10. 3105245	\$7.90

Novelty Items

Micro Bible. 89585	\$15.00
View with a 10X objective all 1,245 pages of the King James version of the Holy Bible on this approximately 1" square microfilm.	



Ruler, 1 Decimeter, stainless steel, Germany. 605-060	\$7.40
Ruler is 1 dm long, 1 cm wide, and 1 mm thick. Its volume is 1 cm ³ and its weight in grams is the specific gravity. One edge of one side is graduated in centimeters with 1 mm divisions.	

Reduce waste & save money with this heavy-duty wringer by completely extracting contents from even the most rigid tube. Features 3" aluminum rollers with steel handle.



Tube-Wringer®, USA. 295401	\$29.95
-----------------------------------	---------