



Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. 1596..... \$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene

10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**..... \$6.95

18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**..... \$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. DEM-345.00..... \$49.95 ↓



Tools not included.

Dissecting Cases

6-Loop with snap fastener, black leatherette. **1257370**..... \$1.30

12-Loop with snap fastener, black leatherette. **1257371**..... \$1.40

13-Loop with zipper, brown leatherette. **36900**..... \$7.80

Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

Inquire about customizing a kit for your surgical procedure.

Kit components are discounted 10%.

Basic kit. 38501..... \$5.30

Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551..... \$17.85

Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)

20 7/8 x 12 13/16 x 2 1/2". **AM20029**..... \$21.90

25 3/4 x 17 3/4 x 2 1/4". **AM68357**..... \$64.90



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic

Black Enameled, Pk/100..... \$7.00

#000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**

#00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**

#0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**

Stainless Steel, Pk/100..... \$8.75

#000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**

#00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**

#0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. 150-1005..... \$43.00




Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel

Small, 7 x 3 x 1 1/2". **35-406**..... \$26.60

Large, 10 x 4 x 2". **35-408**..... \$35.00

Stainless Steel



Rulers, 15 cm/6" long, scaled in both metric and English units

Vinyl, flexible

Transparent. **4461**..... \$0.95

White. **4460**..... \$0.55

Stainless Steel, 12 mm wide, **USA. 605-065**..... \$5.50



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, USA. 5-1450..... \$419.00

Glass Beads (replacement), 1.5 mm, Bg/10 oz. **5-1455**..... \$25.00

Continued



Fuses

¼ x 1¼"

Fast-Acting, AGC

0.1A/250V. 312-100	\$2.00
0.2A/250V. 312-200	\$2.00
1A/250V. 312-001	\$2.00
2A/250V. 312-002	\$2.00
3A/250V. 312-003	\$2.00
4A/250V. 312-004	\$2.00

Time-Delay (Slow-Blow), MDL

0.7A/250V. 313-700	\$2.00
1A/250V. 313-001	\$2.00
2.5A/250V	\$2.00

5 x 20 mm

Fast-Acting, GMA


0.5A/250V. 235-500	\$2.00
1A/250V. 235-001	\$2.00
3A/250V. 235-003	\$2.00
5A/125V. 235-005	\$2.00

Time-Delay (Slow-Blow), Glass tube


0.5A/250V. 239-500	\$2.00
1A/250V. 239-001	\$2.00
2A/250V. 239-002	\$2.00
2.5A/250V	\$2.00
3A/250V. 239-003	\$2.00
5A/250V. 239-005	\$2.00

Time-Delay (Slow-Blow), Ceramic

1A/250V. 215-001	\$2.00
-------------------------------	--------



Cone Point



Cup Point

Set Screws, steel

Cone Point, DIN 914

M4 x 0.7 x 5 mm length, 2 mm hex key size. 77702	\$0.50
M4 x 0.7 x 10 mm length, 2 mm hex key size. 77704	\$0.50
M4 x 0.7 x 12 mm length, 2.5 mm hex key size. 77706	\$0.50
M5 x 0.7 x 10 mm length, 2.5 mm hex key size. 77708	\$0.50

Cup Point, DIN 916

M3 x 0.5 x 3 mm length, 1.5 mm hex key size. 77722	\$0.50
M3 x 0.5 x 5 mm length, 1.5 mm hex key size. 77724	\$0.50

Tools




Cutter, diagonal, 5". **PLR-766.00** \$10.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1½ x 3" opening. **DEM-345.00** \$49.95 ↓

Gauges



Vernier Caliper measures in millimeters (0-150 mm with 0.01 mm precision) and inches (0-6"). Provided with a plastic case.

Vernier Caliper, digital, 3-button. **600VCD001** \$33.35

Battery (replacement), Silver Oxide. **357** \$4.00



Reproducibility enhanced by rapid stop. Provided with padded, plastic case.

Micrometer, 0-25 mm in 0.01 mm div, Germany. **601A** \$40.70