

Tubes will withstand centrifugal forces over 30,000 x g and temperatures from -80° to 121°C. Flat, textured caps are pierceable. Free of RNase/DNase.

**Tubes, Microcentrifuge, Premium, PP, G-Tube™, Bio Plas, USA**

Amber color protects light sensitive specimens.

- 0.6 mL, Pk/500 .....\$26.25  
Amber. **4048** Lavender. **4047** Natural. **4040**  
Mixed (blue, green, orange, pink, & yellow). **4049** .....\$31.65
- 1.5 mL, 0.5 mL graduations, frosted area on side, Pk/500.....\$26.25  
Amber. **4038** Blue. **4033** Lavender. **4037** Natural. **4030**  
Mixed (blue, green, orange, pink, & yellow). **4039** .....\$31.65
- 2.0 mL, 0.5 mL graduations, frosted area on side  
Natural, Pk/500. **4050**.....\$28.20

**Tubes, Microcentrifuge, Siliconized, G-Tube™, Bio Plas, USA**

- 1.5 mL, natural color, Pk/500. **4165SL**.....\$52.55  
*Siliconized products should NOT be autoclaved*



**Tubes, PCR, 0.2 mL thin-wall, polypropylene, Simport**

- Individual, flat caps, natural, Pk/1000. **T325-2VN** .....\$49.00



**Vials, Cryogenic, round bottom, internal thread**


- O-Ring Seal, 4 mL, Cs/1000. **T301-4** .....\$297.00
- Washer Seal, 4 mL, Cs/1000. **T311-4** .....\$297.00



Vinyl grip, **USA**

- Vial Decapping Pliers, for 13 mm Crimp Seals. D-13** .....\$72.00

**Tube Storage**



**Baskets, test tube, aluminum, Heathrow**

- Small, 15.2 x 15.2 x 15.2 cm(H). **HS20341A** .....\$27.00
- Medium, 25.6 x 15.2 x 15.2 cm(H). **HS20341B** .....\$31.00

**Cryogenic Storage Boxes**



SB2CC-64 Cryogenic  
Tube Storage Box

**Tube Storage Boxes, Cardboard, vented**

- 2" Box Only, 5¼ x 5¼ x 2", each. **SB2CC** .....\$1.80
- 2" Box w/dividers, white
  - 100-Cell (10 x 10), Pk/12. **SB2CC-100** .....\$28.45
  - ★ 81-Cell (9 x 9), Pk/12. **SB2CC-81** .....\$28.20
  - 64-Cell (8 x 8), Pk/12. **SB2CC-64** .....\$28.20
- 2" Mini box w/dividers, 3 x 3 x 2"
  - 25-Cell (5 x 5), Pk/12. **MB2C-25** .....\$30.65  
Holds 13 mm vials and 4.0-5.0 mL cryovials.
- 3" Box Only, 5¼ x 5¼ x 3", each. **SB3C** .....\$1.80
- 3" Box w/dividers
  - 100-Cell (10 x 10), Pk/12. **SB3C-100** .....\$28.45
  - 81-Cell (9 x 9), Pk/12. **SB3C-81** .....\$28.20
  - 64-Cell (8 x 8), Pk/12. **SB3C-64** .....\$28.20
- Dividers only, Fit 2 & 3" Boxes
  - 100-Cell (10 x 10), Pk/12. **CD-100** .....\$11.85  
Holds 1.5-2.0 mL tubes and 1.2-2.0 mL cryovials.
  - 81-Cell (9 x 9), Pk/12. **CD-81** .....\$11.55  
Holds 13 mm vials and 4.0-5.0 mL cryovials.
  - 64-Cell (8 x 8), Pk/12. **CD-64** .....\$11.55  
Holds 14 mm vials.
  - 49-Cell (7 x 7), Pk/12. **CD-49** .....\$11.85  
Holds 14 mm vials.
  - 25-Cell (5 x 5), 25 mm cell, Pk/12. **CD-25** .....\$11.55
  - 16-Cell (4 x 4), 31 mm cell, Pk/12. **CD-16** .....\$11.55
- 5" Box Only, 5¾ x 5¾ x 5", each. **SB4C** .....\$3.30
- 5" Box w/dividers
  - 64-Cell (8 x 8), Pk/12. **SB4C-64** .....\$43.50
  - 36-Cell (6 x 6), 15 mL tubes, Pk/12. **SB4C-36** .....\$43.50
  - 16-Cell (4 x 4), 50 mL tubes, Pk/12. **SB4C-16** .....\$43.50
- Dividers only, Fit 5" Box
  - 64-Cell (8 x 8), Pk/12. **CD4-64** .....\$11.55
  - 36-Cell (6 x 6), Pk/12. **CD4-36** .....\$11.55
  - 16-Cell (4 x 4), Pk/12. **CD4-16** .....\$11.55



Holds 1.0-2.0 mL tubes and cryovials.  
Features a locking hinged lid, and base with molded alphanumeric coordinates.  
Dimensions: 12.8 x 12.8 x 5.4 cm.

**Tube Storage Boxes, Polypropylene**

81-Well (9x9), 2 blue/2 yellow boxes, Pk/4. **60066501**.....\$30.25

Holds 1.5-2.0 mL tubes in a two-piece rack.



**Tube Storage Boxes, Polystyrene Foam (styrofoam), USA**

50-Well, Pk/10. **HS4283P**.....\$47.80  
100-Well, Pk/10. **HS4282**.....\$59.15

**Tube Racks**




**Tube Racks, Grip**, polypropylene, No-Wire™, Bel-Art

15-16 mm Tubes, 60-place, each. **187490001**.....\$32.65  
18-20 mm Tubes, 40-place, each. **187490002**.....\$32.65



4-Way Rack

HS29022G



4-Way Rotatable Rack


HS29040A

**Tube Racks, Interlocking**, polypropylene, Heathrow

4-Way, mixed colors, Pk/5. **HS29022G**.....\$43.90  
Holds (4) 50 mL tubes, (12) 15 mL tubes, (32) 1.5 mL tubes, or (32) 0.5 mL tubes.  
4-Way Rotatable, Rota-Rack™, four connected modules  
Large Tubes, 22.4 x 11 x 11 cm, each. **HS29040A** .....\$45.20  
Each module holds (2) 50 mL tubes, (8) 15 mL tubes, (6) 20 mm tubes, or (10) 10-13 mm tubes. Mixed colors.




HS24322A/E



HS120112

**Tube Racks, Mega Rack**, polypropylene, fully-assembled

10-13 mm Tubes  
216-Place, 33.7 x 26.7 x 9.1 cm(H)  
Blue, each. **HS24322A**.....\$47.30  
White, each. **HS24322E**.....\$47.30  
432-Place, 67.4 x 26.7 x 9.1 cm(H)  
Blue, each. **HS120113**.....\$84.00  
White, each. **HS120112**.....\$84.00



80-Well



96-Well Reversible

**Tube Racks, Microtube**, polypropylene, Heathrow

80-Well (16 x 5), molded-in coordinates, 22.5 x 6.7 x 2.8 cm  
Holds 1.5-2.0 mL tubes in conical wells.  
Blue, Pk/5. **HS29025B**.....\$26.90  
Pink, Pk/5. **HS29025D**.....\$26.90  
Mixed, Pk/5. **HS29025G**.....\$26.90  
96-Well reversible (12 x 8), snap-on cover, 24.6 x 12.1 x 5 cm  
Holds 1.5-2.0 mL tubes on one side or 0.5 mL on reverse.  
Blue, Pk/5. **HS2345C**.....\$53.75  
Green, Pk/5. **HS2345D**.....\$53.75  
Pink, Pk/5. **HS2345E**.....\$53.75

Available in blue, green, orange, purple & yellow.



**Tube Rack, PCR**, polypropylene, mixed colors, Heathrow

Reversible, Pk/5. **HS2344A**.....\$58.15  
Holds (168) 0.2, (40) 0.5 and (12) 1.5 mL tubes.



**Tube Racks, Peg**, 14-17 mm Tubes, 50-place, polypropylene

White, Pk/2. **HS24311A**.....\$24.80

Continued

Continued