



Tilt your test tube rack 25°, 35°, or 45°. Create agar slant tubes, or adjust the angle of tube racks when pipetting to reduce hand strain.

Tube Rack Slanter, Bel-Art. **188900100**\$32.05
Platform size: 19.6 cm(W) x 12.7 cm(D)



Tube Racks, Teardrop, OneRack™, POM, Heathrow
Half-Size
16 mm Tubes, 36-place, 12.7 x 12.7 x 7 cm(H)
Magenta, Pk/8. **HS27502D**\$109.00
Full-Size
13 mm Tubes, 72-place, 20.2 x 10.4 x 5.9 cm(H)
Blue, Pk/8. **HS27511B**\$143.00
Magenta, Pk/8. **HS27511D**\$143.00
16 mm Tubes, 72-place, 25 x 12.7 x 7 cm(H)
Blue, Pk/8. **HS27512B**\$143.00
Magenta, Pk/8. **HS27512D**\$143.00



Polypropylene cube and lid has two snap-in, gridded inserts to hold either (16) 15 mL tubes or (9) 50 mL tubes.

Tube Racks, Transformer Cube, Pk/5. **HS23271A**\$60.40



Tube Racks, Wire, epoxy-coated, blue
10-13 mm Tubes, 108-place (9 x 12), each. **345120759**\$32.00
13-16 mm Tubes
108-Place (9 x 12), each. **345120763**\$32.00
15-16 mm Tubes
100-Place (5 x 20), **USA**, ea. **187860780**\$69.65
Dimensions: 38.6 x 10.8 x 6.4 cm(H)
180-Place (9 x 20), **USA**, ea. **187860770**\$68.45
Dimensions: 38.6 x 18.3 x 6.4 cm(H)

17515ML-25



Tube Racks, 15 mL

25-Place (5 x 5), stainless steel, **USA**, ea. **17515ML-25**\$14.80
80-Place (5 x 16), stainless steel, **USA**, ea. **17515ML-80**\$33.40
Dimensions: 40.7 x 12.7 x 10.2 cm(H)



17550ML-9

Tube Racks, 50 mL

9-Place (3 x 3), stainless steel, **USA**, ea. **17550ML-9**\$14.80
16-Place (4 x 4), epoxy-coated wire, ea. **187940001**\$65.55

Weighing Dishes



Hexagonal



Standard



Pour

Weigh Boats, polystyrene, Heathrow

Standard
★ Small, 4.6 x 4.6 x 0.8 cm(D), Bx/500. **HS1420A**\$23.45
★ Medium, 8.5 x 8.5 x 2.4 cm(D), Bx/500. **HS1420B**\$38.15
★ Large, 14 x 14 x 2.2 cm(D), Bx/500. **HS1420C**\$62.15
Hexagonal, anti-static
★ Small, 3.8 cm top dia., Bx/500. **HS14251A**\$17.25
★ Medium, 6.4 cm top dia., Bx/500. **HS14251B**\$29.75
★ Large, 12 x 14 cm top OD, Bx/500. **HS14251C**\$77.85
(25433-104)
★ X-Large, 14 x 16 cm top OD, Bx/500. **HS14251D**\$98.75
(25433-106)
Pour (Canoes)
★ Small, 5.8 x 4.3 x 1.3 cm(D), Bx/500. **HS1419A**\$38.90
★ Medium, 13.2 x 8.3 x 2.6 cm(D), Bx/500. **HS1419B**\$76.50
★ Large, 18.3 x 10.8 x 2.6 cm(D), Bx/500. **HS1419C**\$111.00