

Student Microscope



National

131-RLED Student Microscope, Monocular, DIN..... \$259.00

- Widefield 10X eyepiece w/ pointer accepts optional reticle
- Monocular, 45° inclined viewing head rotates 360°
- DIN achromat 4X, 10X, and 40XR objectives
- Stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- LED illuminator, rechargeable batteries (equivalent to a 20 watt tungsten bulb)
- 36 cm (14") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K

Stereo Microscopes



440-440PLL

National

440-440PLL Stereoscope, Binocular, 7.5-45X Zoom \$879.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-45X total magnification; working distance of 110 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (LED with rheostat control and directional beam), bottom only (LED), or both top and bottom
- Post stand with locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Limited lifetime mechanical and optical/3-year electrical

440T-440PLL Stereoscope, Trinocular, 7.5-45X Zoom \$999.00

- Trinocular version of 440-440PLL-10

M30TZ-SM99CL



Swift

M30TZ-SM99CL Stereoscope, Trinocular, 7.5-45X Zoom.... \$914.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Trinocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-45X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators with rheostat control offer cool light from top, bottom, or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warrant

Fiber Optic Illuminator

Please call regarding availability of items not listed



Fiber optic light source with dual gooseneck and two focusing lenses is shown. Items are sold separately.

Fiber Optics, Fiber-Lite™, Dolan-Jenner

Light Source, 150W EKE lamp, rheostat. **MI-150** \$660.00

(FL150/115) Warranty: 3 years

Gooseneck, Dual, chrome plated, 23" long. **EEG2823M** \$420.00

(FL150/05) Warranty: 1 year

Focusing Lens, each. **LH-759** \$75.00

(FL150/50) Warranty: 1 year