

# SLIDE & CASSETTE PROCESSING

## Microscope Slides

### Slides, general, 75 x 25 x 1 mm

Plain, ground edges

Economy, safety corners, C&A

Cs/20 Bx, bx/72. **9101-E**.....\$90.00

Premium, 90° corners, "Globe"

Cs/20 Bx, bx/72. **1301**.....\$134.00

Frosted, beveled edges

Economy, sandblasted, safety corners, C&A

Cs/20 Bx, bx/72. **8205**.....\$146.00

Premium, etched, clipped corners, "Globe"

Cs/20 Bx, bx/72. **1334**.....\$177.00

White End, clipped corners

Economy, C&A

Cs/20 Bx, bx/72. **9108W-E**.....\$148.00

Superfrost®, "SPECIMEN". **4951-001T**.....\$294.00

(12-550-143, 48311-600)

Superfrost® Plus



Concavity, 2-Well



### Slides, specialty

Charged (silanated), 75 x 25 x 1 mm

Recommended for adhering tissue sections  
and cell preparations.

Diamond™, white end, Globe

Cs/20 Bx, bx/72. **1358W**.....\$391.00

Superfrost® Plus, white end, Erie (Epredia)

Cs/20 Bx, bx/72. **4951PLUS-001**.....\$660.00

Preferred slide for Ventana instrumentation.

(12-550-15, 48311-703)

Concavity

75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy

1-Well, Bx/12. **60CS3X11**.....\$7.40

2-Well, Bx/12. **60CS3X12**.....\$10.20

76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie (Epredia)

1-Well, Bx/12. **1526**.....\$97.35

Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, **USA**

**QMS3X1X1**, each.....\$20.00

Fused Silica, Corning 7980, 3" x 1" x 1 mm, **USA**

**FSMS3X1X1**, each.....\$29.00

## Coverslips (Cover Glasses)

Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)

Circles, Germany\*



25 mm, ~222 pcs/oz, Cs/10 oz. **0-25**.....\$214.00

#1 Thickness (0.13-0.17 mm)

Squares & Rectangles, Erie (Epredia) 1 oz cs/10 oz



22 x 22 mm, ~162 pcs/oz. **1-2222**.....\$66.15

22 x 30 mm, ~119 pcs/oz. **1-2230**.....\$59.70

22 x 50 mm, ~71 pcs/oz. **1-2250**.....\$59.70

24 x 50 mm, ~65 pcs/oz. **1-2450**.....\$59.70

24 x 60 mm, ~54 pcs/oz. **1-2460**.....\$59.70

#1 Thickness (0.13-0.17 mm)

Circles, Germany\*



12 mm, ~663 pcs/oz. **1-12**.....cs/5 oz cs/10 oz

18 mm, ~296 pcs/oz. **1-18**.....\$249.00 —

18 mm, ~296 pcs/oz. **1-18**.....\$150.00 —



#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*

Squares & Rectangles, Erie (Epredia) 1 oz cs/10 oz

18 x 18 mm, ~209 pcs/oz. **1.5-1818**.....\$63.90

22 x 22 mm, ~140 pcs/oz. **1.5-2222**.....\$7.80 \$65.40

22 x 50 mm, ~62 pcs/oz. **1.5-2250**.....\$59.70

24 x 40 mm, ~71 pcs/oz. **1.5-2440**.....\$59.70

24 x 50 mm, ~56 pcs/oz. **1.5-2450**.....\$59.70

24 x 60 mm, ~47 pcs/oz. **1.5-2460**.....\$59.70

#1.5 Thickness (0.16-0.19 mm). *Recommended for most users*


Circles, Germany\*

12 mm, ~552 pcs/oz. **1.5-12**.....cs/5 oz cs/10 oz

18 mm, ~245 pcs/oz. **1.5-18**.....\$198.00 —

22 mm, ~164 pcs/oz. **1.5-22**.....\$250.00 —

25 mm, ~127 pcs/oz. **1.5-25**.....\$252.00 —



\*German glass is preferred for cell growth/adhesion

### Coverslips, specialty

Fused Quartz, **USA**

0.15-0.25 mm Thick, NSG Type OZ®

1" Square, each. **QCS1X1**.....\$19.00

1" Circle, each. **QCS1X0**.....\$19.00

0.16 ± 0.01 mm Thick, A1

25 mm Square, each. **R525025**.....\$32.40

25 mm Circle, each. **R525000**.....\$23.40

Fused Silica, **USA**

0.15-0.25 mm Thick, Corning 7980

1" Square, each. **FSCS1X1**.....\$21.50

1" Circle, each. **FSCS1X0**.....\$21.50

0.16 ± 0.01 mm Thick, S1-UV

25 mm Square, each. **R425025**.....\$32.00

25 mm Circle, each. **R425000**.....\$28.80

Hemocytometer Coverslips

20 x 26 x 0.4 mm, Erie, Pk/2 **2026H**.....\$11.05

*See additional Counting Chamber Coverslips in  
MICROSCOPE Accessories & Parts section*

## Accessories

### Slide & Coverslip Handling

#### Forceps

##### Spade Tips



- Angled (Kuehne), 11 cm (4 1/4")
  - Stainless Steel, economy. **33-1041** ..... \$4.15
  - Teflon® Coated, Bel-Art. **379420000** ..... \$26.85
- Straight, stainless steel, 10.5 cm, economy. **57.926** ..... \$8.45

##### Square Point Tips, straight



- Stainless Steel, 12 cm (4 3/4"), economy. **1244** ..... \$2.25
- Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** ..... \$29.70



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

- Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** ..... \$1.95

### Cytology Funnels

Single,  
Metal Clips  
Required



Double,  
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

#### Cytology Funnels, Simport CytoSep™

##### Shandon (Thermo) Cytospin® Cytocentrifuges

##### Single. *Metal slide clips required*

- White Filter & Cap, tray/50, Cs/500. **M964-10FW** ..... \$426.00  
(59 91 040, 110 25 48)
- White Filter & Cap, Pk/1, Cs/500. **M964-10FW1** ..... \$563.00
- Tan Filter & Cap, tray/50, Cs/500. **M964-10FT** ..... \$426.00
- Cytology Funnel Clips, metal, Pk/6. **M964B** ..... \$176.00  
*Usable with the above three items*

##### Single. *All plastic - metal slide clips not required*

- White Filter & Cap, Cs/480. **M965-10FW** ..... \$628.00  
(A787 100003)
- Tan Filter & Cap, Cs/480. **M965-10FT** ..... \$628.00

##### Double. *All plastic - metal slide clips not required*

- White Filter & Cap, Cs/480. **M965-20FW** ..... \$835.00  
(A787 100005)

### Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.

**USA**

#### Immersion Oil

##### Type A, low viscosity (150 cSt). *Short W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18101** ..... \$10.00
- 20 mL in HDPE dispensing bottle. **18102** ..... \$13.90
- 480 mL in amber bottle. **18104** ..... \$94.50

##### Type B, high viscosity (1250 cSt). *Long W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18105** ..... \$10.00
- 20 mL in HDPE dispensing bottle. **18106** ..... \$13.90
- 480 mL in amber bottle. **18108** ..... \$94.50

##### Type NVH, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** ..... \$67.00


##### Type HF, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** ..... \$144.00 ↓

##### Type FF, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** ..... \$140.00


A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



- Bottle, Dispensing, HDPE, 20 mL, each. **18145** ..... \$8.10

**Slide Storage**

**PLASTIC SLIDE BOXES**



- Slide Boxes**, hold 75 x 25 mm (3 x 1") slides
- 25-Slide, squeeze-release cover, cork bottom, 9.5 x 8.5 x 3 cm  
White, ABS plastic, China. **97-0025** ..... \$2.25
  - 50-Slide, hinged cover, cork bottom, 20 x 9.5 x 3.5 cm  
Blue, ABS plastic, China. **97-0056** ..... \$6.30  
Red, ABS plastic, China. **97-0058** ..... \$6.30
  - 100-Slide, steel-hinged cover, cork bottom  
Blue, ABS plastic, China. **97-0126** ..... \$8.10  
Red, ABS plastic, China. **97-0128** ..... \$8.10  
White, ABS plastic, China. **97-0125** ..... \$8.10



Store 4500 75 x 25 mm slides in 6 divided drawers.  
Each cabinet comes with a sample of 6 compressor follow blocks and 6 drawer labels. Safely stack up to 12 cabinet units on a single base.  
18¾"(D) x 15¾"(W) x 5"(H)

- Slide Cabinets**, steel
- 6-Drawer, stackable
  - Green. **513500G** ..... \$219.00
  - Tan. **513500T** ..... \$219.00
  - White. **513500W** ..... \$219.00
  - Memorandum Markers, paperboard, Pk/100. **C108** ... \$12.75
  - Compressor Follow Block, rubber, each. **C109** ..... \$1.90
  - Base, raises cabinets 4¾" off floor
  - Green. **513503G** ..... \$141.00
  - Tan. **513503T** ..... \$141.00
  - White. **513503W** ..... \$141.00





Cardboard Folder



Plastic Folder

- Slide Folders**, 20-slide, hold 75 x 25 mm (3 x 1") slides
- Cardboard, black tape, economy. **HS9902** ..... \$6.00  
(12-587-10)
  - Plastic, gray tray w/2 clear doors. **M750-20GY** ..... \$11.10  
Transparent covers hold visible slides securely in place, but can fold flat against back of tray when opened.  
Slides are slightly elevated for easy retrieval.




- Slide Holders**, perfect for slide drying or staining
- 52-Slide (approx.), 12 x 4", HDPE, **USA. L500-C** ..... \$30.95
  - 91-Slide (approx.), 12 x 7", HDPE, **USA. L500-D** ..... \$41.35




5-Slide mailer can also be used as a staining jar.

- Slide Mailer**, 5-Slide, plastic, end-opening, Cs/200. **62706** ..... \$87.30




SL-3-264



SL-3-182

- Slide Sleeves**, 10 gauge vinyl
- 4 x 5½" Card Holder Plus 4 Slide Pockets. **SL-3-264** ..... \$0.45
  - 8½ x 11" Page Holds 18 Slides. **SL-3-182** ..... \$0.75



Aluminum



Plastic

- Slide Trays**, 20-slide
- Aluminum, Eberbach, **USA. E2870** ..... \$21.05
  - Plastic, white. **445810000** ..... \$19.05

## Slide & Cassette Pens, Marker & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.  
Pilot pens are ideal for cytologists and pathologists.



### Pens, Pilot®, permanent, extra fine point

Black, Bx/12. <b>44102</b> .....	\$17.45
Blue, Bx/12. <b>44103</b> .....	\$17.45
Green, Bx/12. <b>44105</b> .....	\$17.45
Red, Bx/12. <b>44104</b> .....	\$17.45



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

**Marker, Diamond PerfectMark™, black, Bx/10. 1000**..... \$28.90



### Scribes, retractable

Carbide-Tipped (sharpenable). <b>RSC</b> .....	\$24.30
Carbide Tip (replacement). <b>RCT-1</b> .....	\$8.45
Diamond-Tipped. <b>RSD</b> .....	\$28.35
Diamond Tip (replacement). <b>RDT-1</b> .....	\$14.35

## Slide Warmers




Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

### Slide Warmers

Small, 26 x 18 cm (10 x 7"). <b>XH-2002</b> .....	\$328.00
Large, 63 x 20 cm (25 x 8"). <b>XH-2001</b> .....	\$431.00

## Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

### Staining/Moisture Chambers

10-Slide, polystyrene, <b>USA</b>	<u>ea</u>	<u>pk/3</u>
Clear. <b>240-9000-010</b> .....	\$54.95	154.00
Black. <b>240-9020-Z10</b> .....	\$56.10	157.00



### 20-Slide, Simport

Clear Lid (PETG), black base (ABS), each. <b>M920-1</b> .....	\$186.00
Black Lid (ABS), black base (ABS), each. <b>M920-2</b> .....	\$186.00




### Staining Dishes

#### Glass (flint)


20-Slide, holds slides horizontally in metal rack	
Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)	
Set (dish, rack & lid). <b>SDS-20</b> .....	\$29.50
Dish & Lid. <b>SD-20</b> .....	\$18.00
Rack only, stainless steel. <b>SR-20</b> .....	\$14.80
30-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)	
Set (dish, rack & lid). <b>SDS-30</b> .....	\$40.30
Dish & Lid. <b>SD-30</b> .....	\$20.05
Rack only, stainless steel. <b>SR-30</b> .....	\$24.75
60-Slide, holds slides vertically in metal rack	
Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)	
Set (dish, rack & lid). <b>SDS-60</b> .....	\$56.40
Dish & Lid. <b>SD-60</b> .....	\$30.85
Rack only, stainless steel. <b>SR-60</b> .....	\$31.80

Continued



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides  
Set (dish, rack & lid). **0990-99** ..... \$50.90



**Staining Jars**, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, USA. <b>673</b> .....	\$12.50	142.00
Glass Cover, Fits 673 Jar. <b>673-LID</b> .....	\$4.70	—
	<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, <b>USA 674</b> .....	\$12.50	142.00
Phenolic Cap, Fits 674 Jar, black. <b>674-LID</b> .....	\$4.75	—
	<u>ea</u>	<u>cs/12</u>
PP w/Screw-On Cap. <b>60062101</b> .....	\$4.90	47.40
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
	<u>ea</u>	<u>cs/12</u>
Wide Top, Glass w/Square Lid, <b>USA 675</b> .....	\$16.95	189.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. <b>675-LID</b> .....	\$5.25	—



Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)  
12-Slide, EasyDip™, Simport

Jar only, Blue, Cs/6. <b>M900-12B</b> .....	\$64.95
Jar only, Green, Cs/6. <b>M900-12G</b> .....	\$64.95
Jar only, Pink, Cs/6. <b>M900-12P</b> .....	\$64.95
Jar only, White, Cs/6. <b>M900-12W</b> .....	\$64.95
Jar only, Yellow, Cs/6. <b>M900-12Y</b> .....	\$64.95
Rack/Lid only, gray, Cs/6. <b>M905-12DGY</b> .....	\$71.75




Anodized aluminum staining jar holder

(Jars sold separately)  
Staining Jar Holder, 6-Place. **M906** ..... \$79.45



24-Slide, light green, microwaveable, acetal  
Jar & Lid. **4501030** ..... \$11.55  
Rack only, gray. **4501042** ..... \$21.45  
Staining Jar Station, 12-Place. **4501011** ..... \$221.00  
Includes: metal frame, (12) 24-slide dishes w/ lids (4501030) & (1) rack (4501042).  
Dimensions: 29(L) x 5(W) x 5"(H)

### Tissue (Embedding) & Biopsy Cassettes



**Base Molds**



Plastic, disposable, Simport

15 x 15 x 5 mm, Cs/1000. <b>M475-2</b> .....	\$148.00
24 x 24 x 5 mm, Cs/1000. <b>M475-3</b> .....	\$148.00
★ 30 x 24 x 5 mm, Cs/1000. <b>M475-4</b> .....	\$148.00
37 x 24 x 5 mm, Cs/1000. <b>M475-5</b> .....	\$148.00
37 x 24 x 10 mm, Cs/500. <b>M475-10</b> .....	\$76.45
Designed for Simport Macrosette™.	

Metal, reusable

15 x 15 x 5 mm, Pk/12. <b>HB-15</b> .....	\$27.40
24 x 24 x 5 mm, Pk/12. <b>HB-24</b> .....	\$27.40
30 x 24 x 5 mm, Pk/12. <b>HB-30</b> .....	\$27.40
37 x 24 x 5 mm, Pk/12. <b>HB-37</b> .....	\$27.40

Solvent resistant,  
0.2 x 0.2 mm mesh.




**Biopsy Bags**

Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. <b>4502011</b> .....	\$235.00
Paper (Tea Bags), 2¼ x 2¼", Pk/500. <b>56100</b> .....	\$45.15

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



M476-2

**Biopsy Pads**


Rectangular

Black Foam, 30.2 x 25.4 mm, Bg/1000. <b>339014</b> .....	\$42.00
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. <b>M476-1</b> .....	\$213.00

Circular, Simport

Blue Foam, 25.4 x 2.7 mm, Cs/10,000. <b>M476-2</b> .....	\$260.00
--	----------

**Biopsy Cassettes**



Simport  
Histosette I




Simport  
Micromesh

- Biopsy Cassettes, Attached Lids**, acetyl polymer, Simport
- Histosette™ I, 1 mm square pores, Pk/1500 ..... \$147.00  
30° front labeling area.
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M491-6</b>  | Pink. <b>M491-3</b>  |
| Green. <b>M491-4</b> | White. <b>M491-2</b> |
- Unisette™, hinged lids, Pk/1500 ..... \$157.00  
35° front labeling area. *(Alternative to Sakura's Uni-Cassette®)*
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M506-6</b>  | Pink. <b>M506-3</b>  |
| Green. <b>M506-4</b> | White. <b>M506-2</b> |
- Micromesh™, 0.38 mm pores, Pk/1000 ..... \$241.00  
45° front labeling area.
- |                      |                      |                       |
|----------------------|----------------------|-----------------------|
| Blue. <b>M507-6</b>  | Pink. <b>M507-3</b>  | Yellow. <b>M507-5</b> |
| Green. <b>M507-4</b> | White. <b>M507-2</b> |                       |

**Tissue Cassettes**

Tissue  
Cassettes



- Tissue Cassettes, Detached Lids**, slotted pores
- Histosette™ II, Simport, Pk/1000 ..... \$100.00  
45° front labeling area.
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M485-6</b>  | Pink. <b>M485-3</b>  |
| Green. <b>M485-4</b> | White. <b>M485-2</b> |
- Uni-Cassettes™, **USA**, Pk/500 ..... \$74.45
- |                           |                     |
|---------------------------|---------------------|
| Aqua. <b>6059</b>         | Lime. <b>6062</b>   |
| Blue. <b>6054</b>         | Peach. <b>6064</b>  |
| Burnt Orange. <b>6056</b> | Pink. <b>6052</b>   |
| Chartreuse. <b>6061</b>   | Red. <b>6060</b>    |
| Gold. <b>6058</b>         | Rose. <b>6063</b>   |
| Green. <b>6053</b>        | Tan. <b>6055</b>    |
| Lavender. <b>6065</b>     | White. <b>6050</b>  |
| Lilac. <b>6057</b>        | Yellow. <b>6051</b> |



Simport



Simport  
Macrosette

- Tissue Cassettes, Attached Lids**, slotted pores
- Histosette™ I, snap-apart lids, Pk/1500 ..... \$147.00  
30° front labeling area.
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M490-6</b>  | Pink. <b>M490-3</b>  |
| Green. <b>M490-4</b> | White. <b>M490-2</b> |
- Histosette™ II, integral lids, Pk/1500 ..... \$199.00  
45° front labeling area.
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M492-6</b>  | Pink. <b>M492-3</b>  |
| Green. <b>M492-4</b> | White. <b>M492-2</b> |
- Unisette™, hinged lids, Pk/1500 ..... \$157.00  
35° front labeling area. *(Alternative to Sakura's Uni-Cassette®)*
- |                      |                      |
|----------------------|----------------------|
| Blue. <b>M505-6</b>  | Pink. <b>M505-3</b>  |
| Green. <b>M505-4</b> | White. <b>M505-2</b> |
- Macrosette™, 2X normal depth, white, Pk/750. **M512** ..... \$204.00  
Dimensions: 40.1 x 28.5 x 13 mm(H)

Continued

## Cassette Storage



Matchbox Style

**Box, Matchbox Style (Pill Box),** cardboard, white  
 $2\frac{1}{8} \times 1\frac{1}{8} \times \frac{5}{8}$ ", Pk/72. **20047** ..... \$33.80

**Box, Pizza Box Style,** 14 x 14 x 2", white, Pk/50. **8342** ..... \$121.00  
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks  
 in 4 drawers, each divided into 12 rows.  
 Cabinets can be stacked up to 8 high.  
 Dimensions: 49.5 x 38 x 21 cm  
 ( $19\frac{1}{2} \times 15 \times 8\frac{1}{4}$ ")



**Tissue Cassette Cabinet, Steel,** 4-drawer. **485400** ..... \$386.00

Store up to 2100 cassettes in 6 drawers.  
 Each drawer includes sponge block  
 spacers and an identification label.  
 Cabinets can be stacked up to 4 high.  
 Dimensions: 44.5 x 24 x 39.5 cm(H)  
 ( $17\frac{1}{2} \times 9\frac{1}{2} \times 15\frac{1}{2}$ ")



**Tissue Cassette Cabinet, ABS.** **97-0600** ..... \$183.00




Simport

Store 250 cassettes or 165 embedding  
 rings in each stackable drawer.  
 High-impact polystyrene resists  
 temperatures from -80°C to +80°C.  
 Drawer dimensions: 40.5 x 23 x 5.1 cm.

**Tissue Cassette Drawers, Polystyrene**  
 Pk/2. **M495-2** ..... \$75.60  
 Pk/6. **M495-6** ..... \$172.00


## Brushes



### Brushes

Flat Bristles (paint style)  
 Camel Hair, 1" wide, 6½" total length. **BR665B (LB-37)** ..... \$6.65  
 Round Bristles, camel hair  
 $\frac{3}{16}$ " Diameter (artist style), 7½" total length. **LB-31** ..... \$1.90  
 $\frac{3}{4}$ " Diameter, 8½" total length. **BR661 (LB-30)** ..... \$11.15

## Microtome Blades



### Microtome Blades

Accu-Edge®, Sakura Finetek  
 Low Profile, 80 x 8.3 mm, Pk/50. **4689** ..... \$189.00  
 High Profile, 75.7 x 14 mm, Pk/50. **4685** ..... \$208.00  
 Encore™, PTFE coated, **USA**  
 Low Profile, 83 x 8 mm, Pk/50. **7380** ..... \$84.00  
 High Profile, 83 x 14 mm, Pk/50. **7390** ..... \$84.00

## Scrapers

Remove paraffin, ice, etc., safely  
 from your laboratory equipment  
 with these plastic scrapers.



### Scrapers, plastic

5 cm (2") Wide Blade. **311-245** ..... \$3.70  
 10 cm (4") Wide Blade. **311-250** ..... \$4.70  
 15 cm (6") Wide Blade. **311-255** ..... \$5.70

## Tissue Flotation Bath



**Tissue Flotation Bath. XH-1001** ..... \$298.00