

SLIDE & CASSETTE PROCESSING

Microscope Slides

Slides, general, 75 x 25 x 1 mm

Plain, ground edges	
Premium, 90° corners, "Globe"	
Cs/20 Bx, bx/72. 1301	\$137.00
Frosted, beveled edges	
Economy, sandblasted, safety corners, C&A	
Cs/20 Bx, bx/72. 8205	\$150.00
Premium, etched, clipped corners, "Globe"	
Cs/20 Bx, bx/72. 1334	\$181.00
White End, clipped corners	
Economy, C&A	
Cs/20 Bx, bx/72. 9108W-E	\$152.00
Superfrost®, "SPECIMEN". 4951-001T	\$317.00
(12-550-143, 48311-600)	

Superfrost® Plus



Concavity, 2-Well



Slides, specialty

Charged (silanated), 75 x 25 x 1 mm	
Recommended for adhering tissue sections and cell preparations.	
Diamond™, white end, Globe	
Cs/20 Bx, bx/72. 1358W	\$400.00
Superfrost® Plus, white end, Erie (EpreDia)	
Cs/20 Bx, bx/72. 4951PLUS-001	\$748.00
Preferred slide for Ventana instrumentation.	
(12-550-15, 48311-703)	
Concavity	
75 x 25 x 1.3 mm, 16 mm dia., 0.8 mm deep, economy	
1-Well, Bx/12. 600CS3X11	\$7.55
2-Well, Bx/12. 600CS3X12	\$10.35
76 x 26 x 3.2 mm, 18 mm dia., 1.75 mm deep, Erie (EpreDia)	
1-Well, Bx/12. 1526	\$97.35
Fused Quartz, NSG Type OZ®, 3" x 1" x 1 mm, USA	
QMS3X1X1 , each	\$20.00
Fused Silica, Corning 7980, 3" x 1" x 1 mm, USA	
FSMS3X1X1 , each	\$29.00

Coverslips (Cover Glasses)

Coverslips, premium, borosilicate glass, 1 oz/pk, Cs/10 oz

#0 Thickness (0.08-0.12 mm)	
Circles, Germany*	
25 mm, ~222 pcs/oz, Cs/10 oz. 0-25	\$214.00

#1 Thickness (0.13-0.17 mm)	
Squares & Rectangles, Erie (EpreDia)	
22 x 22 mm, ~162 pcs/oz, Cs/10 oz. 1-2222	\$70.15
22 x 30 mm, ~119 pcs/oz, Cs/10 oz. 1-2230	\$63.30
22 x 50 mm, ~71 pcs/oz, Cs/10 oz. 1-2250	\$63.30
24 x 50 mm, ~65 pcs/oz, Cs/10 oz. 1-2450	\$63.30
24 x 60 mm, ~54 pcs/oz, Cs/10 oz. 1-2460	\$63.30

#1 Thickness (0.13-0.17 mm)	
Circles, Germany*	
12 mm, ~663 pcs/oz, Cs/5 oz. 1-12	\$249.00
18 mm, ~296 pcs/oz, Cs/5 oz. 1-18	\$150.00

#1.5 Thickness (0.16-0.19 mm). Recommended for most users	
Squares & Rectangles, Erie (EpreDia)	
18 x 18 mm, ~209 pcs/oz, Cs/10 oz. 1.5-1818	\$67.75
22 x 22 mm, ~140 pcs/oz, Cs/10 oz. 1.5-2222	\$69.35
22 x 50 mm, ~62 pcs/oz, Cs/10 oz. 1.5-2250	\$63.30
24 x 40 mm, ~71 pcs/oz, Cs/10 oz. 1.5-2440	\$63.30
24 x 50 mm, ~56 pcs/oz, Cs/10 oz. 1.5-2450	\$63.30
24 x 60 mm, ~47 pcs/oz, Cs/10 oz. 1.5-2460	\$63.30

#1.5 Thickness (0.16-0.19 mm). Recommended for most users	
Circles, Germany*	
12 mm, ~552 pcs/oz. 1.5-12	cs/5 oz \$198.00 cs/10 oz —
18 mm, ~245 pcs/oz. 1.5-18	\$250.00 —
22 mm, ~164 pcs/oz. 1.5-22	— \$252.00
25 mm, ~127 pcs/oz. 1.5-25	— \$252.00

*German glass is preferred for cell growth/adhesion

Coverslips, specialty

Fused Quartz, USA	
0.15-0.25 mm Thick, NSG Type OZ®	
1" Square, each. QCS1X1	\$19.00
1" Circle, each. QCS1X0	\$19.00
0.16 ± 0.01 mm Thick, A1	
25 mm Square, each. R525025	\$32.40
25 mm Circle, each. R525000	\$23.40
Fused Silica, USA	
0.15-0.25 mm Thick, Corning 7980	
1" Square, each. FSCS1X1	\$21.50
1" Circle, each. FSCS1X0	\$21.50
0.16 ± 0.01 mm Thick, S1-UV	
25 mm Square, each. R425025	\$32.00
25 mm Circle, each. R425000	\$28.80