



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).

Plastic (polyoxymethylene [POM]), holds 25 slides  
Set (dish, rack & lid). **0990-99** .....\$50.90



674



60062101



675

**Staining Jars**, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, USA. <b>673</b> .....	\$13.10	146.00
Glass Cover, Fits 673 Jar. <b>673-LID</b> .....	\$4.85	—
	<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, <b>USA 674</b> .....	\$13.10	146.00
Phenolic Cap, Fits 674 Jar, black. <b>674-LID</b> .....	\$4.95	—
	<u>ea</u>	<u>cs/12</u>
PP w/Screw-On Cap. <b>60062101</b> .....	\$4.90	47.40
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
	<u>ea</u>	<u>cs/12</u>
Wide Top, Glass w/Square Lid, <b>USA 675</b> .....	\$17.65	193.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. <b>675-LID</b> .....	\$5.40	—



12-Slide

Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)  
12-Slide, EasyDip™, Simport

Jar only, Blue, Cs/6. <b>M900-12B</b> .....	\$66.55
Jar only, Green, Cs/6. <b>M900-12G</b> .....	\$66.55
Jar only, Pink, Cs/6. <b>M900-12P</b> .....	\$66.55
Jar only, White, Cs/6. <b>M900-12W</b> .....	\$66.55
Jar only, Yellow, Cs/6. <b>M900-12Y</b> .....	\$66.55
Rack/Lid only, gray, Cs/6. <b>M905-12DGY</b> .....	\$73.55



Anodized aluminum staining jar holder

(Jars sold separately)

Staining Jar Holder, 6-Place. **M906**.....\$81.45




4501030



4501042

24-Slide, light green, microwaveable, acetal  
Jar & Lid. **4501030** .....\$11.55  
Rack only, gray. **4501042** .....\$21.45  
Staining Jar Station, 12-Place. **4501011** .....\$221.00  
Includes: metal frame, (12) 24-slide dishes w/ lids (4501030) & (1) rack (4501042).  
Dimensions: 29(L) x 5(W) x 5"(H)

**Tissue (Embedding) & Biopsy Cassettes**




**Base Molds**, Plastic, disposable, Simport


15 x 15 x 5 mm, Cs/1000. <b>M475-2</b> .....	\$152.00
24 x 24 x 5 mm, Cs/1000. <b>M475-3</b> .....	\$152.00
30 x 24 x 5 mm, Cs/1000. <b>M475-4</b> .....	\$152.00
37 x 24 x 5 mm, Cs/1000. <b>M475-5</b> .....	\$152.00
37 x 24 x 10 mm, Cs/500. <b>M475-10</b> .....	\$78.35

Designed for Simport Macrosette™.

Solvent resistant,  
0.2 x 0.2 mm mesh.



Nylon




Paper

**Biopsy Bags**  
Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011** .....\$235.00  
Paper (Tea Bags), 2¼ x 2¼", Pk/500. **56100** .....\$46.50

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



M476-2

**Biopsy Pads**  
Rectangular  
Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014** .....\$42.00  
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1** .....\$217.00  
Circular, Simport  
Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2** .....\$266.00