


Biopsy Cassettes



Simport
Histosette I



Simport
Micromesh

- Biopsy Cassettes, Attached Lids, acetyl polymer, Simport**
- Histosette™ I, 1 mm square pores, Pk/1500\$150.00
30° front labeling area.
- Blue. **M491-6** Pink. **M491-3**
Green. **M491-4** White. **M491-2**
- Unisette™, hinged lids, Pk/1500\$160.00
35° front labeling area. (*Alternative to Sakura's Uni-Cassette®*)
- Blue. **M506-6** Pink. **M506-3**
Green. **M506-4** White. **M506-2**
- Micromesh™, 0.38 mm pores, Pk/1000\$245.00
45° front labeling area.
- Blue. **M507-6** Pink. **M507-3** Yellow. **M507-5**
Green. **M507-4** White. **M507-2**

Tissue Cassettes


Tissue
Cassettes



- Tissue Cassettes, Detached Lids, slotted pores**
- Histosette™ II, Simport, Pk/1000\$104.00
45° front labeling area.
- Blue. **M485-6** Pink. **M485-3**
Green. **M485-4** White. **M485-2**
- Uni-Cassettes™, **USA**, Pk/500\$75.50
- Aqua. **6059** Lime. **6062**
Blue. **6054** Peach. **6064**
Burnt Orange. **6056** Pink. **6052**
Chartreuse. **6061** Red. **6060**
Gold. **6058** Rose. **6063**
Green. **6053** Tan. **6055**
Lavender. **6065** White. **6050**
Lilac. **6057** Yellow. **6051**



Simport



Simport
Macrosette

- Tissue Cassettes, Attached Lids, slotted pores**
- Histosette™ I, snap-apart lids, Pk/1500\$150.00
30° front labeling area.
- Blue. **M490-6** Pink. **M490-3**
Green. **M490-4** White. **M490-2**
- Histosette™ II, integral lids, Pk/1500\$199.00
45° front labeling area.
- Blue. **M492-6** Pink. **M492-3**
Green. **M492-4** White. **M492-2**
- Unisette™, hinged lids, Pk/1500\$160.00
35° front labeling area. (*Alternative to Sakura's Uni-Cassette®*)
- Blue. **M505-6** Pink. **M505-3**
Green. **M505-4** White. **M505-2**
- Macrosette™, 2X normal depth, white, Pk/750. **M512**\$215.00
Dimensions: 40.1 x 28.5 x 13 mm(H)

Continued