

GENERAL LAB SUPPLIES & EQUIPMENT


Applicators



Applicators, 6" long

- Cotton Tipped, wooden shaft
- Nonsterile, 10 bg/bx, Bx/1000. **95-8702**\$14.05
- Sterile Pk/1, Bx/100. **44024-106-1S**\$7.85
- Sterile Pk/2, Bx/200. **44024-106-2S**\$9.80
- Plain, wood, Bg/1000. **95-8701**\$8.60

Aprons



Protect yourself with this disposable 0.8-mil thick apron.

Aprons, disposable, polyethylene, 46"(L) x 28"(W)

- 0.8 mil Thick (0.02 mm), Bx/100. **105DA2846**\$16.65
- 1.75 mil Thick (0.044 mm), Bx/100. **59010485**\$38.00

Autoclaves

B4000-28-120



B4000-M

B4000-18


Autoclaves, benchtop, Benchmark (1-year warranty)

- 8 L, BioClave™ Mini. **B4000-M**\$4,399.00↓
Chamber Dimensions: 6.7" ID x 12.5" (D)
- 18 L, BioClave™ 18. **B4000-18**\$6,019.00↓
Chamber Dimensions: 9.8" ID x 13.75" (D)
- 28 L, BioClave™ 28. **B4000-28-120**\$11,080.00
Chamber Dimensions: 12.6" ID x 16.4" (D)

Bags



Collect biohazardous waste in this space-saving, autoclavable bag. Printed with a biohazard symbol and sterilization indicator patch. Measures 20.3 x 30.5 cm (8 x 12") and is 1.57 mil (0.04 mm) thick polypropylene.



HS10044H

- Bags, Biohazard, Benchtop, 8 x 12", Bx/500. HS10320**\$102.00
- Biohazard Bag Holder, Benchtop. **HS10044H**\$27.00
Epoxy-coated wire holder measures 9.1 x 12.8 x 21.4 cm.



Bags display a biohazard symbol and precautionary procedures. Sterilization indicator patch darkens when autoclaved. Polypropylene bags are bright red. Durable, extra heavy-duty bags provide maximum protection against dangerous tears and leaks.

Bags, Biohazard, Standard, 1.5 mil (0.038 mm) thick

- 14 x 19" (36 x 48 cm), Bx/200. **HS10321**\$101.00
- 19 x 23" (48 x 58 cm), Bx/200. **131641923**\$144.00
- 25 x 35" (64 x 89 cm), Bx/200. **131642535**\$309.00

Bags, Biohazard, Heavy-Duty, 2.2 mil (0.055 mm) thick

- 12 x 18"(L), 9 L, Bx/200. **600UNU40105**\$62.00
- 12 x 24"(L), 12 L, Bx/200. **600UNU40106**\$74.25
- 18 x 22"(L), 16 L, Bx/200. **600UNU40103**\$89.00
- 25 x 35"(L), 37 L, Bx/200. **600UNU40104**\$155.00



Bag, Waste Disposal, HDPE, Cleanware™

- 8 x 10", Pk/50. **131741008**\$29.35



Seal polyethylene or polypropylene bags, sheets, or flat tubing from 0.01 to 0.8 mm thick. Sliding blade cuts excess material 19 mm from the seal.

- Bag Sealer, electric, 12" (30.5 cm) Wide. HS28272**\$248.00
- Bag Sealer Replacement Part Kit. **HS28280**\$37.80