

## Balances

### Portable



<b>Model</b>	<b>DWP-1001</b>
<b>Capacity (g)</b>	1000
<b>Readability (g)</b>	0.1
<b>Modes</b>	g, oz, ct, gn
<b>Pan</b>	14.6 cm (5.75") stainless steel
<b>Power</b>	4 AA/AC adapter
<b>Price</b>	\$72.00

Warranty: 1 year

### Toploading

Benchmark  
W3200-1200



Model	W3200-500	W3200-1200
<b>Capacity (g)</b>	500	1200
<b>Readability (g)</b>	0.001	0.01
<b>Modes</b>	g, oz, ct	g, oz, ct
<b>Pan</b>	11.5 cm (4.5")	16 cm (6.3")
<b>Draft Shield</b>	Yes, 3-door glass	No
<b>Price</b>	\$1,299.00	\$1,053.00

Warranty: 2 years

### Analytical

Benchmark



Model	★ W3100-120	W3100A-120
<b>Calibration</b>	External (100 g weight provided)	Internal (Quick-Cal™ software)
<b>Capacity (g)</b>	120	120
<b>Readability (g)</b>	0.0001	0.0001
<b>Modes</b>	g, oz, lb, ct	g, oz, lb, ct
<b>Pan</b>	9 cm (3.5")	9 cm (3.5")
<b>Price</b>	\$1,948.00 ↓	\$2,067.00 ↓

Warranty: 2 years

## Baths



BSH1002

Aluminum  
Blocks


**Dry Baths**, digital, Benchmark (2-year warranty)

1-Block Capacity. <b>BSH1001</b> .....	\$527.00
2-Block Capacity. <b>BSH1002</b> .....	\$649.00
4-Block Capacity. <b>BSH1004</b> .....	\$1,029.00

*Blocks sold separately – please specify configuration*

Blocks, aluminum

Holds (48) 0.2 mL Tubes. <b>BSW02</b> .....	\$152.00
Holds (24) 0.5 mL Tubes. <b>BSW05</b> .....	\$152.00
Holds (24) 1.5 mL Tubes. <b>BSW1500</b> .....	\$152.00
Holds (24) 1.5/2.0 mL Tubes. <b>BSW1520</b> .....	\$152.00
Holds (20) 12/13 mL Tubes. <b>BSW13</b> .....	\$152.00
Holds (12) 15 mL Tubes. <b>BSW15</b> .....	\$152.00
Holds (5) 50 mL Tubes. <b>BSW50</b> .....	\$152.00



40 x 1.5/2.0 mL



15 x 50 mL



41 x 15 mL

Racks are only included with model B2000-4.

**Water Baths**, digital, MyBath™, Benchmark (2-year warranty)

4 L, 3 tube racks. <b>B2000-4</b> .....	\$1,009.00
8 L. <b>B2000-8</b> .....	\$1,115.00
12 L. <b>B2000-12</b> .....	\$1,329.00