

## Beakers

Polypropylene Beakers



### Beakers

Glass (borosilicate), low form, double scale

Economy

10 mL, no graduations, Pk/12 <b>600BG10</b> .....	\$32.40
50 mL, Pk/12. <b>600BG50</b> .....	\$39.10
100 mL, Pk/12. <b>600BG100</b> .....	\$39.45
150 mL, Pk/12. <b>600BG150</b> .....	\$42.25
250 mL, Pk/12. <b>600BG250</b> .....	\$45.45
400 mL, Pk/12. <b>600BG400</b> .....	\$57.25
600 mL, Pk/6. <b>600BG600</b> .....	\$41.90
800 mL, Pk/6. <b>600BG800</b> .....	\$75.00
1000 mL, Pk/6. <b>600BG1000</b> .....	\$67.50
2000 mL, each. <b>600BG2000</b> .....	\$21.20
4000 mL, each. <b>600BG4000</b> .....	\$63.75

Heavy-Duty

50 mL, Pk/12. <b>600BG100350</b> .....	\$40.60
100 mL, Pk/12. <b>600BG1003100</b> .....	\$73.00
150 mL, Pk/12. <b>600BG1003150</b> .....	\$81.00
250 mL, Pk/12. <b>600BG1003250</b> .....	\$93.50
400 mL, Pk/12. <b>600BG1003400</b> .....	\$109.00
600 mL, Pk/6. <b>600BG1003600</b> .....	\$81.00
1000 mL, each. <b>600BG100310000</b> .....	\$19.65

Polypropylene, premium, low form, single scale

100 mL, 25 mL graduations, each. <b>262110000</b> .....	\$3.95
250 mL, 50 mL graduations, each. <b>262130000</b> .....	\$6.55
600 mL, 100 mL graduations, each. <b>262150000</b> .....	\$10.20
1000 mL, 200 mL graduations, each. <b>262160000</b> .....	\$10.35
2000 mL, 500 mL graduations, each. <b>262170000</b> .....	\$31.35
4000 mL, 1000 mL graduations, each. <b>262190000</b> .....	\$52.25

**NEW**

## Bench Protector



Protect your bench surface using this 40 mil absorbent liner backed with 0.5 mil polyethylene film.

**Bench Liner**, Roll, 50' x 20", white, Labmat™. **246750000** ..... \$89.50

## Bottles




Deliver controlled drops from standard 0.030" tip opening with these LDPE bottles. Sealer cap snaps on easily.

### Bottles, Dispensing, squeeze bottles

LDPE (low-density polyethylene), Yorker cap

30 mL (1 oz), Bg/12. <b>116370001</b> .....	\$27.75
185 mL (6 oz), Bg/12. <b>116370006</b> .....	\$49.70

A long nozzle, textured sides, and slim profile make this chemically inert, HDPE container ideal for dispensing immersion oil.



HDPE (high-density polyethylene)

20 mL, each. <b>18145</b> .....	\$8.10
---------------------------------	--------

Amber glass protects light-sensitive liquids. Screw-cap dropper features a black rubber bulb and glass pipet.



Boston Round

### Bottles, Dropping, Amber Glass (flint)

30 mL (1 oz), Bx/12. <b>3104053</b> .....	\$18.35
Bottle Dropper Assembly (replacement)	
20/400, 7 x 76 mm Pipet, Pk/12. <b>3104053D</b> .....	\$10.75
60 mL (2 oz), Bx/12. <b>3104054</b> .....	\$26.85
Bottle Dropper Assembly (replacement)	
20/400, 7 x 89 mm Pipet, Pk/12. <b>3104054D</b> .....	\$11.00

French Square



### Bottles, Dropping, Clear Glass (flint)

15 mL (1/2 oz), Bx/12. <b>3104062</b> .....	\$33.80
Bottle Dropper Assembly (replacement)	
20/400, 7 x 60 mm Pipet, Pk/12. <b>3104074</b> .....	\$11.55
30 mL (1 oz), Bx/12. <b>3104063</b> .....	\$38.50
Bottle Dropper Assembly (replacement)	
24/400, 7 x 75 mm Pipet, Pk/12. <b>3104075</b> .....	\$14.70

Continued