


Hinged
Polystyrene



Boxes, Plastic, (see following chart)

Plastic Boxes				
Catalog #	Inside Dimensions	Foam Insert	Suggested Application	Price
H-1	1 x 1 x 1/4"	-	Reticles, Pointers	\$0.40
H-1/11B	1 x 1 x 1/4"	1	Reticles, Pointers	\$0.50
H-54	2 9/16 x 2 9/16 x 1/2"	-	Filters	\$1.20
H-54/2916G	2 9/16 x 2 9/16 x 1/2"	1	Filters	\$1.35
H-28/334B2	3 3/4 x 1 3/16 x 5/8"	2	Slides, Stage Micrometers	\$1.40
H-10	2 7/8 x 1 3/16 x 3/8"	-	Small Parts	\$0.60
H-21/512B	5 1/2 x 1 3/4 x 1/2"	1	Forceps & Scissors	\$1.80

Brushes



Brushes

- Bench, horse hair, 10" long. **BRS-969.10**\$9.95
Bristle head: 4 3/4"(L) x 1 1/2"(W); bristle length: 1"
- Flat Bristles (paint style)
Camel Hair, 1" wide, 6 1/2" total length. **BR665B (LB-37)**\$6.65
- Horsehair, heavy-duty bristles, metal handle. **33761880**\$0.50
- Round Bristles, camel hair, 8 1/2" long. **BR661 (LB-30)**\$11.25
- Swivel-Head, black nylon bristles. **BR565 (LB-14)**\$4.70
Bristle head: 4 1/4"(L) x 1 1/2"(W); total length: 12"

Continued

Remove debris, rust, or corrosion in hard to reach places.



- Pen Style (scratch), retractable, 11.5 cm (4 1/2") long, Germany
Brass. **BRS-295.00**\$12.20
- Nylon. **BRS-293.00**\$9.70
- Steel. **BRS-292.00**\$12.80

Burners



- Burner, Alcohol**, 120 mL, glass body, metal cap. **97-5313**\$7.50
- Wicks (replacement), 10-14 cm x 8 mm, Pk/12. **97-5310**\$5.45



H-5870 vertical mixing tube is equipped with a flame stabilizer and adjustable air supply ports. Serrated gas inlet accepts 5/16" ID tubing. Height is 15.6 cm (6 1/8")

H-5550 gridded top, 3.8 cm (1 1/2") diameter, produces a uniform blue cone of light. Height is 21.6 cm (8 1/2")



Burners, Bunsen, natural gas, Humboldt Mfg, **USA**

- Standard, flame stabilizer. **H-5870**\$32.65
- High Temperature, Meeker-type. **H-5550**\$98.80
- Accessory
Bunsen Burner Natural Gas Connector
3' Long. **H-2087**\$39.15



- Burner Lighter**, replaceable single flint, round file. **SOL-290.00**\$3.95
- Flints (replacement), Pk/5. **SOL-290.01**\$2.85

Carboys



- Carboys with Spigot**, polypropylene, round
10 L (2.5 gal). **60034051**\$80.95
- 20 L (5 gal). **60034052**\$112.00