

Drying Racks



188050000

189330023



HS23243A

Drying Racks, Benchtop, Bel-Art

Single-sided. **188050000**\$170.00

Dry labware on this epoxy-coated steel wire rack with polystyrene tray. Holds up to 26 pieces on 12 straight pins and 14 hairpin loops. Dimensions: 43 x 17 x 48 cm(H).

Double-sided

2-Tier. **189330022**\$391.00

3-Tier. **189330023**\$539.00

NEW

Drying Rack, Benchtop/Wallmount, Heathrow

Single-sided. **HS23243A**\$146.00

Dry up to 90 pieces of labware with this self-standing or wall-mounted rack. Constructed of vinyl-coated wire hairpin loops with a polystyrene tray and rack. Dimensions: 46.2 x 18.2 x 52.5 cm(H).

Electrophoresis



Electrophoresis System, MyGel™ Mini. E1101\$607.00

Eyewear

Absorbs
99.9% UV



Glasses, Safety, integral side shield. **43574002**\$6.00



Flexible frame and adjustable elastic strap provides a tight seal.

Goggles, Safety, vented (chemical splash). **97-2001**\$4.00

Filter & Blotting Papers



Ahlstrom - Munksjo
Filter Paper

Qualitative, Standard Grades (for gravity filtration)

Slow flow (4 mL/min), 1.5 µm retention, Grade 610

(Whatman 5; S&S 602)

9 cm, Pk/100. **6100-0900**\$12.20

12.5 cm, Pk/100. **6100-1250**\$15.25

15 cm, Pk/100. **6100-1500**\$18.40

20.5 cm, Pk/100. **6100-2050**\$29.05

24 cm, Pk/100. **6100-2400**\$40.45

Medium flow (30 mL/min), 2 µm retention, Grade 642

(Whatman 2, 2CHR, 20CHR; S&S 2043B, 2045A/B)

9 cm, Pk/100. **6420-0900**\$9.50

15 cm, Pk/100. **6420-1500**\$17.05

Medium flow (55 mL/min), 2.5 µm retention, Grade 601

(Whatman 1, 1CHR; S&S 593, 597, 2043A)

4.25 cm, Pk/100. **6010-0425**\$7.25

9 cm, Pk/100. **6010-0900**\$9.80

12.5 cm, Pk/100. **6010-1250**\$14.10

15 cm, Pk/100. **6010-1500**\$17.15

18.5 cm, Pk/100. **6010-1850**\$24.85

20.5 cm, Pk/100. **6010-2050**\$26.80

33 cm, Pk/100. **6010-3300**\$72.40



Medium flow (60 mL/min), 6 µm retention, Grade 613

(Whatman Student; S&S 595, 0903)

9 cm, Pk/100. **6130-0900**\$8.35

15 cm, Pk/100. **6130-1500**\$14.50

20.5 cm, Pk/100. **6130-2050**\$21.80

33 cm, Pk/100. **6130-3300**\$48.45

Medium flow (200 mL/min), 10 µm retention, Grade 631

(Whatman 4, 4CHR; S&S 604, 2040A/B; Fisher P8)

20.5 cm, Pk/100. **6310-2050**\$27.00

33 cm, Pk/100. **6310-3300**\$73.15

Fast flow (235 mL/min), 25 µm retention, creped, Grade 615

9 cm, Pk/100. **6150-0900**\$8.60

12.5 cm, Pk/100. **6150-1250**\$12.50

15 cm, Pk/100. **6150-1500**\$15.00

20.5 cm, Pk/100. **6150-2050**\$22.50

33 cm, Pk/100. **6150-3300**\$49.90

Very Fast flow (360 mL/min), 35 µm retention, creped, Grade 617

9 cm, Pk/50. **6170-0900**\$7.85

12.5 cm, Pk/50. **6170-1250**\$10.45

15 cm, Pk/50. **6170-1500**\$12.15

20.5 cm, Pk/50. **6170-2050**\$18.55

33 cm, Pk/50. **6170-3300**\$39.80

Qualitative, Pleated (Fluted)

Slow flow (30 mL/min), 2 µm retention, Grade 542

(Whatman 2V)

12.5 cm, Pk/100. **5420-1250**\$78.15

18.5 cm, Pk/100. **5420-1850**\$88.45

Continued