

Qualitative, Wet-Strengthened Grades (for vacuum filtration)

- Medium flow** (65 mL/min), 3 µm retention, **Grade 909**
 (Whatman 91; S&S 0860, 1574, 3003; Reeve-Angel 201)
 32 cm, Pk/100. **9090-3200**\$78.85
- Fast flow** (245 mL/min), 27 µm retention, creped, **Grade 962**
 (Whatman 114; S&S 0856, 0905, 3002; Reeve-Angel 202)
 11 cm, Pk/100. **9620-1100**\$11.70

Quantitative, Ashless

(Designed for use in analytical and gravimetric analysis)

- Slow flow** (4 mL/min), 1.5 µm retention, **Grade 94**
 (Whatman 42; S&S 589/3)
 9 cm, Pk/100. **0940-0900**\$33.00
 15 cm, Pk/100. **0940-1500**\$56.85
- Medium flow** (30 mL/min), 2 µm retention, **Grade 74**
 (Whatman 40; S&S 589/2, 589/4)
 9 cm, Pk/100. **0740-0900**\$33.00
 15 cm, Pk/100. **0740-1500**\$56.85
- Fast flow** (220 mL/min), 15 µm retention, **Grade 55**
 (Whatman 54, 541; S&S 1505)
 9 cm, Pk/100. **0550-0900**\$37.90
 15 cm, Pk/100. **0550-1500**\$59.95

Blotting/Chromatography/Electrophoresis Papers


- 0.19 mm thick**, 125 mm/30 min absorption, **Grade 601**
 (Whatman 1, 1CHR; S&S 593, 597, 2043A)
 46 x 57 cm, Pk/100. **6018-4657**\$147.00
- 0.34 mm thick**, 140 mm/30 min absorption, **Grade 238**
 (Whatman 3MM, S&S 598, 231)
 20 x 20 cm, Pk/100. **2388-2020**\$75.25
 46 x 57 cm, Pk/100. **2388-4657**\$286.00
- 0.42 mm thick**, 150 mm/30 min absorption, **Grade 237**
 (Whatman 3, 3CHR; S&S 22/2, 591, 2469, 5703;
 Albet 2040B, FP591; Munktell 3S/H, FN7a)
 46 x 57 cm, Pk/100. **2378-4657**\$329.00
- 0.95 mm thick**, 180 mm/30 min absorption, **Grade 243**
 (Whatman 17CHR; S&S 2668)
 11 x 14 cm, Pk/100. **2438-1114**\$31.75
 20 x 20 cm, Pk/100. **2438-2020**\$62.70
- 2.48 mm thick**, 240 mm/30 min absorption, **Grade 320**
 (S&S 2727)
 11 x 14 cm, Pk/50. **3208-1114**\$60.10
 20 x 20 cm, Pk/50. **3208-2020**\$87.45



Smooth tips & polypropylene grips

- Forceps, Filter, 11.5 cm. 300-025**\$4.40

Freezer Racks



Freezer Racks, stainless steel, **USA**

175V-13-2L

Restocking fee of 20% on all returned freezer racks listed below.
 No returns on custom racks.

Horizontal (fits upright freezers)

- Horizontal, Holds 2" Boxes**
- | | | |
|---|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). 175HS-9-2 | \$59.00 | 53.00 |
| Dimensions: 16.31(L) x 5.5(W) x 6.5"(H) | | |
| 12-Place (3L x 4H). 175HS-12-2 | \$78.00 | 71.00 |
| Dimensions: 16.31(L) x 5.5(W) x 8.75"(H) | | |
| 15-Place (3L x 5H). 175HS-15-2 | \$98.00 | 89.00 |
| Dimensions: 16.31(L) x 5.5(W) x 11"(H) | | |
| 16-Place (4L x 4H). 175HL-16-2 | \$105.00 | 95.00 |
| Dimensions: 22(L) x 5.5(W) x 8.75"(H) | | |
| 20-Place (4L x 5H). 175HL-20-2 | \$132.00 | 119.00 |
| Dimensions: 22(L) x 5.5(W) x 11"(H) | | |
- Horizontal, Holds 3" Boxes**
- | | | |
|---|---------|-------|
| 9-Place (3L x 3H). 175HS-9-3 | \$61.00 | 55.00 |
| Dimensions: 16.31(L) x 5.5(W) x 9.5"(H) | | |
| 12-Place (4L x 3H). 175HL-12-3 | \$80.00 | 73.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |



Horizontal Drawer (fits upright freezers)

- Horizontal Drawer, Holds 2" Boxes**
- | | | |
|--|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 9-Place (3L x 3H). 175HD-332 | \$85.00 | 77.00 |
| Dimensions: 16.5(L) x 5.5(W) x 7.125"(H) | | |
| 12-Place (3L x 4H). 175HD-342 | \$111.00 | 101.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H) | | |
| 15-Place (3L x 5H). 175HD-352 | \$141.00 | 128.00 |
| Dimensions: 16.5(L) x 5.5(W) x 11.812"(H) | | |
| 16-Place (4L x 4H). 175HD-442 | \$151.00 | 137.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |
| 20-Place (5L x 4H). 175HD-542 | \$188.00 | 171.00 |
| Dimensions: 27.25(L) x 5.5(W) x 9.5"(H) | | |
- Horizontal Drawer, Holds 5" Boxes**
- | | | |
|---|------------|------------|
| | <u>1-9</u> | <u>≥10</u> |
| 6-Place (3L x 2H). 175RD-325 | \$74.00 | 67.00 |
| Dimensions: 18.125(L) x 6(W) x 10.75"(H) | | |
| 8-Place (4L x 2H). 175RD-425 | \$99.00 | 90.00 |
| Dimensions: 24(L) x 6(W) x 10.75"(H) | | |
- Horizontal Drawer, Holds 15 mL tubes**
- | | | |
|---|----------|--------|
| 160-Place (80L x 2H). 175HD-32ML-15 | \$99.00 | 90.00 |
| Dimensions: 16.5(L) x 5.5(W) x 9.5"(H) | | |
| 220-Place (110L x 2H). 175HD-42ML-15 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |
- Horizontal Drawer, Holds 50 mL tubes**
- | | | |
|---|----------|--------|
| 84-Place (42L x 2H). 175HD-42ML-50 | \$129.00 | 118.00 |
| Dimensions: 22(L) x 5.5(W) x 9.5"(H) | | |

Continued