

Horizontal Drawer-Tray (fits upright freezers)


Drawer-Tray, Holds 2" Boxes , aluminum trays	1-9	≥10
12-Place (3L x 4H). 175HDT-342	\$96.00	86.00
Dimensions: 16.63(L) x 5.5(W) x 8.5"(H)		
18-Place (3L x 6H). 175HDT-362	\$143.00	129.00
Dimensions: 16.63(L) x 5.5(W) x 12.75"(H)		
24-Place (4L x 6H). 175HDT-462	\$189.00	172.00
Dimensions: 22(L) x 5.5(W) x 12.75"(H)		

Vertical (fits chest freezers)

Vertical, Holds 2" Boxes		
13-Place. 175V-13-2L	\$88.00	79.00
Dimensions: 5.5(L) x 5.625(W) x 28.562"(H)		

Tube Storage Boxes (cryogenic) – see in GENERAL LAB SUPPLIES & EQUIPMENT section

Funnels



Funnels

Disposable, polypropylene, short stems, Simport		
★ 55 mm Top ID, takes 11 cm filter, Cs/100. F490-3	\$40.25	
★ 65 mm Top ID, takes 12.5 cm filter, Cs/100. F490-4	\$40.25	
Glass (borosilicate), long stem (150 mm)		
50 mm Top ID, Pk/6. 600GF614050	\$45.60	
65 mm Top ID, Pk/6. 600GF614065	\$54.90	
75 mm Top ID, Pk/6. 600GF614075	\$56.75	
Heavy-Duty, Molded-In Handle, polypropylene		
250 mm Top ID, 2.5 L capacity, each. 60057232	\$25.70	
Micro, polypropylene, Bel-Art		
24 mm Top ID, 2.6 mL capacity, Pk/12. 1146850024	\$41.05	
28 mm Top ID, 5.2 mL capacity, Pk/12. 146850028	\$33.60	
43 mm Top ID, 15.7 mL capacity, Pk/12. 146850043	\$70.80	
Powder, polypropylene		
Parallel Stem		
78 mm Top ID, Pk/12. 60057222	\$16.45	
100 mm Top ID, Pk/6. 146600100	\$15.95	
145 mm Top ID, Pk/6. 60057224	\$20.70	
Tapered Stem		
152 mm Top ID, each. 146820000	\$34.80	

Continued

Funnels (cont)

Standard Stem, Premium, polypropylene, Bel-Art	
Exclusive filter-step internal design provides speedier filtration; external ribs prevent air-lock.	
55 mm Top ID, takes 9 cm filter, Pk/12. 146760055	\$47.75
65 mm Top ID, takes 11 cm filter, Pk/12. 146760065	\$49.60
75 mm Top ID, takes 12.5 cm filter, Pk/6. 146760075	\$25.65
100 mm Top ID, takes 18.5 cm filter, Pk/4. 146760100	\$32.60
Standard Stem, Economy, polypropylene	
50 mm Top ID, each. 310FUNPP50	\$0.80
62 mm Top ID, each. 310FUNPP62	\$0.95
75 mm Top ID, each. 310FUNPP75	\$1.15
100 mm Top ID, each. 310FUNPP100	\$1.65

Gauze Pads

These 4 x 4" nonsterile gauze sponges are 12-ply.



Gauze Pads , 10 pk/200, Cs/2000. 44122000	\$88.45
---	---------

Germicidal Wipes


Sani-Cloth Plus	
Germicidal Wipes , Tub/160 440Q89072	\$10.00



Gloves & Hand Protection

Cotton Gloves	
Protect optics and photographs from fingerprints. Can also be used as glove liners.	
Small, Pk/12. GLV-190.10	\$5.25
Large, Pk/12. GLV-190.20	\$5.75
Purple Nitrile, powder-free, 5.5 mil, textured fingers, Halyard	
Nitrile gloves are hypoallergenic and provide greater puncture and chemical resistance than latex or vinyl gloves.	
Small, bx/100, Cs/10 Bx. 55081	\$139.00
Medium, bx/100, Cs/10 Bx. 55082	\$139.00
Large, bx/100, Cs/10 Bx. 55083	\$139.00
X-Large, bx/90, Cs/10 Bx. 55084	\$139.00

Nitrile Gloves



Gloves

Cotton

Continued