

## Student Microscope




### National

**131-RLED** Student Microscope, Monocular, DIN.....\$267.00

- Widefield 10X eyepiece w/ pointer accepts optional reticle
- Monocular, 45° inclined viewing head rotates 360°
- DIN achromat 4X, 10X, and 40XR objectives
- Stage provided with safety stop to protect slides & objectives
- Low positioned coarse and fine focus; tension adjustment
- In-stage 0.65 N.A. condenser with 5-place disc diaphragm
- LED illuminator, rechargeable batteries (equivalent to a 20 watt tungsten bulb)
- 36 cm (14") tall
- Lifetime mechanical and optical/3-year electrical warranty
- Grades K

## Stereo Microscopes



440-440PLL

### National

**440-440PLL** Stereoscope, Bino cular, 7.5-45X Zoom .....\$899.00

- Widefield 10X/20 eyepieces; accept 23 mm reticle
- Binocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-45X total magnification; working distance of 110 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual illuminators offer light from top only (3W LED with rheostat control and directional beam), bottom only (3W LED), or both top and bottom
- Post stand with locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- 24.2 cm(D) x 17.1 cm(W) x 43.3 cm(H)
- Limited lifetime mechanical and optical/3-year electrical

**440T-440PLL** Stereoscope, Trinocular, 7.5-45X Zoom .....\$1,029.00

- Trinocular version of 440-440PLL-10

M30TZ-SM99CL



### Swift

**M30TZ-SM99CL** Stereoscope, Trinocular, 7.5-45X Zoom \$1,095.00

- Widefield 10X/23 eyepieces; accept 24.95 mm reticle
- Trinocular, 45° inclined head with dual diopters rotates 360°; 55-75 mm interpupillary adjustment
- Paired achromat objectives provide continuous 7.5-45X total magnification; working distance of 113 mm with supplied 1X objective
- Heavy-duty rack and pinion focusing features a slip clutch and tension adjustment
- Dual 3W LED illuminators offer cool light from top only (directional beam), bottom only (variable control), or both top and bottom
- Post stand has a large ergonomic base and locking support collar
- Frosted glass and reversible black/white plastic stage plates included
- Lifetime mechanical and optical/3-year electrical warrant

## Fiber Optic Illuminator

*Please call regarding availability of items not listed*



Fiber optic light source with dual gooseneck and two focusing lenses is shown. Items are sold separately.

### Fiber Optics, Fiber-Lite™, Dolan-Jenner

Light Source, 150W EKE lamp, rheostat. **MI-150** .....\$700.00  
(FL150/115) Warranty: 3 years

Gooseneck, Dual, chrome plated, 23" long. **EEG2823M** .....\$450.00  
(FL150/05) Warranty: 1 year

Focusing Lens, each. **LH-759** .....\$80.00  
(FL150/50) Warranty: 1 year