

Pipets, Transfer, LDPE (low-density polyethylene)

Graduated, Nonsterile, 15 cm long

- ★ 5 mL capacity, 1 mL graduated, 25 drops/mL
 - Bx/500. **PTP-02**.....\$21.30
 - Cs/10 Bx. **PTP-02**.....\$196.00
- 5 mL capacity, thick-walled
 - Cs/10 Bx. **188BP1200**.....\$164.00
- 7 mL capacity, 3 mL graduated, 14 drops/mL
 - Bx/500. **PTP-01**.....\$21.30
 - Cs/10 Bx. **PTP-01**.....\$196.00
- Graduated, Sterile Pk/1 (individually wrapped), 15 cm long
 - 5 mL capacity, 1 mL graduated, Simport
 - Bx/400. **P200-521S**.....\$49.85



These mini, air displacement pipettors provide an economical and accurate method for pipetting fixed volumes. Pipettors feature a two-step plunger and are fully autoclavable.

Pipettors, Single Channel

Fixed Volume, pocket-size, Diamond® Jr.

- 5 µL, red. **3322-5**.....\$27.75
- 10 µL, pink. **3322-10**.....\$27.75
- 20 µL, green. **3322-20**.....\$27.75
- 50 µL, purple. **3322-50**.....\$27.75
- 100 µL, orange. **3322-100**.....\$27.75

Single channel, adjustable volume pipettors are fully autoclavable.



Adjustable Volume, Diamond® Apex™ (3-year warranty)

- 0.5-10 µL, 0.02 µL increments, red. **3351-10**.....\$198.00
- 5-50 µL, 0.1 µL increments, yellow. **3351-50**.....\$198.00
- 10-100 µL, 0.1 µL increments, blue. **3351-100**.....\$198.00
- 100-1000 µL, 1 µL increments, orange. **3351-1000**.....\$198.00



Pipettor Holders, Carousel, Flip & Grip™, Heathrow

- Blue. **HS120026**.....\$92.45
- Lime. **HS120056**.....\$92.45

Pitchers



Polypropylene




Stainless Steel

Pitchers

- Polypropylene, graduated (mL), tall form, tapered
 - 1000 mL, 50 mL graduations, each. **289910000**.....\$7.00
 - 2000 mL, 100 mL graduations, each. **289920000**.....\$8.05
 - 3000 mL, 100 mL graduations, each. **289930000**.....\$13.80
- Stainless Steel, graduated (mL/oz), tapered, Polar Ware
 - 1000 mL/32 oz, each. **T1063**.....\$73.20

Rockers, Shakers & Mixers



Rockers, Fixed Speed, Benchmark (2-year warranty)

- 2-D (seesaw), 12 x 12" platform w/flat mat
 - BenchBlotter™. **BR1000**.....\$543.00
 - Fixed speed of 12 rocks/minute
 - Platform angle is ±8°
 - Cold room & incubator safe (4-45°C)
 - Load capacity is 2.3 kg (5 lbs)
- 3-D (nutating)
 - Small, 10.5 x 7.5" platform w/dimpled mat
 - MiniMixer™. **B3D1020**.....\$449.00
 - Large, 12 x 12" platform w/dimpled mat
 - BioMixer™. **B3D1320**.....\$558.00
 - Fixed speed of 24 rocks/minute
 - Platform angle is ±20°
 - Cold room & incubator safe (4-45°C)
 - Load capacity is 1.8 kg (4 lbs)

Continued