



Rockers, Variable Speed

Benchmark (2-year warranty)

- 2-D (seesaw), 14 x 12" platform w/flat mat
 BenchRocker™ 2D. **BR2000**\$830.00
 Stacking Platform. **BR2000-STACK**.....\$155.00
 14 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP**\$55.00
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 3.6 kg (8 lbs)

- 3-D (nutating), 12 x 12" platform w/dimpled and flat mats
 BenchRocker™ 3D. **B3D2300**\$1,029.00
 Stacking Platform. **BR1000-STACK**.....\$142.00
 12 x 12" platform w/flat mat has 3" of clearance.
 Stacking Platform Extenders, Pk/4. **BR2000-SP**\$55.00
 Adds 1.25" to platform separation.

- Variable speed of 2-30 rocks/minute
- Platform tilt angle adjusts from 0-30°
- Cold room & incubator safe (4-45°C)
- Load capacity is 2.3 kg (5 lbs)



Rockers, Variable Speed

Reliable Scientific (5-year warranty), **USA**

- 2-D (seesaw)
 Standard, 11 x 14" platform w/rubber mat
 Single Platform. **55S1114**\$603.00
 Double Platform, 3¼" apart. **55D1114**.....\$665.00
 Large, 12 x 16" platform w/rubber mat
 Single Platform. **55S1216**\$647.00
 Double Platform, 3¼" apart. **55D1216**.....\$704.00

- Variable speed of 40-160 rocks/minute
- Platform tilt angle adjusts from 0-20°
- Cold room & incubator safe (0-65°C)
- Load capacity is 4.5 kg (10 lbs)



Rocker, Tube, 12 x 4.5" grooved mat, Benchmark. **M2100**\$405.00

- Fixed speed of 24 rocks/minute
- Platform tilt angle is ±30°
- Coldroom & incubator safe (4-45°C)
- Accommodates 1.5 to 50 mL tubes up to 2.3 kg (5 lbs)
- 2-year warranty

Orbital Shakers



Aerate and mix with these low to high capacity orbital shakers.

Orbital Shakers, Benchmark (2-year warranty)

- Low Speed, 14 x 12" platform w/flat rubber mat
 Orbi-Blotter™. **BT30**\$933.00

- Speed adjusts from 3-70 rpm
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 3 x 2 L, or 4 kg (9 lbs)

- Small, 9.5 x 8" platform w/flat rubber mat
 Orbi-Shaker™ JR. **BT302**\$1,126.00

- Speed adjusts from 50-300 rpm in 10 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 4 x 1 L, or 4 kg (9 lbs)
- Timer is programmable from 1 min - 99 hrs, 59 min



- Large, 13 x 12" platform w/flat rubber mat
 Orbi-Shaker™. **BT3001**\$1,730.00

- Speed adjusts from 30-300 rpm in 1 rpm increments
- Orbit diameter measures 19 mm
- Cold room & incubator safe (4-65°C)
- Load capacity is 9 x 500 mL, or 5 x 1 L, or 2 x 2 L, or 6.4 kg (14 lbs)
- Timer is programmable from 1 min - 48 hrs/continuous

Optional Magic Clamp™ platform & accessories sold separately - specify configuration

- Steel Platform, Magic Clamp™, 14 x 12". **BT3000-MR**....\$327.00

- Clamps, magnetic
 250 mL Flask. **H1000-MR-250**\$46.10
 500 mL Flask. **H1000-MR-500**\$46.10
 1000 mL Flask. **H1000-MR-1000**\$70.50