


Rotating Mixers



Tube Rotator, programmable, Benchmark (2-year warranty)
 Roto-Bot Mini. **R4040**\$679.00
 Roto-Bot. **R4045**\$1,019.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



Vortex Mixers, Benchmark (2-year warranty)
 Vornado™, mini vortexer, fixed 2800 rpm. **BV101**_*\$229.00
 *Specify **Gray**, **Blue** or **Purple cup head**
 BenchMixer™, single cup head, 200-3200 rpm. **BV1000**\$296.00
 Flat Head Attachment, 3". **BV1000-FLAT**\$44.75

Rulers, Meter Stick & Tape Measure



Rulers, 15 cm/6" long, scaled in both metric and English units
 Vinyl, flexible
 Transparent. **4461**.....\$0.95
 White. **4460**\$0.55
 Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065**.....\$6.25

Meter Stick, lacquered wood, plain ends. **605-050**\$4.90
 One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036**.....\$3.70
 One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



Sharps Containers, red polypropylene
 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269**\$9.90
 Dimensions: 12" x 4.75" x 11"(H)
 12 Gallon, hinged vertical entry lid, each. **440321391**\$30.75
 Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers



Specimen Containers
 Screw-Cap, polypropylene (PP)
 90 mL, sterile, Cs/400. **840175**\$150.00
 480 mL (16 oz), Cs/24. **59014510**\$71.95
 83 mm (3 1/4") opening
 960 mL (32 oz), Cs/12. **59014511**\$59.05
 111 mm (4 3/8") opening
 Snap-Cap, white HDPE bottom/white LDPE lid
 480 mL (16 oz), Cs/100. **271016W**\$86.20
 9.9 cm opening x 11.1 cm high
 1000 mL (34 oz), Cs/100. **271032W**\$129.00
 11.7 cm opening x 13.7 cm high

Spot Plate



Spot Plate, glazed white porcelain, 21 mm diameter wells
 12-Well, 115 x 95 x 7 mm, each. **3104187701**\$4.45