


Tubes & Vials



- Tubes, Centrifuge, Conical** bottom, sterile, white writing patch
RNase-/DNase-/Pyrogen-free
- 15 mL, PP (polypropylene) tube/HDPE cap, natural
Rack/25, 20 racks, Cs/500. **BP1501**\$116.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - 15 mL, PS (polystyrene)
Rack/50, 10 racks, Cs/500. **6274**\$246.00
Maximum RCF is 3,000 x g.
 - 50 mL, natural, polypropylene (PP) tube/HDPE cap
Bag/25, 20 bags, Cs/500. **BP5002**\$128.00
Maximum RCF is 12,000 x g. Frosted writing patch.
 - Bag/25, 20 bags, Cs/500. **6297**\$200.00
Maximum RCF is 20,000 x g
 - Rack/25, 20 racks, Cs/500. **BP5001**\$135.00
Maximum RCF is 12,000 x g. Frosted writing patch.

Self-Standing



- Tubes, Centrifuge, Self-Standing** (skirted)
- 30 mL
Nonsterile, PP tube/PP cap, Cs/500. **109220**\$171.00
Dimensions: 23(OD) x 90 mm
 - Sterile, PP tube/HDPE cap
Cs/500, 20 bags/25. **C572-1**\$156.00
Molded-in graduations from 5 to 30 mL in 5 mL increments. Dimensions: 25.3(OD) x 111 mm

Plastic culture tubes with dual-position caps.



Tubes, Culture,

- 12 x 75 mm (5 mL)
Polystyrene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110428**\$89.50
- Polystyrene, Non sterile, Cs/1000. **110410**\$62.90
- Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110438**\$89.70
- 13 x 100 mm (10 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1510**\$97.20
- 17 x 100 mm (16 mL)
Polypropylene, Sterile, dual-position caps
Cs/500, 20 sterile bags/25. **110178**\$143.00
- 16 x 125 mm (19 mL)
Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1515**\$157.00
- 16 x 150 mm (23 mL)
★ #1 Glass (Type 1, Class A borosilicate), rimless
Cs/1000. **1517**\$161.00

Culture Tube Caps, polypropylene
Up to 16 mm Tubes



- Blue, Cs/1000. **118152B**\$99.65
- Green, Cs/1000. **118152G**\$99.65
- Red, Cs/1000. **118152R**\$99.65
- Yellow, Cs/1000. **118152Y**\$99.65



RNase-/DNase-free

Tubes, Microcentrifuge, Economy, polypropylene

- 0.5 mL, 0.1 mL graduations, natural, Pk/1000. **HS4422**\$38.45
- 1.5 mL, 0.25 mL graduations, Heathrow
Natural, max RCF is 16,400 x g
Pk/500. **HS4323**\$18.70
Pk/1000. **HS4325**\$35.35
- Black, max RCF is 18,000 x g, Pk/500. **HS4323K**\$20.80
- 5 mL, natural, max RCF is 25,000 x g, Benchmark
Nonsterile, 4 bags/50, Pk/200. **C1005T5**\$65.40
Sterile, 10 bags/20, Pk/200. **C1005T5ST**\$71.45

Continued