

Metal Halide

60V/24W MR11, 300 hrs, Welch Allyn. **09500-U**\$270.00
(AL-1824, 8-61257)

Mercury Short Arc

50W Mercury Arc, 100 Hrs, Osram. **HBO50W/AC/L1**\$115.00
(78564, 78588, 8-B199, 38-16-19, 414140, 2055,
3040-0090, 500-137, SHL-50AC1, BA-001, 5000487)
50W Mercury Arc, 100 Hrs, Osram. **HBO50W/AC/L2**\$115.00
(31-31-65, SHL-50AC2, BA-002, 5000857)
50W Mercury Arc, 200 Hrs, Osram. **HBO50W/3**\$140.00
(BA-003)
100W Mercury Arc, 200 Hrs, Ushio. **USH-102D**\$109.00
(78589, 8-B192U, 500-138, 3800-18-4060/D,2087-001,
3040-0380, BA-005A, HBO100W/2, M/5000/133, 5000274)
100W Mercury Arc, 300 Hrs, Ushio. **USH-103D**\$135.00
(78590, 8-B1920, 380301-9350, USH-103OL, 5001329,
HBO103W/2)
200W Mercury Arc, 1000 hrs, Ushio. **USH-200DP**\$135.80
(HBO200W/2, 31-01-01, 649, 860014, 78580, 8-B190,
3800-18-4070, 38-16-15, 050-320, MA885,
91-00-50, 505-620, 78590, 5000276)

Xenon Short Arc

75W Xenon Arc, 400 hrs, Ushio. **UXL-75XE**\$195.00
(XBO75W/2, 8-B198, 910881-2013, 38-16-02, 500-139,
91-01-42, 3800-79-9190, 78682, 620-125, 903-050,
BA-013, 5000371)
150W Xenon Arc, Short Gap. **XBO150W/CROFR**\$260.00
(8-B195, BA-018, UXL-150S-0, 3801-22-2450,
500-008, 5000326)
180W Xenon Arc, Osram. **XBOR180W45COFR**\$450.00
(69183, 407.309)

Fluorescent

4W Fluorescent, Tubular, Eiko. **F4T5/D**\$7.00
(31-31-03, MA309/05, 3000108)
5W Fluorescent, G23, Sylvania. **CF5DS/841**\$11.70
(21278, 800-138, MA68, 48920-81, 20303)
9W Fluorescent, G23, Sylvania. **CF9DS/841**\$4.95
(20305)
22W Fluorescent, 8" OD Circle, Sylvania. **FC8T9/D/RS**\$19.00
(29980)
32W Fluorescent, 12" OD Circle, Sylvania. **FC12T9/D/RS**\$15.90
(20030)

Germicidal

4W Germicidal (UV-C), G5, Osram. **G4T5/OF**\$10.35
(23397)
15W Germicidal (UV-C), G13, Sylvania. **G15T8**\$17.00
(21612, 3000007)
25W Germicidal (UV-C), G13, Ushio. **G25T8**\$19.80
(3000008)
39W Germicidal, single-pin, Ushio. **G36T5L**\$35.80
(23443, 00206, 3000312)
65W Germicidal, single-pin, 155.4 cm, Higuichi. **G64T5L**\$83.00

Heat (Infrared)

250W, R40, shatter resistant, GE. **250BR40/HI/TC**\$14.00

* Recommend Bulb Removal/Installation Tool (shown below)



Remove and install prefocused/flanged bulbs safely with this simple tool. Push rubber-lined tube directly onto top of bulb and turn clockwise to install, or counterclockwise to remove.

Bulb Removal/Installation Tool. 3283.....\$16.00

Batteries

Batteries, button cell, 11.6 mm diameter
Silver Oxide, 1.55V, 5.4 mm high, Energizer. **357**\$4.00
(303, A76, AG13, EPX76, LR1154, LR44)
Alkaline, 1.5V, 3.05 mm high. **389**\$2.25
(LR1130, LR54, SR1130, AG10, 189)

Continued