





Microscope Filters, glass

- 29.8 mm Diameter
 - Blue, clear. **MA289/05**\$50.00
 - Blue Frosted. **MA289**\$50.00
 - Diffusion (frosted). **MA288**\$50.00
- 32 mm Diameter
 - Blue, clear. **923-160**\$15.00
 - Blue Frosted. **3236**\$7.50
 - Diffusion (frosted). **3235-W**\$7.50
 - Green, clear. **CM-0321-G**\$7.00
- 40 mm Diameter
 - Blue Frosted. **MA563/05**\$70.00
 - Diffusion (frosted). **MA566**\$70.00
- 45 mm Diameter
 - Blue, clear. **14671**\$10.40
 - Diffusion (frosted). **14673**\$35.00
 - Green Interference, 546/60, Nikon. **MBN11200**\$52.00
 - Neutral Density, Nikon
 - 50% Transmittance, ND2A. **MBN21802**\$42.00
 - 12.5% Transmittance, ND8A. **MBN21808**\$42.00



Renew and improve stage surface. Cleans easily.
 Light surface provides contrast for previewing specimens.
 Additionally, PTFE coating aids moving slides on the stage.

- Microscope Stage Laminate, USA. 16251**\$8.00
- 6.25 x 7.75" sheet provided with a #10 scalpel.
- White Teflon® (PTFE) on self-adhesive vinyl is 0.010" thick.



- Nosepiece Dust Plug, RMS thread (20.32 mm Dia.). 15420** \$0.50



Mark slides, dishes, flasks, and chambers with a 1.8 mm ring for later referral. Works on upright and inverted microscopes.

- Object Marker, RMS thread (20.32 mm Dia.)**
- Refill Ink, black, 10 mL, Shachihata. **22112**\$7.00
- Adapters
 - RMS/45 to CFI60, 15 mm long. **14111**\$43.00
 - RMS ID to 25 mm OD Thread, brass. **14115**\$52.00

Microscopes with Nikon CFI60 infinity-corrected (∞) objectives require a 14111 adapter.

Microscopes with Leica 25 mm thread/45 mm parfocal length objectives require a 14115 adapter.

Microscopes with 160 mm tube length objectives do not require an adapter.



Objectives, 160 mm tube length, DIN, Meiji, Motic, National
 Achromats, Motic/National

- 4X, N.A. 0.10. **704-160**\$32.00
- 10X, N.A. 0.25. **710-160**\$33.00
- 40X, N.A. 0.65. **740-160**\$42.00
- 100X Oil, N.A. 1.25. **799-160**\$69.00

Super High Contrast Achromats, National

- 4X, N.A. 0.10. **704-160ASC**\$34.00
- 10X, N.A. 0.25. **710-160ASC**\$43.00
- 40X, N.A. 0.65. **740-160ASC**\$79.00
- 100X Oil, N.A. 1.25. **799-160ASC**\$117.00

Semi-Plan Achromats, National (unless noted)

- 2X, N.A. 0.05, Accu-Scope. **3193SP-00**\$70.00
- 4X, N.A. 0.10. **704-160SP**\$39.00
- 10X, N.A. 0.25. **710-160SP**\$40.00
- 40X, N.A. 0.65. **740-160SP**\$51.00
- 50X Oil, N.A. 0.95, Meiji. **MA694**\$330.00
- 60X, Strain Free, N.A. 0.85. **MA645**\$360.00
- 100X Oil, N.A. 1.25. **799-160SP**\$82.00

Plan Achromats, National (unless noted)

- 4X, N.A. 0.10. **704-160P**\$40.00
- 10X, N.A. 0.25. **710-160P**\$65.00
- 20X, N.A. 0.40. **720-160P**\$147.00
- 40X, N.A. 0.65. **740-160P**\$123.00
- 50X Oil, N.A. 0.95, Accu-Scope. **3190PL-00**\$315.00
- 60X, N.A. 0.85. **760-160P**\$232.00
- 100X Oil, N.A. 1.25. **799-160P**\$197.00

Protect spare microscope objectives.



- Objective Case, DIN, polystyrene. 15601**\$0.50



Extender tubes do not contain lenses.

Objective Extender Tubes, RMS thread, brass, USA

- 10 mm Long. **14105**\$44.20
- 18 mm Long. **14106**\$42.75