


**Spade Tips**

- Angled (Kuehne), 11 cm (4¼")
- Stainless Steel, economy. **33-1041** .....\$4.15
- Teflon® Coated, Bel-Art. **379420000** .....\$35.60
- Straight, stainless steel, 10.5 cm, economy. **57.926** .....\$8.45



**Square Point Tips, straight**

- Teflon® Coated, 11 cm (4¼"), Bel-Art. **379430000** .....\$39.40

**Locking**



**Forceps, stainless steel**

- Self-Closing (cross-locking), squeeze handles to open tips
- Fine Point, smooth tips
- Economy, 13 cm (5⅞"). **300-066** .....\$4.95
- Peer-Vigor, 11 cm (4¼"), Swiss. **57.324** .....\$42.40
- Medium Point, serrated tips, 17 cm (6⅝"). **TWZ-948.00** .....\$3.95



**Slide-Locking, grid-serrated tips**


- Fine Point
- Economy, 14 cm (5½"). **57.0357A** .....\$9.00
- Dumont, 16 cm (6¼"). **TWZ-101.72** .....\$45.75
- Medium Point
- Economy, 16 cm (6¼"). **203** .....\$9.00
- Dumont, 16 cm (6¼"). **TWZ-101.76** .....\$43.60
- Snap-Closing, angled, fine points, 15.5 cm (6"). **3610** .....\$5.80
- serrated tips & handles, guide pin

**Tissue**



**Forceps, stainless steel**

- Adson-Brown, 9 x 9 side grasping teeth, 12 cm. **17-574** .....\$16.05
- Adson Tissue, 1 x 2 teeth, 12 cm (4¾"). **17-560** .....\$8.55
- Iris Tissue, 1 x 2 teeth
- Straight, economy, 10 cm (4"). **26-295** .....\$10.75
- Angled, economy, 11.5 cm (4½"). **19-350ANG** .....\$5.00
- Half Curve, standard, 10 cm (4"). **45-822** .....\$17.60





- Russian Tissue, circular serrated tips, 15.5 cm (6"). **2275** .....\$5.25
- Tissue, 1 x 2 teeth, 14 cm (5½"). **17-145** .....\$8.55

**Tongs**



- Tongs, Crucible, serrated jaws, stainless steel**
- 24.5 cm (9¾"). **320-101** .....\$11.65
- 51 cm (20"). **320-120** .....\$28.95

**Hemostats & Clamps**



**Hemostats, serrated jaws**


- Hartman Mosquito, stainless steel
- Straight
- Economy, 9 cm (3½"). **12-180S** .....\$6.05
- Standard, 9 cm (3½"). **20-352 (525-025)** .....\$21.30
- Curved
- Economy, 9 cm (3½"). **306-012 (14-030)** .....\$9.10
- Standard, 10 cm (4"). **20-360** .....\$21.70
- Halstead Mosquito, 12.5 cm (5"), stainless steel
- Straight
- Economy. **14-065** .....\$10.00
- Standard. **20-501** .....\$22.45
- Curved
- Economy. **14-070** .....\$10.00
- Standard. **20-519** .....\$22.85
- Kelly, stainless steel
- Straight
- Economy, 14.5 cm (5¾"). **306-016 (14-225)** .....\$8.20
- Standard, 14.5 cm (5¾"). **20-550 (525-225)** .....\$23.60
- Curved
- Economy, 14 cm (5½"). **306-017 (14-230)** .....\$9.90
- Standard, 14.5 cm (5¾"). **20-568** .....\$23.95
- Rochester-Pean, 20.5 cm (8"), stainless steel
- Straight
- Economy. **1920** .....\$6.50
- Standard. **21-977 (525-435)** .....\$33.95
- Curved
- Economy. **1930** .....\$6.50
- Standard. **21-979 (525-440)** .....\$33.95
- Polypropylene, 12.5 cm (5"). **HEM-300.00** .....\$1.95

Continued

**Hemostats (cont)**



Specimen Forceps, 25 cm (9<sup>3</sup>/<sub>4</sub>" ), serrated jaws. **300-500** ..... \$17.25



Sponge Forceps, Foerster, serrated jaws, stainless steel  
 Straight  
 17.5 cm (7"), economy. **14-660** ..... \$17.50  
 24 cm (9<sup>1</sup>/<sub>2</sub>" ), economy. **14-670** ..... \$21.10  
 Curved, 17.5 cm (7"), economy. **14-665** ..... \$17.50

Straight  300-120

Curved  300-121

**Clamps**


Dieffenbach Serrefine, stainless steel, economy  
 Straight, serrated jaws and handles  
 60 mm (2<sup>3</sup>/<sub>8</sub>" ). **300-120** ..... \$6.80  
 Curved, serrated jaws and handles  
 48 mm (1<sup>7</sup>/<sub>8</sub>" ). **12-248** ..... \$9.15  
 60 mm (2<sup>3</sup>/<sub>8</sub>" ). **300-121** ..... \$6.80

 Riveted handle opens for cleaning.

Johns Hopkins Bulldog, Roboz  
 Straight, 38 mm (1<sup>1</sup>/<sub>2</sub>" ). **RS-7440** ..... \$85.15

 RS-5452

Schwartz Vessel Clips, 1.7 x 8 mm jaws, 2.5 cm long, Germany  
 Straight Jaws, Roboz. **RS-5452** ..... \$34.00  
 Slightly Angled Jaws, Roboz. **RS-5457** ..... \$34.00



  
 8.6 cm shaft; 5 mm serrated jaw

Micro Alligator Forceps, stainless steel. **2769** ..... \$50.00

**Needle Holders & Wound Closure**

 Crile-Wood


 Olsen-Hegar

 Webster  


**Needle Holders, stainless steel**

Baumgartner, cross-serrated jaws, 14 cm (5<sup>1</sup>/<sub>2</sub>" ). **16-285S** ..... \$11.05  
 Castroviejo, cross-serrated jaws w/lock, 14 cm. **15-302** ..... \$50.00  
 Crile-Wood, cross-serrated jaws, 15 cm (6")  
 Economy, central groove. **16-293S** ..... \$13.45  
 Standard. **32-485** ..... \$33.75  
 German. **320-366** ..... \$60.80  
 Derf, 12 cm (4<sup>3</sup>/<sub>4</sub>" ), economy. **16-238** ..... \$9.70  
 serrated jaws with central groove  
 Olsen-Hegar, cross-serrated jaws with scissors, 14 cm (5<sup>1</sup>/<sub>2</sub>" )  
 Economy, central groove. **16-296** ..... \$14.30  
 Economy, No Groove. **16-296A** ..... \$14.45  
 Standard, central groove. **32-500 (555-160)** ..... \$35.25  
 Webster, smooth jaws, 12 cm (4<sup>3</sup>/<sub>4</sub>" )  
 Economy. **15-035** ..... \$17.45  
 Standard. **32-101 (555-035)** ..... \$30.75  
 Tungsten Carbide. **0817** ..... \$45.00  
 German. **320-101** ..... \$60.50

**Michel Wound Closure System**



Michel Clip Applier/Remover, double-ended. **3700** ..... \$13.75  
 11.5 cm (4<sup>1</sup>/<sub>2</sub>" ), stainless steel



Michel Wound Clips, stainless steel  
 11 mm, Bx/1000. **380-519** ..... \$74.00  
 14 mm, Bx/1000. **380-543** ..... \$89.00  
 16 mm, Bx/1000. **380-568** ..... \$97.00  
 18 mm, Bx/1000. **380-584** ..... \$98.00  
 22 mm, Bx/1000. **380-626** ..... \$116.00

**Reflex Wound Closure System**

CellPoint Scientific, USA

Reflex Clip Applier  
 For 7 mm Clips. **204-1000** ..... \$162.00  
 For 9 mm Clips. **202-1000** ..... \$162.00

Reflex Wound Clips, nonsterile, stainless steel  
 7 mm, 5 packs of 20 clips, Bx/100. **203-1000** ..... \$54.00  
 9 mm, 5 packs of 20 clips, Bx/100. **201-1000** ..... \$44.00

Reflex Clip Remover, fits 7 mm & 9 mm clips. **205-1000** ..... \$99.00

## Accessories

### Slide & Coverslip Handling

#### Forceps

##### Spade Tips



- Angled (Kuehne), 11 cm (4 1/4")  
Stainless Steel, economy. **33-1041** .....\$4.15
- Teflon® Coated, Bel-Art. **379420000** .....\$35.60
- Straight, stainless steel, 10.5 cm, economy. **57.926** .....\$8.45

##### Square Point Tips, straight



- Stainless Steel, 12 cm (4 3/4"), economy. **1244** .....\$2.25
- Teflon® Coated, 11 cm (4 1/4"), Bel-Art. **379430000** .....\$39.40



Minimize contamination during slide & coverslip handling. Will not scratch; corrosion and heat resistant.

- Hemostat, Polypropylene, 12.5 cm. **HEM-300.00** .....\$1.95

### Cytology Funnels

Single,  
Metal Clips  
Required



Double,  
All Plastic



White filter cards will accept up to 0.5 mL of sample, while tan filters will accept volumes up to 0.4 mL.

#### Cytology Funnels, Simport CytoSep™

##### Shandon (Thermo) Cytospin® CytoCentrifuges

Single. *Metal slide clips required*

- White Filter & Cap, tray/50, Cs/500. **M964-10FW** .....\$433.00  
(59 91 040, 110 25 48)
- White Filter & Cap, Pk/1, Cs/500. **M964-10FW1** .....\$574.00
- Tan Filter & Cap, tray/50, Cs/500. **M964-10FT** .....\$433.00
- Cytology Funnel Clips, metal, Pk/6. **M964B** .....\$178.00  
*Usable with the above three items*

Single. *All plastic - metal slide clips not required*

- White Filter & Cap, Cs/480. **M965-10FW** .....\$639.00  
(A787 100003)
- Tan Filter & Cap, Cs/480. **M965-10FT** .....\$639.00

Double. *All plastic - metal slide clips not required*

- White Filter & Cap, Cs/480. **M965-20FW** .....\$850.00  
(A787 100005)

### Immersion Oil & Bottles



Cargille immersion oils A, B, and NVH can be mixed with each other for custom viscosities.  
**USA**

#### Immersion Oil

**Type A**, low viscosity (150 cSt). *Short W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18101** .....\$10.00
- 20 mL in HDPE dispensing bottle. **18102** .....\$14.50
- 480 mL in amber bottle. **18104** .....\$104.00

**Type B**, high viscosity (1250 cSt). *Long W.D. optics*

- 7.5 mL in LDPE dispensing bottle. **18105** .....\$10.00
- 20 mL in HDPE dispensing bottle. **18106** .....\$14.50
- 480 mL in amber bottle. **18108** .....\$104.00

**Type NVH**, very high viscosity (21,000 cSt)

- 480 mL in amber bottle. **18113** .....\$67.00


**Type HF**, low fluorescence (700 cSt)

- 480 mL in amber bottle. **18118** .....\$159.00

**Type FF**, fluorescence free (170 cSt)

- 480 mL in amber bottle. **18124** .....\$154.00

A long nozzle, textured sides, and slim profile make this chemically inert, high-density polyethylene (HDPE) container ideal for dispensing immersion oil.



- Bottle, Dispensing**, HDPE, 20 mL, each. **18145** .....\$8.10