


Bone Cutters, stainless steel, economy

Liston

- 14 cm (5½"), 20 mm blades. **320-086**\$24.30
- 21.5 cm (8½"), 30 mm blades. **320-087 (1095)**.....\$29.85
- Scissor Style, 10 cm (4"), 18 mm blades. **320-025**\$15.15
straight, strong tips



Bone Nipper, 10 cm (4") long, 6 mm cutting edge. **B286**\$12.00
angled offset jaws




Bone Shears, Constable, 23 cm (9") long, 50 mm (2") cutting jaws

- Economy. **33-1001**\$45.00
- German. **965-620**\$160.00



Postmortem Hammer, hook handle, 25 cm (10") long

- Economy. **33-1032**\$35.70
- Standard. **74-043**\$80.00



16 mm (5/8") cutting blade;
9 cm (3½") long


Skull Breaker, Virchow, economy. **33-1011**\$42.40

Scalpels



Scalpels

Disposable, economy, sterile, gray plastic handle	bx/10	cs/50 bx
★ #10 carbon steel blade. 9410	\$10.90	480.00
#11 carbon steel blade. 9411	\$10.90	480.00
#22 carbon steel blade. 9422	\$10.90	480.00



Miniature, plastic thumb handle, nonsterile, England
#11P stainless steel blade, each. **2201111**\$1.15




Smooth Handle



Fluted Handle

Reusable, one-piece stainless steel, 1½" blade

- Smooth Handle, 15.5 cm (6"). **14-350**\$3.15
(1445)
- Fluted Handle, narrow blade, 15.5 cm (6"). **14-360**\$3.15
(1460)



Scalpel Blades, individually wrapped

- Carbon Steel, Economy, Bx/100\$21.75
 - #10 blade. **9310** #21 blade. **9321**
 - #11 blade. **9311** #22 blade. **9322**
 - #15 blade. **9315** #23 blade. **9323**
 - #20 blade. **9320** #24 blade. **9324**
- Carbon Steel, Germany, nonsterile
 - #60 blade, 25 strips of 4, Bx/100. **2200260**\$57.00
 - #70 blade, 25 strips of 4, Bx/100. **2200270**\$69.50
- Carbon Steel, Swann-Morton, nonsterile
 - #60 blade, Bx/50. **22002SM60**\$37.20
- Stainless Steel, Economy, Bx/100.....\$21.95
 - #10 blade. **36-010S** #21 blade. **052-610**
 - #11 blade. **052-602** #22 blade. **052-612**
 - #20 blade. **36-020S**
- Stainless Steel, Lance, Bx/100\$45.50
 - #11 blade. **92111** #15 blade. **92115**
- Stainless steel, Swann-Morton, sterile
 - #22 blade, Bx/100. **22001SM22**\$51.00
(RS-9861-19)



Safely remove and store up to 100 used blades in this disposable box.

Scalpel Blade Remover. 22048\$4.00



Scalpel Handles, reusable

- Stainless Steel, solid
 - #3 Style, accepts blades #9-16
 - Economy. **02-490** (36-405)\$3.45
 - Standard. **10-305**\$11.50
 - #4 Style, accepts blades #20-26
 - Economy. **02-492** (320-060, 36-410).....\$3.45
 - Standard. **10-310** (565-410)\$11.50
 - #7 Style, accepts blades #9-16
 - Economy. **02-494** (36-415)\$5.35
 - Standard. **10-315** (565-415)\$18.50
- Stainless Steel, Hollow
 - #8 Style, accepts blades #60 & #70. **HAN8**\$11.00

PVC Coated Stainless Steel. *Not autoclavable*

- #3 Style, accepts blades #10-17. **KNF-225.00**\$5.75
- Plastic, contoured to fit palm of hand
 - #5 Style, accepts blades #10-17. **KNF-221.00**\$5.00
 - #6 Style, accepts blades #20-26. **1605**\$5.00
 - #8 Style, accepts blades #60 & #70
- Bard-Parker™. **371080**\$44.95

Scissors



Scissors, Spring, ambidextrous, stainless steel

- Barraquer (butterfly), 6 cm long, economy. **3737**\$34.50
- 10 mm angled-on-flat blades

36205/5-254/5-855



M3720




Castroviejo, straight, sharp tips, economy

- 7 cm (2¾"), 5 mm blades. **36205**\$45.00
- Effective cutting edge is 3 mm; tip-to-tip is 3 mm.
- 8 cm (3"), 7 mm blades. **5-254**\$45.00
- Effective cutting edge is 6 mm; tip-to-tip is 4 mm.
- 8.5 cm (3¼"), 5 mm blades. **5-855**\$45.00
- Effective cutting edge is 4 mm; tip-to-tip is 3 mm.
- 11.5 cm (4½"), 17 mm blades. **M3720**\$26.00
- Effective cutting edge is 12 mm.



Uniband™

Lifetime guarantee against rust and defects.



Featherlite™

Uniband™, sharp tips, Laschal Surgical, **USA**

- Straight, 12 mm Blades. **LA-4**\$32.00
- Straight, 20 mm Blades. **LA-1**\$32.00
- Straight, 20 mm Serrated Blades. **LA-2**\$32.00
- Curved, 12 mm Blades. **LA-4C**\$32.00
- Curved, 20 mm Blades. **LA-1C**\$32.00

Featherlite™, Castroviejo, sharp tips, Laschal Surgical, **USA**

- Greater tactile sensitivity than Uniband™ scissors.
- Straight, 12 mm Blades. **MPF-4**\$80.80
- Curved, 12 mm blades. **MPF-4C**\$80.80

Continued

Continued