


Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**.....\$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**.....\$6.95
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**\$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**.....\$49.95



Tools not included.

Dissecting Cases
 6-Loop with snap fastener, black leatherette. **1257370**\$1.30
 12-Loop with snap fastener, black leatherette. **1257371**\$1.40
 13-Loop with zipper, brown leatherette. **36900**\$7.80
 Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
 Kit components are discounted 10%.*

Basic kit. 38501\$5.30
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551\$17.85
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)
 207/8 x 12¹³/₁₆ x 2 1/2". **AM20029**\$21.90
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**\$64.90



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic
 Black Enameled, Pk/100\$7.00
 #000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
 #00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
 #0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**
 Stainless Steel, Pk/100\$8.75
 #000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
 #00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
 #0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**\$55.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel
 Small, 7 x 3 x 1 1/2". **35-406**\$26.60
 Large, 10 x 4 x 2". **35-408**\$35.00



Rulers, 15 cm/6" long, scaled in both metric and English units
 Vinyl, flexible
 Transparent. **4461**\$0.95
 White. **4460**\$0.55
 Stainless Steel, 12 mm wide, **USA. 605-065**\$6.25



Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, **USA. 5-1450**\$419.00
 Glass Beads (replacement), 1.5 mm, Bg/12 oz. **5-1455**\$28.00



Continued



Safely incinerate microorganisms in 5-7 seconds at 815°C (1500°F) within this ceramic, stainless steel-encased chamber. Compact and flameless operation allows this unit to be used in anaerobic and aerobic chambers. Ideal for sterilizing loops and needles, glass tube/pipette mouths, and various glass and metal instruments. 1-year warranty.

Sterilizer, Infrared Heat, BactiZapper™. B1000\$445.00

6100A



Sterilizing Trays, plastic, autoclavable, PST, USA

B/L/M (base, lid & mat) included.

6 x 2½ x ¾" (15.2 x 6.4 x 1.9 cm). **2600A**\$23.00

(500252)



6½ x 4 x ¾" (16.2 x 10.2 x 1.9 cm). **4600A**\$33.00

10 x 6 x ¾" (25.4 x 15.2 x 1.9 cm). **6100A**\$67.00

(500254)



60081736

Instrument Trays

Polypropylene, deep w/flat cover, autoclavable

8½ x 6 x 2¾" (22 x 15 x 7 cm), Pk/2. **60081736**\$22.95



Shallow Tray



Deep Trays

Stainless Steel

Shallow (Mayo), **USA**

15½ x 10½ x ¾" (38.4 x 27 x 1.9 cm). **15F**\$51.65


Deep w/Dome Cover

8¾ x 3¼ x 1½" (22.2 x 8.4 x 3.8 cm). **35-309**\$31.25

Deep - Order trays and covers separately, **USA**

8⅞ x 5 x 2" (22.7 x 12.9 x 5.1 cm). **952**\$28.40

Cover, Fits 952 Tray, dome w/handle. **952D-2**\$44.90



This natural Arkansas stone outlasts and produces a better finish than a man-made stone. Mounted in a wooden box. Dimensions: 4 x 2 x ½"

Sharpening Stone, fine grit (hard). ABR-800.10\$31.95



Includes a plastic cleaning basket. Requires only 42 kHz of power. Dimensions: 6 x 3½ x 2"(D).

Ultrasonic Cleaner, 600 mL, mechanical timer. CLN-520.20\$45.95

Accessories

Cleaning Basket, fine mesh, 5 x 4 x 1¾". **CLN-650.10**\$14.95

Cleaning Solution, 1 Quart. **CLN-850.10**\$9.95

Biodegradable concentrate mixes 1:30.