


Trocar



Trocar, 4-in-1 set, 13 cm, stainless steel. **1596**.....\$17.50

Accessories/Maintenance/Storage



Cutting Boards, white high-density polyethylene
 10 x 6 x 1/2" (25.4 x 15.2 x 1.3 cm). **1054319**.....\$6.95
 18 x 12 x 1/2" (45.7 x 30.5 x 1.3 cm). **105CB1218**\$16.95



Demagnetize or magnetize parts and tools.

Demagnetizer/Magnetizer, 1 1/2 x 3" opening. **DEM-345.00**.....\$49.95



Tools not included.

Dissecting Cases
 6-Loop with snap fastener, black leatherette. **1257370**\$1.30
 12-Loop with snap fastener, black leatherette. **1257371**\$1.40
 13-Loop with zipper, brown leatherette. **36900**\$7.80
 Features red velveteen interior, protector flap, and ID pouch.

Dissecting Kits

*Inquire about customizing a kit for your surgical procedure.
 Kit components are discounted 10%.*

Basic kit. 38501\$5.30
 Contents: Dissecting Case, 6-loop (1257370); Forceps, broad point (EA-11); Pipet, plastic transfer, 5 mL (PTP-02); Pipet, plastic transfer, 7 mL (PTP-01); Scalpel, disposable, #10 blade (9410); Scissors, steel wire (22); Teasing Needle, straight (97-1301); Teasing Needle, angular (97-1302); and Ruler, 15 cm/6" (4460).

Mouse kit. 38551\$17.85
 Contents: Dissecting Case, 6-loop (1257370); Forceps, fine point, straight, 11.5 cm (300-130); Forceps, micro, #3 (TWZ-736.16); Forceps, micro, #7 (TWZ-736.28); Scissors, dissecting, straight, 30 mm blades (310-040); Scissors, operating, straight, S/S, 14 cm (12-055).

Dissecting Pans, stainless steel (without wax)
 207/8 x 12¹³/₁₆ x 2 1/2". **AM20029**\$21.90
 25 3/4 x 17 3/4 x 2 1/4". **AM68357**\$64.90



Insect Pins, 38 mm (1 1/2") long, nylon head, Czech Republic
 Black Enameled, Pk/100\$7.00
 #000, 0.25 mm dia. **355425** #1, 0.40 mm dia. **355440**
 #00, 0.30 mm dia. **355430** #2, 0.45 mm dia. **355445**
 #0, 0.35 mm dia. **355435** #3, 0.50 mm dia. **355450**
 Stainless Steel, Pk/100\$8.75
 #000, 0.25 mm dia. **355525** #1, 0.40 mm dia. **355540**
 #00, 0.30 mm dia. **355530** #2, 0.45 mm dia. **355545**
 #0, 0.35 mm dia. **355535** #3, 0.50 mm dia. **355550**



Instrument Cup, silicone, 3 1/2" opening, 2 1/2" deep. **150-1005**\$55.00



Instrument Storage Boxes, 10 x 4 x 2", 18-8 stainless steel
 Small, 7 x 3 x 1 1/2". **35-406**\$26.60
 Large, 10 x 4 x 2". **35-408**\$35.00



Rulers, 15 cm/6" long, scaled in both metric and English units
 Vinyl, flexible
 Transparent. **4461**\$0.95
 White. **4460**\$0.55
 Stainless Steel, 12 mm wide, **USA. 605-065**\$6.25




Rapidly decontaminate metal dissecting instruments with this dry heat sterilizer. Glass beads are heated to a constant 500°F (260°C) in a 2" diameter x 4" deep stainless steel well. Sterilize instrument tips in 15 seconds. 1-year warranty.

Sterilizer, Hot Bead, Germinator™ 500, **USA. 5-1450**\$419.00
 Glass Beads (replacement), 1.5 mm, Bg/12 oz. **5-1455**\$28.00

Continued

Rotating Mixers



Tube Rotator, programmable, Benchmark (2-year warranty)
 Roto-Bot Mini. **R4040**\$679.00
 Roto-Bot. **R4045**\$1,019.00

- Speed adjusts from 1-80 rpm
- Tube capacity (starter set)
 - R4040: 12 x 1.5 mL, 2 x 15 mL & 2 x 50 mL
 - R4045: 30 x 1.5 mL, 8 x 15 mL & 2 x 50 mL
- Timer is programmable from 1 sec – 99 hrs, 59 min / Cont.

Vortexers



Vortex Mixers, Benchmark (2-year warranty)
 Vornado™, mini vortexer, fixed 2800 rpm. **BV101**_*\$229.00
 *Specify **Gray**, **Blue** or **Purple cup head**
 BenchMixer™, single cup head, 200-3200 rpm. **BV1000**\$296.00
 Flat Head Attachment, 3". **BV1000-FLAT**\$44.75

Rulers, Meter Stick & Tape Measure



Rulers, 15 cm/6" long, scaled in both metric and English units
 Vinyl, flexible
 Transparent. **4461**.....\$0.95
 White. **4460**\$0.55
 Stainless Steel, 12 mm wide, 0.5 mm thick, **USA. 605-065**.....\$6.25

Meter Stick, lacquered wood, plain ends. **605-050**\$4.90
 One side is graduated on both edges in centimeters with 1 mm divisions, while the other side is graduated on one edge in inches with 1/8" divisions.

Tape Measure, 152.4 cm/60", reinforced fiberglass. **11036**.....\$3.70
 One side is graduated on the top edge in inches with 1/16" divisions and on the bottom edge in cm with 1 mm divisions. The opposite side is graduated in inches with 1/16" divisions.

Sharps Containers



Sharps Containers, red polypropylene
 5.4 Quart, horizontal entry lid, Prevent™, each. **4402269**\$9.90
 Dimensions: 12" x 4.75" x 11"(H)
 12 Gallon, hinged vertical entry lid, each. **440321391**\$30.75
 Dimensions: 18.25" x 12.75" x 18.75"(H)

Specimen Containers



Specimen Containers
 Screw-Cap, polypropylene (PP)
 90 mL, sterile, Cs/400. **840175**\$150.00
 480 mL (16 oz), Cs/24. **59014510**\$71.95
 83 mm (3 1/4") opening
 960 mL (32 oz), Cs/12. **59014511**\$59.05
 111 mm (4 3/8") opening
 Snap-Cap, white HDPE bottom/white LDPE lid
 480 mL (16 oz), Cs/100. **271016W**\$86.20
 9.9 cm opening x 11.1 cm high
 1000 mL (34 oz), Cs/100. **271032W**\$129.00
 11.7 cm opening x 13.7 cm high

Spot Plate



Spot Plate, glazed white porcelain, 21 mm diameter wells
 12-Well, 115 x 95 x 7 mm, each. **3104187701**\$4.45