

Slide & Cassette Pens, Marker & Scribes

Write a 0.4 mm stroke on glass slides with this permanent marker.
Pilot pens are ideal for cytologists and pathologists.



Pens, Pilot®, permanent, extra fine point

- Black, Bx/12. **44102**\$17.45
- ★ Blue, Bx/12. **44103**\$17.45
- Green, Bx/12. **44105**\$17.45
- Red, Bx/12. **44104**\$17.45



The 0.75 mm tip is ideal for use on small surfaces, e.g., histology cassettes, microscope slides and cryogenic vials. The permanent ink is waterproof and resistant to most laboratory chemicals and solvents.

Marker, Diamond PerfectMark™, black, Bx/10. **1000**\$29.80



Scribes, retractable

- Carbide-Tipped (sharpenable). **RSC**\$24.30
- Carbide Tip (replacement). **RCT-1**\$8.45
- Diamond-Tipped. **RSD**\$28.35
- Diamond Tip (replacement). **RDT-1**\$14.35

Slide Warmers




Evenly warm your slides from room temperature to 70°C. Black anodized surface provides excellent contrast.

Slide Warmers

- Small, 26 x 18 cm (10 x 7"). **XH-2002**\$337.00
- Large, 63 x 20 cm (25 x 8"). **XH-2001**\$443.00

Staining



Ten slide compartments, separated by 1/2" of space, are additionally isolated by eight barrier dividers when the lid is closed. Individual slides are suspended 1/2" above the chamber floor by four pedestal posts and four corner posts.

Staining/Moisture Chambers

- 10-Slide, polystyrene, **USA**

	<u>ea</u>	<u>pk/3</u>
Clear. 240-9000-010	\$54.95	154.00
Black. 240-9020-Z10	\$56.10	157.00



20-Slide, Simport

- ★ Clear Lid (PETG), black base (ABS), each. **M920-1**\$191.00
- Black Lid (ABS), black base (ABS), each. **M920-2**\$191.00



Staining Dishes

Glass (flint)

- 20-Slide, holds slides horizontally in metal rack
 - Approx. inside dimensions: 9.6 x 7.6 x 5.7 cm(H)
 - Set (dish, rack & lid). **SDS-20**\$29.90
 - Dish & Lid. **SD-20**\$18.20
 - Rack only, stainless steel. **SR-20**\$15.00
- 30-Slide, holds slides vertically in metal rack
 - Approx. inside dimensions: 9.6 x 8.5 x 8.3 cm(H)
 - Set (dish, rack & lid). **SDS-30**\$40.80
 - Dish & Lid. **SD-30**\$20.35
 - Rack only, stainless steel. **SR-30**\$25.00
- 60-Slide, holds slides vertically in metal rack
 - Approx. inside dimensions: 15.5 x 8.3 x 8.5 cm(H)
 - Set (dish, rack & lid). **SDS-60**\$57.10
 - Dish & Lid. **SD-60**\$31.20
 - Rack only, stainless steel. **SR-60**\$32.25

Continued



Plastic rack and dish resist alcohol and xylene; can be microwaved or frozen. Lightproof staining dish has a tight lid with a perforated slide partition for supporting a drying agent. Small size of dish saves on reagents, since only 150 mL is needed (holds 250 mL).


Plastic (polyoxymethylene [POM]), holds 25 slides
Set (dish, rack & lid). **0990-99**\$50.90



674



60062101



675

Staining Jars, hold slides (up to 76 x 26 mm) vertically

Coplin, holds 10 slides back-to-back	<u>ea</u>	<u>cs/12</u>
Glass w/Ground Glass Cover, USA. 673	\$13.10	146.00
Glass Cover, Fits 673 Jar. 673-LID	\$4.85	—
	<u>ea</u>	<u>cs/12</u>
Glass w/Screw-On Phenolic Cap, USA 674	\$13.10	146.00
Phenolic Cap, Fits 674 Jar, black. 674-LID	\$4.95	—
	<u>ea</u>	<u>cs/12</u>
PP w/Screw-On Cap. 60062101	\$4.90	47.40
Polypropylene jar is unbreakable, microwave proof, stain resistant, and chemically inert. Cap is domed.		
	<u>ea</u>	<u>cs/12</u>
Wide Top, Glass w/Square Lid, USA 675	\$17.65	193.00
Holds 16 slides back-to-back		
Glass Lid, Fits 675 Jar. 675-LID	\$5.40	—



12-Slide

Move slides, without forceps, between jars linked together in a workstation. Swing out the rack's integral lid for easy accessibility to slides. Stable, wide-based jar requires only 80 mL of reagent. Polyoxymethylene (POM) resists alcohol and xylene; can be microwaved or frozen.

Wide Top, Plastic (POM)
12-Slide, EasyDip™, Simport

Jar only, Blue, Cs/6. M900-12B	\$66.55
Jar only, Green, Cs/6. M900-12G	\$66.55
Jar only, Pink, Cs/6. M900-12P	\$66.55
Jar only, White, Cs/6. M900-12W	\$66.55
Jar only, Yellow, Cs/6. M900-12Y	\$66.55
Rack/Lid only, gray, Cs/6. M905-12DGY	\$73.55



Anodized aluminum staining jar holder

(Jars sold separately)
Staining Jar Holder, 6-Place. **M906**.....\$81.45




4501030



4501042

24-Slide, light green, microwaveable, acetal
Jar & Lid. **4501030**\$11.55
Rack only, gray. **4501042**\$21.45
Staining Jar Station, 12-Place. **4501011**\$221.00
Includes: metal frame, (12) 24-slide dishes w/ lids (4501030) & (1) rack (4501042).
Dimensions: 29(L) x 5(W) x 5"(H)

Tissue (Embedding) & Biopsy Cassettes




Base Molds, Plastic, disposable, Simport


15 x 15 x 5 mm, Cs/1000. M475-2	\$152.00
24 x 24 x 5 mm, Cs/1000. M475-3	\$152.00
30 x 24 x 5 mm, Cs/1000. M475-4	\$152.00
37 x 24 x 5 mm, Cs/1000. M475-5	\$152.00
37 x 24 x 10 mm, Cs/500. M475-10	\$78.35

Designed for Simport Macrosette™.

Solvent resistant,
0.2 x 0.2 mm mesh.



Nylon




Paper

Biopsy Bags
Nylon, 30 x 45 mm, 5 bags/200, Pk/1000. **4502011**\$235.00
Paper (Tea Bags), 2¼ x 2¼", Pk/500. **56100**\$46.50

M476-1

Made of 60 PPI cellular polyester urethane foam. Use rectangles with cassettes, circles with capsules.



M476-2

Biopsy Pads
Rectangular
Black Foam, 30.2 x 25.4 mm, Bg/1000. **339014**\$42.00
Blue Foam, 30.2 x 25.4 x 2 mm, Cs/10,000. **M476-1**\$217.00
Circular, Simport
Blue Foam, 25.4 x 2.7 mm, Cs/10,000. **M476-2**\$266.00