

Cassette Storage



Box, Matchbox Style (Pill Box), cardboard, white
 2 1/8 x 1 1/8 x 5/8", Pk/72. **20047**\$35.45

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**\$121.00
 Each box can store approximately 540 paraffin blocks.

Store 3,150 cassettes or 1,950 blocks
 in 4 drawers, each divided into 12 rows.
 Cabinets can be stacked up to 8 high.
 Dimensions: 49.5 x 38 x 21 cm
 (19 1/2 x 15 x 8 1/4")



Tissue Cassette Cabinet, Steel, 4-drawer. **485400**.....\$386.00

Store up to 2100 cassettes in 6 drawers.
 Each drawer includes sponge block
 spacers and an identification label.
 Cabinets can be stacked up to 4 high.
 Dimensions: 44.5 x 24 x 39.5 cm(H)
 (17 1/2 x 9 1/2 x 15 1/2")



Tissue Cassette Cabinet, ABS. **97-0600**.....\$186.00




Store 250 cassettes or 165 embedding
 rings in each stackable drawer.
 High-impact polystyrene resists
 temperatures from -80°C to +80°C.
 Drawer dimensions: 40.5 x 23 x 5.1 cm.

Simport

Tissue Cassette Drawers, Polystyrene
 Pk/2. **M495-2**.....\$75.60
 Pk/6. **M495-6**.....\$172.00

Brushes



Brushes
 Flat Bristles (paint style)
 Camel Hair, 1" wide, 6 1/2" total length. **BR665B (LB-37)**\$6.65
 Round Bristles, camel hair
 3/4" Diameter, 8 1/2" total length. **BR661 (LB-30)**\$11.15

Microtome Blades



Microtome Blades

Accu-Edge®, Sakura Finetek
 Low Profile, 80 x 8.3 mm, Pk/50. **4689**\$196.00
 High Profile, 75.7 x 14 mm, Pk/50. **4685**\$208.00
 Encore™, PTFE coated, **USA**
 Low Profile, 83 x 8 mm, Pk/50. **7380**\$84.00
 High Profile, 83 x 14 mm, Pk/50. **7390**\$84.00

Scrapers

Remove paraffin, ice, etc., safely
 from your laboratory equipment
 with these plastic scrapers.



Scrapers, plastic

5 cm (2") Wide Blade. **311-245**.....\$3.70
 10 cm (4") Wide Blade. **311-250**.....\$4.70
 15 cm (6") Wide Blade. **311-255**.....\$5.70

Tissue Flotation Bath



Tissue Flotation Bath. XH-1001.....\$307.00

Borosilicate glass bottle;
blue PP cap (GL45)



Bottles, Media, graduated, autoclavable

100 mL, Cs/10. 600BM0100	\$90.90
250 mL, Cs/10. 600BM0250	\$103.00
500 mL, Cs/10. 600BM0500	\$118.00
1000 mL, Cs/10. 600BM1000	\$168.00
Replacement Items	
Bottle Caps, GL45, blue, Pk/10. B3000-CAP	\$44.15
Bottle Drip-Free Rings, GL45, Pk/10. B3000-RIN	\$23.90



Globally Harmonized System (GHS) wash bottles are vented to prevent solvent drips. Polyethylene bottles (LDPE) with polypropylene (PP) caps.

Bottles, Wash, GHS Labeled, 500 mL (empty)

Acetone, Pk/6. 600370051	\$39.25
Bleach, Pk/6. 600370059	\$39.25
Deionized Water, Pk/4. 124160003	\$39.25
Ethanol, Pk/4. 600370052	\$39.25
Methanol, Pk/6. 600370053	\$39.25
Water, Pk/6. 600370061	\$39.25

Bottles, Wash, Safety Labeled, 500 mL (empty)

4-Color diamond, LDPE bottle/PP cap	
Soap, Pk/4. 117160014	\$39.25
Universal (label your own bottle), Pk/4. 117160009	\$39.25↓



Bottles, Wash, Integral Spout, LDPE bottle/PP cap

250 mL (8 oz), Cs/10. 601633	\$70.70
500 mL (16 oz), Cs/10. 601634	\$77.90↓



Bottles, Narrow-Mouth, HDPE bottle/polypropylene cap

125 mL, Pk/12. 106200015	\$24.50
250 mL, Pk/12. 106200016	\$35.20
500 mL, 28 mm cap, Pk/12. 106200007	\$40.95
1000 mL, Pk/6. 106200008	\$32.30



Bottles, Wide-Mouth

HDPE bottle/polypropylene cap	
125 mL, 38 mm cap, Pk/12. 60033407	\$24.50
250 mL, Pk/12. 106250006	\$35.20
500 mL, 53 mm cap, Pk/12. 60033409	\$40.95
1000 mL, Pk/6. 60033410	\$34.55
Polypropylene bottle and cap, autoclavable	
125 mL, 38 mm cap, Pk/12. 60033307	\$24.20
250 mL, 43 mm cap, Pk/12. 60033308	\$35.20
500 mL, 53 mm cap, Pk/12. 106320007	\$40.70
1000 mL, 63 mm cap, Pk/6. 60033310	\$35.95



Bottle Top Dispensers

0.5-5 mL. 600BTDR2	\$293.00
5-60 mL. 600BTDR5	\$349.00

Boxes



Matchbox Style

Box, Matchbox Style (Pill Box), cardboard, white

2 1/8 x 1 1/8 x 5/8", Pk/72. 20047	\$35.45
---	---------

Box, Pizza Box Style, 14 x 14 x 2", white, Pk/50. **8342**

Each cardboard box can store approximately 540 paraffin blocks.



246530002

246530001



132040001

Boxes, Glass Disposal, cardboard, polyethylene bag included

Bench Model, 8 x 8 x 10", Pk/6. 246530002	\$116.00
Floor Model, 12 x 12 x 27", Pk/6. 246530001	\$160.00
Glass Disposal Box Cover, reusable, PP. 132040001	\$52.65
<i>Replacement item for the floor model. Nonrefundable.</i>	
<i>Polypropylene cover is autoclavable.</i>	

Continued